



PRODUCT CATALOGUE

# SMALL ANIMALS



# 

DIAGNOSTIC	C
Otoscopes	8-10
Videootoscopes	10-12
Telescopes	13
Face shields	13
Sthetoscopes	14-19
Thermometers	20-21
Laryngoscopes	24
ECG	22-23
Blood Pressure Devices	25-26



# NJECTION AND INFUSION

Syringes	30
IV catheters	30-31
Needles	31
Infusion equipement	32
Stands	33-34
Compresion bands	34

35
38-39
40
40
40
41
41
41
42
43-44
44
45
46
47
48
49-50



51
54
55
56
57-58
59-64
65
66
66
67
67
68
68
68
69
69
69
70-76

# DRESSING MATERIALS

SUTURE PGA sutures NYLON sutures

PDO sutures

Suture needles Skin stapler

SURGERY

Pean Rochester Needle forceps

Forceps

Clamps

Collars

Clippers Lasers

Claw pincers Catheters for dogs

Catheters for cats

Weight scales

Absorbable sutures in casette

Absorbable sutures in casette

Disposable skin staple remover

Non-absorbable sutures in casettes 96-107

Adhesive tape	79
Viscose knit bands	79
Elastic bands	79-80
Wadding and rolls	81-82
Compresses	82
Sterile gauze	83
Plaster bandage	83
Synthetic plaster lining	84
Bands	84
Paw stabiliser	84

DIAGNOSTIC IMAGING	
Protection	1
Positioning	
Generator Ecotron	
	5

Generator Ecotron	159
Intermediate radiography system	160-165
Direct radiography	166-171
Dental radiography	172
USG	173-179
Endoscopes	180-182

ANAESTHETICS	
Anaesthetic machines	186-187
Ambu bags	188
Breathing bags	188
Anaesthetic masks	189
Endotracheal tubes	189-190
Monitors	191
Pumps	192

DENTISTRY	193
Stomatology units	196-198
Scalers	199, 210-2
Dental luxators	200
Bone ripper	200-201
Dental pincers	202
Extractor	202
Dental forceps	203-207
Prophy paste	208
Prophy angle	208
Doxiprobe	208-209
Dental mirror	211
Mouth gags	211-215
Dental model	216



Little lights	220-221
Lights	222-224
Instrument tables	225-226
Medical cabinets	226-227
Stools	228-229
Tables	230-233
Cages and accessories	234-238

# 156 157-159

101

a synthese

Standards stra



Weitlander wound dilator 139 Scalpel accessories 139 140-141 Electrocauters Adhesive cloth spray 141 Electrocauter accessories Operating field 142 Absorbent lining

CARE AND MEDICAL 143 PRODUCTS

88-89

92-95

108-122

128-133 134-135

135

136

146

147

147

148-149

149-151

137-138



# polnet.pl

# Diagnostic

# Videootoscope MD SCOPE

# Littmann Elektronic Stethoscope





MORE INFORMATION ON PAGE 9

01

MORE INFORMATION ON PAGE 15

### OTOSCOPES

01

### Videoskope MD-Scope

- invaluable for everyday ear nose and throat use;
  easily adjustable Neutral 2.5V vacuum light (made in Germany);
- | quality 3 times magnifying lens;
- ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations;
- supplied with 3 autoclavable ear funnels in a colour coordinated zip bag. All latex free;
- | require 2 x C batteries (not included).

Cat. No. 31442



### Xenon-halogen diagnostic set: otoscope/ophthalmoscope (plastic cannulae)

### | xenon-halogen light;

otoscope + ophthalmoscope + handle + 2 heads;
5 plastic cannulae (2.5, 3.5, 4.5, and 5.5-9 mm) - set 31451;
metal spatula providing light during throat examination;
shape of metal nasal speculum fits fingers.

Cat. No. 31424



### Xenon-halogen otoscope

xenon-halogen light;

- otoscope + handle and head;
- 5 plastic cannulae (2.5, 3.5, 4.5, 5.5, and 9 mm
- in diameter) 31451;
- | 3 plastic cannulae (4, 5, and 7 mm in diameter) 31447.

Cat. No. 31422



Xenon-halogen diagnostic set: otoscope/ophthalmoscope (metal cannulae)

| xenon-halogen light;

- | otoscope and ophthalmoscope + handle and 2 heads;
- | 3 metal cannulae (4, 5, and 6 mm in diameter);
- | 2 metal specula (4.5-7 mm in diameter).

Cat. No. 31425



### ACCESSORIES

Cat. No.	Name
31451	Set of 5 plastic cannulae (2.5, 3.5, 4.5, 5.5, 9 mm in diameter)
31447	Set of 3 plastic cannulae (4, 5, 7 mm in diameter)
31453	Set of 3 metal cannulae + 2 metal specula (cannulae:4, 5, 6 mm in diameter 2 metal specula (4.5-7 mm in diameter).
31454	Nasal speculum
31427	Xenon-halogen lamp for otoscopes
31428	Standard otoscope lamp
31449	Lampa otoskopowa standardowa

# Diagnostic



OTOSCOPES









### VIDEOOTOSCOPES

### Videootoscope MD SCOPE

videootoscope with image freezing function;
2.4 inch colour digital screen;
image resolutions 320 x 240 pixels;
ending made of stainless steel (5 mm in diameter);
sets of 2 alkaline batteries;
external video port .

### The set includes:

- | 2 plastic camera covers;
- 24 disposable specula;
- | 1-1.5 m video wire;
- standard speculum, 45 mm, for ear examination (M5101-002);
- | speculum, 75 mm, for nose examination (M5101-012).

Cat. No. 270175S



### MD-Scope videootoscope – Deluxe 2

### **Product description**

- | videootoscope with image freezing function;
- 2.4 inch colour digital screen;
- A clear and sharp image can be acquired and displayed momentarily on a builtin screen. Just press a single button and wait for 3 seconds;
- | image resolutions 320 x 240 pixels;
- | ending made of stainless steel (5 mm in diameter);
- sets of 2 alkaline batteries;
- external video port;

### Kit contents

- | videootoscope;
- | 24 disposable specula;
- RCA and USB cable;
- | 1 standard speculum, 45 mm, for ear examination (M5101-002);
- | 1 flexible probe 150 mm;
- carrying case 25.5 x 16 x 7 cm;
- | lightweight: 218 g;
- size: 7.5 x 6.3 x 21.2 cm.
- Cat. No. **32161**



### MD-Scope videootoscope – Elite pack

### Product description

- | videootoscope with image freezing function;
- | 2.4 inch colour digital screen;
- A clear and sharp image can be acquired and displayed momentarily on a builtin screen. Just press a single button and wait for 3 seconds;
- image resolutions 320 x 240 pixels;
- ending made of stainless steel (5 mm in diameter);
- sets of 2 alkaline batteries;
- external video port.

### Kit contents

- videootoscope;
- | 24 disposable specula;
- RCA video cable;
- | 1 standard speculum, 45 mm, for ear examination (M5101-002);
- carrying case 25.5 x 16 x 7 cm;
- | lightweight: 218 g;
- | size: 7.5 x 6.3 x 21.2 cm.

Cat. No. 32160

### MD-Scope videootoscope - Vet 3 probes

### Product description

- | videootoscope with image freezing function;
- | 2.4 inch colour digital screen;
- A clear and sharp image can be acquired and displayed momentarily on a builtin screen. Just press a single button and wait for 3 seconds;
- image resolutions 320 x 240 pixels;
- | ending made of stainless steel (5 mm in diameter);
- | sets of 2 alkaline batteries;
- external video port.

### Kit contents

- videootoscope;
- | 24 disposable specula;
- | RCA and USB cable;
- 3 flexible probe 75 mm, 150 mm, 300 mm;
- Image dongle;
- carrying case 25.5 x 16 x 7 cm;
- lightweight: 218 g;
- | size: 7.5 x 6.3 x 21.2 cm.

# Diagnostic







### VIDEOOTOSCOPES

### MD-Scope videootoscope – DeLuxe 4

### Product description

- | videootoscope with image freezing function;
- | 2.4 inch colour digital screen;
- A clear and sharp image can be acquired and displayed momentarily on a builtin screen. Just press a single button and wait for 3 seconds;
- image resolutions 320 x 240 pixels;
- ending made of stainless steel (5 mm in diameter);
- sets of 2 alkaline batteries;
- external video port.

### Kit contents

videootoscope;
24 disposable specula;
RCA and USB cable
1 standard speculum, 45 mm, for ear examination (M5101-002);

| 3 flexible probe 75 mm, 150 mm , 300 mm;

| Image dongle;

- carrying case 25.5 x 16 x 7 cm;
- lightweight: 218 g;
- | size: 7.5 x 6.3 x 21.2 cm.

Nr kat. 32162

### ACCESSORIES

Cat. No.	Name
32165	Standard speculum, 45 mm, for ear examination
32176	USB-memory key: compatible with Mac and OC computers
32176	Speculum, 75 mm, for nose examination
32168	Flexible speculum, 150 mm, for throat examination
32169	Long, flexible probe, 300 mm
32172	Disposable speculum covers 5 mm, 250 pcs.



# 

### NIKE GALILEIAN TELESCOPE 3x 35 cm + PROJECTOR

- | IUKO VISION MODULAR SYSTEM: unique modular system for all your needs;
- you can combine elements to protect, magnify, light up.
   Any component can be used indipendently or combined with one or more;
- | nIKE Galilean telescope 3X (35cm) + projector;
- this system includes Galilean telescope (CALLIOPE) and LED light projector (IRIDE);
- supplied in a case with accessories.

Nr kat. 30870

### NIKE GALILEIAN TELESCOPE 3x 45 cm + PROJECTOR

- IUKO VISION MODULAR SYSTEM: unique modular system for all your needs;
- you can combine elements to protect, magnify, light up.
   Any component can be used indipendently or combined with one or more;
- | nIKE Galilean telescope 3X (45cm) + projector;
- this system includes Galilean telescope (CALLIOPE) and LED light projector (IRIDE);
- | supplied in a case with accessories.

Nr kat. 30871

### TETI protective shield

- comfortable, ultra light stylish protective shield;
- | it provides high protection from squirts and sprays;
- | folding shield behind the nape;
- | adjustable size, fits majority of healthcare professionals;
- extra wide visor allows to fit over glasses and most of magnification loupes.;
- supplied with 2 face transparent shields and 2 padding with velcro;
- super soft padding with velcro fixing for confortable use;shield thickness 0.5 mm;
- vertificated EN166;
- made in Italy.

12

# Diagnostic



TELESCOPES/GUARDS







### TELESCOPES/STETHOSCOPES

### CALLIOPE GALILEIAN TELESCOPE 3x - dental

- | IUKO VISION MODULAR SYSTEM: unique modular system for all your needs;
- you can combine elements to protect, magnify, light up.
   Any component can be used indipendently or combined with one or more;
- CALLIOPE Galilean telescope 3X dental (working distance 35 cm);
- magnifying optical system with Galilean telescopes;waterproof telescope;
- compatible with protective shield TETI;
- supplied in a case with accessories .

### Nr kat. 30864

### CALLIOPE GALILEIAN TELESCOPE 3x – surgical



3

| IUKO VISION MODULAR SYSTEM: unique modular system for all your needs;

you can combine elements to protect, magnify, light up. Any component can be used indipendently or combined with one or more;

- CALLIOPE Galilean telescope 3X surgical (working distance 45 cm)
- $\cdot$  magnifying optical system with Galilean telescopes
- waterproof telescope,
- · compatible with protective shield TETI,
- supplied in a case with accessories.

### Nr kat. 30865

### JOTARAP stethoscope

- replaceable chromic tips;ensures universal acoustics in both large and small
- animals;
- | detects wide bands of cardiac and pulmonary sounds.

### Cat. No. 32580

### ACCESSORIES


Cat. No.	Name
32590	A set of JOTARAP replacement parts

Cat. No.	JOTARAP stethoscope tube colours
32579	purple
32580	black
32581	blue
32582	orange
32583	green
32584	yellow
32585	light blue
32586	fuchsia
32587	light pink
32588	sea green
32589	black



### Stethoscope Littmann Classic III

- high acoustic sensitivity;
- | dual single-piece tunable diaphragms (pressure sensitive);
- versatile adult and pediatric two-sides chestpiece;
- small-side diaphragm prevents debris in open bell;
- snap tight soft-sealing eartips;
- | next generation longer-life tubing;
- small side converts to open bell;
- lenght: 69 cm;
- head Ø: 43/33 mm.



Cat. No.	Colour sample	Tube colour
32292	0	raspberry
32295		crimson
32290	()	black
32291	()	navy blue
32301		dark green
32299	0	grey
32293		burgundy
32296	0	teal
32303	0	blue
32302		orange
32304	0	pink
32298	0	light green
32300	0	cyan
32297	0	peach
32305	0	pine green

# Diagnostic



STETHOSCOPES





STETHOSCOPES

### Stethoscope Littmann Classic II Pediatric

| intended for small patients, an open funnel on one side and a traditional membrane on the other;

- excellent acoustics;
- 7 points on a scale of 1-10 for acoustic quality;
- | optimal sound transmission due to anatomically designed;
- lyre;
- length: 71 cm;
- head diameter: 33/25 mm.

Nr kat. 32381-32387; 32404-32408



### Stethoscope cardiology

Nr kat. 32606-32613; 32631-32638; 32646-32649; 32371-32373; 32429-32439; 32443-32444; 32472-32475; 32359-32363; 32330-32332





CARDIOLOGY III

- | intended for small patients, an open funnel on one side and a traditional membrane on the other;
- excellent acoustics;
- 7 points on a scale of 1-10 for acoustic quality;
- | optimal sound transmission due to anatomically designed; lyre;
- length: 71 cm;
- head diameter: 27/19 mm.

Nr kat. 32376-32380; 32410-32412





MASTER CARDIOLOGY

# Diagnostic




STETHOSCOPES

### STETHOSCOPE ACCESSORIES

### Diaphragm+retaining ring

big one; | black; box of 5.

Cat. No. 32450



big one; gray; box of 5.

Cat. No. 32464





| 10 pcs in package; colour: black.

Cat. No. 32457, 32459




### Earpiece

| 10 pairs in package; colour: black.

Cat. No. 32451




- | provides exceptional ausculting;
- 85% noise reduction;
- | perfect acoustic thightness;
- | easily identifies difficult cardiac tones, abnormal pulmonary sounds, valve regurgitation;
- | the sound sensor is located at the bottom of the tube, so the sound is more natural on the ear;
- designed to visualise sounds;
- can replay anytime and anywhere;
- | sounds can be stored or transferred for further analysis;
- Bluethooth ready;
- can be connected to a computer;
- prints and plays sounds in the form of phonocardiograms.



# Diagnostic

01

### STETHOSCOPES / BLOOD PRESSURE METER



### THERMOMETERS

### Digital thermometer VisioFocus

- non-contact, able to take temperature on pet's eyes, rectum area, auricle or on the gum;
- | needs few seconds to measure temperature of pets;
- | automatically adjusts to the surrounding conditions;
- | body temperature mesuring range: 32 -43,5 ℃, general temperature measuring range: 1-80℃;
- | suitable for various animal species;
- | supplied by 4 AAA bateries ( not included);

Nr kat. 06620



### Galaxy digital thermometer

| waterproof;| large LCD screen;| measurement in 10 s;| flexible tip.

Cat. No. 25609



### Digiflash digital thermometer

- | a thermometer for large and small animals;
- | polycarbonate packaging;
- | measurement in 10 s;
- | temperatures from 32°C to 42°C;
- | sound signal at the end of measurement;
- | battery charge indication;
- turns off automatically;
- waterproof;
- shock resistant.

Cat. No. 00196000

### Vec Rectal digital thermometer

- | measurement in 8 s;
- | ergonomy suitable for animals;
- | often selected for pigs and cattle;
- | also suitable for other animal species.

Cat. No. 25588

# Diagnostic



### THERMOMETERS





### ECG MACHINES

### ECG 1-channel 100 G

- large LCD display waveform and data;
  high-resolution thermal printing array system;
  f requency response is as high as 150 Hz;
  high resolution digital filter inhibits baseline drift without
  affecting ECG waveform;
  echargeable battery supporting more than 50 patient
  examinations;
  12 leads;
  size : 315 x 215 x 77 mm;
- weight 1,5 kg

### INCLUDED:

- ECG cable 1 set;
  limb electrode clamp 1 set;
  chest electrodes 1 set;
  ground wire;
  power connection;
- | paper roll.

Nr kat. **33220** 

### ECG- 3 channel 300 G

- 3,5 " TFTcolour LCD screen;
- | thermal printing system;
- | frequency response is as high as 150 Hz;
- high resolution digital filter inhibits baseline drift without affecting ECG waveform;
- rechargeable battery supporting more than 150 patient examinations;
- | built-in Database and memory that can store 1000 cases;
- auto-measurement, auto-analusis, auto-interpretation;
- | 12 leads simultaneous display;
- size : 315 x 215 x 77 mm;
- weight 1,6 kg

### INCLUDED:

- ECG cable 1 set;
- | limb electrode clamp 1 set;
- chest electrodes 1 set;
- ground wire;
- | power connection;
- paper roll.

Nr kat. 33221



### VE-100 ECG machine

- | 1 ECG channel;
- LCD screen;
- | configured parameters;
- optimal printing position;
- | built-in battery for 3 h of use;
- digital filter.

Cat. No. 33305

### VE-300 ECG machine

- 3 ECG channels;
- LCD screen;
- | configured parameters;
- | optimal printing position;
- | built-in battery for 3 h of use;
- digital filter;
- | automatic parameters of measurement and analyses;
- can be connected to a computer or printer with a USB cable.

Cat. No. 33306

### ACCESSORIES:

| roll of thermical paper

| 2200 mA batteries;

spare wire.

Cat. No.		Dimensions
32963	Rolls for model 33305	50 x 25 m
32964	Rolls for model 33306	80 x 25 m
33318	Reusable lithium batteries	
33319	ECG cable	

# Diagnostic



ECG MACHINES





### Laryngoscope

adult set (4 blades Mc INTOSH 1-2-3-4);

- great features of the unique GIMA blades and handles;standard 2.5 V vacuum bulb;
- | light ergonomic, satin finish handle to reduce glare;
- Stainless Steel contacts guarantee long working life;easy sterilization ETO or steam;
- | full range of Mc Intosh & Miller blades, latex free;
- set is packaged in an anti-shock case and in cardboard box
- 2.5 V handles work with both alkaline batteries or rechargeable batteries
- good light transmission of 3,500 LUX with 2.5 V handles
- autoclavable up to 134° for 5 minutes (approx. 2,000 times) non magnetic blades.

### Cat. No. 34303

### Laryngoscope tips

Cat. No.	Name	Dimensions
34318	For cats and small dog breeds	90 mm x 14 mm
34319	For medium-sized dog breeds	100 mm x 14 mm
34320	For large dog breeds	130 mm x 18 mm
34321	For giant dog breeds and other large animals	156 mm x 18 mm

### Miller Laryngoscope

- vacuum bulb 2.5 V;
- | ergonomic light;
- easy sterilization without removing the bulb;
- | stainless steel guarantees long life;
- 3 Miller blades (1,2,3)
- | sold in a suitcase

Nr kat. 34306

### Miller Laryngoscope Bladesw

Cat. No.	Name	Dimensions
34310	For very small cats and dogs	65 mm x 12,5 mm
34311	For cats and small breed dags	78 mm x 11,5 mm
34312	For medium breed dogs	102 mm x 11,5 mm
34313	For big breed dogs	155 mm x 16,5 mm
34313	For giant breed dogs	195 mm x 16,5 mm





### Blood Pressure Measurement Device PetMAP graphic II

- specially developed for veterinary purposes;
- PPO technology that allows you to measure systolic and diastolic blood pressure, mean arterial pressure (MAP) and heart rate;
- the set includes 7 cuffs in various sizes (2 cm; 2.5 cm; 3 cm; 3.5 cm; 4 cm; 4.5 cm; 5.5 cm)
- | easy to use, large 4.3 "touch screen;
- | built-in memory;
- blood pressure range: systolic 40-265 mmHg, diastolic
   25-195 mmHg, MAP 30-200 mmHg, heart rate 40-240
   beats per minute;
- | weight: 400 g, size 15 cm x 8 cm x 3 cm
- | powered by 4 AA batteries.


Nr kat. 290083

### Blood Pressure Measurement Device PetMAP g3

- | specially developed for veterinary purposes;
- PPO technology that allows you to measure systolic and diastolic blood pressure, mean arterial pressure (MAP) and heart rate;
- | temperature measurement;
- the set includes 7 cuffs in various sizes (2 cm; 2.5 cm; 3 cm; 3.5 cm; 4 cm; 4.5 cm; 5.5 cm)
- easy to use, large 5 "touch screen (12,7 cm);
- built-in memory;
- blood pressure range: systolic 40-265 mmHg, diastolic
   25-195 mmHg, MAP 30-200 mmHg, heart rate 40-240
   beats per minute;
- weight: 460 g, size 15 cm x 8 cm x 3 cm
- | powered by 4 AA batteries.

Nr kat. 290083\_g3

# Diagnostic







### BLOOD PRESSURE

### Blood Pressure Measurement Device PetMAP + II

- | specially developed for veterinary purposes;
- PPO technology that allows you to measure systolic and diastolic blood pressure, mean arterial pressure (MAP) and heart rate;
- | temperature measurement ( range 29-46 °C);
- | optional ECG and Capno modules
- the set includes 7 cuffs in various sizes (2 cm; 2.5 cm; 3 cm; 3.5 cm; 4 cm; 4.5 cm; 5.5 cm)
- easy to use, large 5 "touch screen (12,7 cm);
- built-in memory;
- blood pressure range: systolic 40-265 mmHg, diastolic
   25-195 mmHg, MAP 30-200 mmHg, heart rate 40-240
   beats per minute;
- | weight: 400 g, size 15 cm x 8 cm x 3 cm
- | powered by 4 AA batteries.

Nr kat. 290083\_+II



### Blood Pressure Measurement Device PetMAP + II Sp02

- | specially developed for veterinary purposes;
- PPO technology that allows you to measure systolic and diastolic blood pressure, mean arterial pressure (MAP) and heart rate;
- | Sp02 measurement
- | temperature measurement ( range 29-46 °C);
- optional ECG and Capno modules
- the set includes 7 cuffs in various sizes (2 cm; 2.5 cm; 3 cm; 3.5 cm; 4 cm; 4.5 cm; 5.5 cm)
- easy to use, large 5 "touch screen (12,7 cm);
- built-in memory;
- blood pressure range: systolic 40-265 mmHg, diastolic
   25-195 mmHg, MAP 30-200 mmHg, heart rate 40-240
   beats per minute;
- | weight: 400 g, size 15 cm x 8 cm x 3 cm
- | powered by 4 AA batteries.

Nr kat. 290083\_+IISp02





# Injection and infusion

# B. BRAUN syringe



IV stand

MORE INFORMATION ON PAGE 30

# Injection and infusion





MORE INFORMATION ON PAGE 33

# 02 Injection and infusion

### SYRINGES / IV CATHETERS

### POLnet disposable syringe 3 components

three-component syringe: corpus, rubber seal, piston;disposable;

| made of top quality materials;

| clear scale.

### Cat. No. Name

112341P	Polnet disposable syringe 3 components 2 ml (100 pcs / package )
112342P	Polnet disposable syringe 3 components 5 ml ( 100 pcs / package )
112343P	Polnet disposable syringe 3 components 10 ml ( 100 pcs / package )
112344P	Polnet disposable syringe 3 components 20 ml ( 50 pcs / package )
112350P	Polnet disposable syringe 3 components 2 ml withoth needle ( 100 pcs / package )

# B. BRAUN disposable syringe 2 components

two-component syringe: corpus, piston;
disposable;
made of top quality materials;
box contains 100 pcs.

Cat. No.	Name
112341M	2 ml syringe, 100 pcs package
112342M	5 ml syringe, 100 pcs package
112343M	10 ml syringe, 100 pcs package
112344M	20 ml syringe, 100 pcs package

### Catheter POLNET

kateter dożylny

- | trójkątne ostrze igły
- | opakowanie zawiera 50 sztuk

Cat. No.	Name
121702P	Catheter POLNET 2.1x45 mm, G-14 – orange
121705P	Catheter POLNET 1.3x45 mm , G-18 – green
121706P	Catheter POLNET 1.1x32 mm , G-20 – pink
121708P	Catheter POLNET 0.9x25 mm , G-22 – blue
121709P	Catheter POLNET 0.7x19mm , G-24 – yellow
121710P	Catheter POLNET 0.62x19 mm , G-26 – purple







### B.BRAUN peripheral venous catheter

### | IV catheter;

- | triangular needle point;
- | flared-tip cannula;
- | box contains 50 pcs.


Cat. No.	
121709-BR	B.BRAUN G24 / 0.7x19 PVC, yellow
121708-BR	B.BRAUN G22 / 0.9x25 PVC, blue
121706-BR	B.BRAUN G20 / 1.1x33 PVC, pink
121705-BR	B.BRAUN G18 / 1.3x33 PVC, green
121704-BR	B.BRAUN G17 / 1.5x45 PVC, white
121703-BR	B.BRAUN G16 / 1.7x50 PVC, gray
121702-BR	B. BRAUN G14 / 2.2x50 PVC, orange

### POLnet disposable needle

- | ultra-sharp edge guarantees minimum resistance of skin and muscles during injection;
- silicon coating of the needle minimises resistance of tissues during injection;
- | ultra-thin needle wall facilitates the administration of highly viscupus fluids;
- | box contains 100 pcs.

Cat. No.	
121112P	Disp. needle 0,5 x 25 mm 100 pcs/pack POLnet
121141P	Disp. needle 0.9 x 13 mm 100 pcs/pack POLnet
121250P	Disp. needle 0.45 x 16 mm 100 pcs/pack POLnet
121253P	Disp. needle 0,6 x 25 mm 100 pcs/pack POLnet
121254P	Disp. needle 0,7 x 30 mm 100 pcs/pack POLnet
121255P	Disp. needle 0,7 x 40 mm 100 pcs/pack POLnet
121256P	Disp. needle 0,8 x 25 mm 100 pcs/pack POLnet
121257P	Disp. needle 0.8 x 22 mm 100 pcs/pack POLnet
121258P	Disp. needle 0,8 x 40 mm 100 pcs/pack POLnet
121142P	Disp. needle 0,9 x 25 mm 100 pcs/pack POLnet
121143P	Disp. needle 0,9 x 40 mm 100 pcs/pack POLnet
121260P	Disp. needle 1,1 x 25 mm 100 pcs/pack POLnet
121262P	Disp. needle 1,1 x 40 mm 100 pcs/pack POLnet
121263P	Disp. needle 1.2 x 25 mm 100 pcs/pack POLnet
121264P	Disp. needle 1,2 x 40 mm 100 pcs/pack POLnet
121266P	Disp. needle 1,6 x 25 mm 100 pcs/pack POLnet
121268P	Disp. needle 1,6 x 40 mm 100 pcs/pack POLnet
121269P	Disp. needle 1,8 x 40 mm 100 pcs/pack POLnet
121276P	Disp. needle 0.8 x 16 mm 100 pcs/pack POLnet
121350P	Disp. needle 2,1 x 40 mm 100 pcs/pack POLnet

# Injection and infusion







# 02 Injection and infusion

### NEEDLES / INFUSION EQUIPMENT / CONTAINERS

# 15 Margomed disposable infusion set

with filter and deaerator;length: 2.16 mm.

Cat. No. 230616

# IS Margomed disposable infusion set with connector

with filter and deaerator;
length: 2.16 mm;
steril;
best quality.

Nr kat. 230613

# IS Margomed disposable infusion set with connector

| with filter and deaerator;

- with needle;
- length: 2.16 mm;
- steril;
- best quality.

Nr kat. 230614

### Used needle container

- containers for used needles are certified by the National Institute of Hygiene;
- carried with an international safety warning label and instructions for use.

Cat. No.	Name	Capacity
25820	Used needle container	0,5 l
25821	Used needle container	11
25823	Used needle container	21
25829	Used needle container	3,5 l
25824	Used needle container	5
25825	Used needle container	10 I
25822	Used needle container	15 l
25827	Used needle container	30 I







### IV stand, 4 hooks

- stem and hook made of chromed stainless steel, the crosspiece
- | and wheels made of durable plastic;
- height: 160-245 cm;
- weight: 3.2 kg;
- 4 hooks.

Cat. No. 27842

### IV stand, 4 hooks

- | corpus made of stainless steel;
- | hook, crosspiece, and wheels made of durable plastic;
- height: 160–220 cm;
- weight: 2.9 kg;
- 4 hooks.

Cat. No. 27844

### IV stand, 4 hooks

- | whole: corpus, hook, crosspiece, and wheels made of stainless steel;
- | height: 160-220 cm;
- weight: 3.3 kg;
- 4 hooks.

Cat. No. 27839

# Injection and infusion

O2 STANDS





# O2 Injection and infusion

### STANDS / COMPRESSION BANDS

### IV stand, 2 hooks

corpus made of chromed steel;
hook and wheels made of durable plastic;
height: 160-245 cm;
weight: 2.9 kg;
2 hooks.

Z HOOKS.

Cat. No. **27841** 

### IV stand, 2 hooks

corpus made of aluminium;
hook and wheels made of durable plastic;
height: 158-248 cm;

| weight: 2.9 kg;

2 hooks.

Cat. No. **27843** 



### Compression banc

- | the band can be easily put on with one hand;
- | top quality materials used for production are durable and light;
- | wide band prevents pulling hair.

Cat. No.	Name	Fastener
25724	Jetpull baby compression band	double locking system; fast and slow release system
25722	Jetpull2 compression band	double locking system; fast and slow release system
25727	Fast compression band	single locking system; slow release system



# Laboratory

# Hematology analyser 3-diff

# Biochemistry analyser





MORE INFORMATION ON PAGE 45

# Laboratory

03



### MORE INFORMATION ON PAGE 46

Laboratory 03

### LABORATORY CONTAINERS

| made of stainless steel;

- | autoclavable at 121°C;
- available in various shapes, sizes, and capacities;
- | containers for instruments or medical waste,

with or without lid;

| bowls for fluids, jars with lid, thermometer stands, forceps stands.









Cat. No.	Name/Capacity
26600	Instrument container without lid /208 x 109 x h 15 mm
26601	Instrument container without lid /210 x 160 x h 25 mm
26599	Instrument container without lid /223 x 126 x h 45 mm
26602	Instrument container without lid /264 x 172 x h 47 mm
26603	Instrument container without lid /306 x196 x h 50 mm
26604	Instrument container without lid /355 x 254 x h 50 mm
26597	Instrument container without lid /380 x304 x h 50 mm
26598	Instrument container without lid /440 x 320 x h 64 mm
26605	Instrument container with lid / 223 x 126 x h 45 mm
26606	Instrument container with lid /264 x 172 x h 47 mm
26607	Instrument container with lid /306 x 196 x h 50 mm
26608	Instrument container with lid /355 x 254 x h 50 mm













26628

Cat. No.	Name/Dimensions/Ccapacity
26616	Fluids bowl Ø 56 x h 32 mm - 50 mL
26617	Fluids bowl Ø 88 x h 46 mm - 170 mL
26618	Fluids bowl Ø 128 x h 60 mm - 450 mL
26619	Fluids bowl Ø 158 x h 75 mm - 940 mL
26621	Fluids bowl Ø 208 x h 84 mm - 1940 mL
26622	Fluids bowl Ø 258 x h 92 mm - 3200 mL
26623	Jar with lid 0.5 i Ø 102 x 65 mm
26624	Jar with lid 1 i Ø 103 x 126 mm
26625	Jar with lid 2 i Ø 127 x 161 mm
26628	Thermometer stand $\phi$ 33 x h 80 mm
26626	Forceps stand $\phi$ 55x h 140 mm
26627	Forceps stand $\phi$ 55x h 180 mm

26627

# Laboratory

### LABORATORY CONTAINERS

03

- 39

### LITTER / URINE STRIPS / PLATES

### Cat litter

enables easy and stress-free sampling of animal urine for analysis..

Cat. No. 275021



### Urine analysis strips, 11 parameters

- strips for urinalysis used for qualitative and semi-qualitative urinalysis;
- sensitive to: leukocytes, nitrites, urobilinogen, protein, pH, blood, specific gravity, ascorbinic acid, ketone bodies, bilirubin,

and glucose;

- | for visual analysis and readers;
- | reading in about 1 minute (leukocytes in 60-120 seconds);
- | store in a dark place, at 10-30°C;
- | 1 box contains 100 strips.

Cat. No. 290798P

### Petri dishes

- culture Petri dishes;
- | glass, 90 mm in diameter;
- made of polystyrene (sterile and aseptic),
   diameter 90 mm, height 14.2 mm,
   without ventilation ribbing.

### Cat. No. Name/Capacity

290560	Glass Petri dishes, diameter 90 mm
610.092	Polystyrene Petri dishes, diameter 90 mm, sterile
611.092	Polystyrene Petri dishes, diameter 90 mm, aseptic



### Laboratory cover glasses

- high-quality microscope cover glasses, indispensible in every laboratory with a miscoscope;
- cover glasses available in 4 sizes;
- box contains 50 pieces (base glasses) and 100 pieces (cover glasses).

Cat. No.	
290284	Cover glasses, 20 x 20 (100 pcs)
290280	Cover glasses, 18 x 18 (100 pcs)
290282	Cover glasses, 22 x 22 (100 pcs)
290284A	Cover glasses, 24 x 24 (100 pcs)
290270P	Base glasses, unpolished, 76 x 26 mm (50 pcs)

### Urine container

- | urinalysis container with screw cap;
- with adhesive label;
- | made of polypropylene (25970) or polystyrene (25972).

Cat. No.	
25970	100 ml, urinalysis container, with screw cap

### Stool container

- | stool analysis container;
- | with screw cap and small spatula;
- | made of polypropylene.

Cat. No.	
25971	30 mL stool analysis container, with screw cap and small spatula

# Laboratory



### COVER GLASSES / CONTAINERS







03 Laboratory

### AGAR BASES

### MÜLLER HINTON agar base

- | used for determining microbial sensitivity to antibiotics and chemotherapeutics;
- indications: E. Coli, Streptococcus, Staphylococcus, Micrococcus;
- box contains 6 pcs.

Cat. No. 290614



### Direct action germicidal lamps

- high disinfection efficiency due to the emission of UV-C radiation with a wavelength of 253.7 nm;
- designed to raise and maintain the level of microbiological cleanliness of the rooms;
- direct germicidal lamp with different power, with one or two radiators;
- irreversibly deactivates bacteria, viruses, molds, fungi and all other microorganisms
- available in wall, ceiling and mobile version
- can be used with a remote control and also have a time counter.

Cat. No. NBV

### SABOURAUD agar base

used for culturing yeasts and fungi form material sampled mainly from skin infections;
indications: Candida;

box contains 6 pcs.

Cat. No. 290612



### MANNITOL agar base with salt

used for isolating and differentiating Staphylococci.
 Staphylococcus aureus (coagulase-positive)

- creates colonies with yellow rim;

- Staphylococcus epidermidis does not change the red colour of the base;
- growth of other microorganisms is minimal or arrested;box contains 6 pcs.

Cat. No. 290616



the growth of Trichophyton or Microsporum changes the base colour from yellow to red;box contains 6 pcs.

Cat. No. **290610** 





### Flow germicidal lamps

- high disinfection efficiency by means of UV-C rays operating inside the disinfection chamber;
- designed to raise and maintain the level of microbiological cleanliness of the rooms;
- the contaminated air is sucked into the disinfection chamber by the fan, where it is cleaned, the air blown out is practically free from germs;
- | possibility of use in the presence of man;
- germicidal lamp with different power, with one or two radiators;
- irreversibly deactivates bacteria, viruses, molds, fungi and all other microorganisms
- available in wall, ceiling and mobile version
- can be used with a remote control and also have a time counter.

Cat. No. NBVE

### Laboratory



UV LAMPS / SWABS





03 Laboratory

UV LAMPS / SWABS

### Flow two functional germicidal lamps

- high disinfection efficiency by means of UV-C rays operating inside the disinfection chamber (function I), as well as direct disinfection of the room (function II)
- designed to raise and maintain the level of microbiological cleanliness of the rooms;
- I the contaminated air is sucked into the disinfection chamber by the fan, where it is cleaned, the air blown out is practically free from germs;
- | possibility of use in the presence of man;
- bactericidal lamp with different power (60/30 or 110/55);
- | irreversibly deactivates bacteria, viruses, molds, fungi and all other microorganisms
- | available in wall, ceiling and mobile version
- can be used with a remote control and also have a time counter.

Nr kat. NBVE



### Swab in a tube

wooden stick;sterile.

Nr kat. wymazówka I

### Swab in a tube

swabs in transport tubes;

| plastic applicator;

sterile.

at. No. Name

581.2D Swab in transport tube

### Hematology Analyzer 3- diff Genrui VH-30

60 samples per hour;

- | 22 parameters and 3 histograms;
- | sample volume 10 µL;
- | 10.4" colour touch screen;
- | species: cat and dogs;
- | memory up to 600,000 measurements;
- | built-in thermal printer;
- advanced device diagnostics and self-cleaning methods;parameters: WBC, Neu#, Lym#, Gran#,
- Neu%, Lym%, Gran%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, MPV, PDW-SD, PDW-CV, PCT, P-LCR, P-LCC.

Nr kat. GEN\_VH30

# Laboratory







### WP21B VET chemistry analyzer

WP21B is a semi-automatic chemistry analyzer with large touch-screen and easy-to-use operation software, which is suitable for small labs, clinics.

### **Product Specifications:**

- 5-inch touch screen and keypad
- Permanent quartz flowcell
- | Long-life halogen lamp
- | Reaction curve real time display
- Support LIS
- | Flowcell & cuvette modes
- | Built-in thermal printer
- | Less and easy maintenance

### Parameters:

ALT, ALP, GGT, TP, Alb, TBIL, DBIL
LDH, CK, CK-MB, AST, A-HBDH
CHO, TG, HDL\_C, LDL\_C
Apo A1, Apo B, Lp(a)
BUN, CRE, UA
CO2-CP, AMS, LIP, GLU, IgA, IgG, IgM, C<sub>3</sub>,C<sub>4</sub>, CRP



### BA 600 VET

- 5-inch touch screen
- Strips: 11 parameters
- | 100 samples/hour
- Dual wavelength reflectance test method
- Auto & manual test mode
- | Thermal printer, auto & manual print
- | Less and easy maintenance

### Cat. No. GEN\_600B

### Electrolyte analyzer

- Animal species: Canine, Feline, Equine, Bovine, Porcine, Others
- | Throughput: 60 samples per hour
- Sample volume:100-200µL
- lon: K+, Na+, Cl-, iCa2+, pH, Li+, customized selection of ions, optional CO<sub>2</sub>
- | Calibration: Automatic or On-demand
- Closed reagent, real time monitoring reagent
- | Printer: Internal thermal printer
- Display: LCD
- Storage: Up to 50,000 results
- | Interface: RS-232, support LIS

Cat. No. GEN\_300G

### Coagulation analyzer

- | Parameters: PT, APTT, TT, FIB, etc.
- | Principle: Clotting, optical method
- | Disposable cuvette, no carry-over
- | Precise electronic pipette
- Automatic start at reagent addition
- Built-in thermal printer, auto&manual print


Cat. No. GEN\_520C

# Laboratory



ANALYZER







03 Laboratory

### CENTRIFUGE

### Centrifuge NF 200

- the NF 200 table centrifuge is designed for separating blood samples, sedimentation of urine particles, and for performing routine applications in microbiological laboratories;
- the use of a microprocessor control system allows the user to program and control the speed and time;
- | operation with 1.5/2 mL eppendorf tube rotors.

### Technical parameters:

- maximum speed: 5000 rpm;
- | maximum RCF: 2822 x g;
- maximum capacity: 12 x 1.5 mL;
- | microprocessor control system;
- | speed range: 1000–5000 rpm;
- | speed step: 10 rpm;
- | time range: 01-99 min + HOLD POSITION;
- time step: 1 min;
- induction motor;
- | power supply: 230 V/50 Hz;
- power: 350 W;
- external dimensions: 274 x 355 x 238;
- | net/gross weight (kg): 11/12.



Cat. No.	Name/Capacity
300005	NF 200 centrifuge
300006	1.5/2 mL tube, rotor reductor set
300007	13 x 100 mm tube, rotor reductor set
300008	13 x 75 mm tube, rotor reductor set

### Mikroscope Binocular Delta Optical Genetic Pro Bino

### | binocular microscope;

- two 10x wide field glasses;
- | 18 mm field of view;
- | achromatic lenses providing sharp image -
- 4x, 10x, 40x, 100x;
- LED lighting source;
- | pupil distance of 48-75 mm;
- | magnification of the microscope 40x, 100 x, 400 x, 1000 x;
- | 160 mm correction;
- | microscope table 142 mm x 132 mm;
- green filter, power cable, immersion oil, dust cover
- 24 months warranty.

Nr kat. DO-3402

### Mikroscope Trinocular Delta Optical Genetic Pro Bino

- | trinocular microscope;
- two 10x wide field glasses;
- 18 mm field of view;
- achromatic lenses providing sharp image
- 4x, 10x, 40x, 100x;
- LED lighting source;
- pupil distance of 48-75 mm;
- | magnification of the microscope 40x, 100 x, 400 x, 1000 x;
- | 160 mm correction;
- | microscope table 142 mm x 132 mm;
- green filter, power cable, immersion oil, dust cover
- 24 months warranty.

Nr kat. DO-3406

# Laboratory



### MICROSCOPES





03 Laboratory

### MICROSCOPES

### Mikroscope Delta Optical ProteOne

| biological microscope;

- | semi-planachromat class optics
- trinocular head, blocked in a way that prevents disassembly without the use of tools
- | large field of view of the glasses 20 mm
- | upper table lock
- solid, stable tripod and special storage for storing the power cord.
- has a mechanical table with vernier a special scale that increases the accuracy of reading the position of the preparation
- coaxial micro-macro knobs and a convenient knob-lever that allows, after obtaining a sharp image, to block the upper position of the table (thanks to this, when lowering the table to change the preparation, we easily place it back in the focus plane)
- 4x, 10x, 40x, 100x;
- LED lighting source;
- | pupil distance of 48-75 mm;
- magnification of the microscope 40x, 100 x, 400 x, 1000 x;infinity correction;
- green filter, power cable, immersion oil, dust cover24 months warranty.

Nr kat. DO-3800



### Temperature control table

- heating table intended for observing preparations and the behavior of objects at different temperatures;
- | useful among others in semen, parasites or bacteria;
- | the table is mounted on a standard microscope table;
- adjustable temperature range 20-50 ° C,
- display accuracy 0.1  $^\circ$  C.

Nr kat. DO-4811



# Hygiene

# Dry Gimette steriliser

# Europa B autoclave





MORE INFORMATION ON PAGE 58

# Hygiene



### MORE INFORMATION ON PAGE 59

### **Sweatshirts**

| professional;

- stylish;
- | 100% cotton;
- all female and male sizes;

| white, green, light blue, patterned.

Size	White colour	Green colour	Light blue colour	Pattern
XS	26124	26126	26128	-
s	26130	26140	26150	26170
М	26131	26141	26151	26171
L	26132	26142	26152	26172
XL	26133	26143	26153	26173
XXL	26134	26144	26154	-



### Trousers

- | professional;
- stylish;
- | 100% cotton;
- | all female and male sizes;
- | white, green, light blue.

Size	White colour	Green colour	Light blue colour
XS	26125	26127	26129
S	26135	26145	26155
М	26136	26146	26156
L	26137	26147	26157
XL	26138	26148	26158
XXL	26139	26149	26159



### Ultrasonic cleaner, mechanical timer

- | capacity: from 2.8 litre;
- | fast and safer than other methods;
- cleans oil and wax;
- | removes all thoroughly;
- | cleans the most stubborn dirt;
- | mechanical timer up to 60 min. with optional extension;
- | weight: from 4.1 kg.

# Ultrasound cleaner, mechanical time + overheating sensor

- | capacity: from 2.8 litre;
- | fast and safer than other methods;
- cleans oil and wax;
- | removes all thoroughly;
- | cleans the most stubborn dirt;
- mechanical timer, heating option with overheating sensor (in case of lack of fluids);
- | weight: from 4.1 kg.

### Digital ultrasound cleaner

- capacity 2.8 litre;
- | fast and safer than other methods;
- | cleans oil and wax;
- | removes all thoroughly;
- | cleans the most stubborn dirt;
- | digital timer, heating option with overheating sensor (in case of lack of fluids);
- | temperature monitor.

Ultrasound cleaners	Branson 2510 2,8 litre	Branson 3510 5,5 litre	Branson 5510 9,5 litre	Branson 8510 20 litre
Mechanical timer	35501	35510	35520	35531
Mechanical times + heating	35502	35511	35521	35532
Digital control, heating, timer	35503	35512	35522	35533
Chamber size [mm]:	240 x 140 x 100	290 x 150 x 150	290 x 240 x 150	495 x 280 x 150
Internal dimensions [mm]:	338 x 295 x 305	399 x 371 x 318	399 x 371 x 401	599 x 376 x 467
Capacity [litres]	2,8	5,5	9,5	20
Seal set + net	nie	tak	tak	tak
Weight [kg]:	4,1	4,9	5,9	11,8
Power consump. without heating [W]:	130	130	185	320
Power consumption with heating	239	335	469	881
Frequency [kHz]:	40	40	40	40

# Hygiene



ULTRASOUND CLEANERS











### ACCESSORIES FOR ULTRASOUND CLEANERS

















35537







Ultrasound cleaners	Branson 2510 2,8 litre	Branson 3510 5,5 litre	Branson 5510 9,5 litre	Branson 8510 20 litre
Branson 2510 cover	35507	35517	35527	35537
Perforated tray	35508	35518	35528	35538
Beaker positioning cover	35509	35519	35529	35539
Glass beaker	35540	35540	-	-
Glass beaker	35541	35541	35541	35541

| meets all requirements of everyday sterilisation;

- | fully automatic;
- durable;
- | guarantees complete safety of operations;
- | double walls of internal isolation;
- | interior made of 18/8 stainless steel;
- | timer up to 120 minutes.

Cat. No. 35600

- | meets all requirements of everyday sterilisation;
- | fully automatic;
- durable;
- | guarantees complete safety of operations;
- double walls of internal isolation;
- | interior made of 18/8 stainless steel;
- | timer up to 120 minutes;
- | set includes: fuses, neon lamps for time control, thermostat, power supply cord.

Cat. No. 35619

35517

# Hygiene



STERILISERS





### STERILISERS

### Dry Gimette steriliser

automatic;

durable;

shock-resistant;

| kills even particularly resistant microorganisms;

| made of stainless steel;

- | four perforated aluminium trays;
- | LED light indicating improper thermoregulation.

Cat. No. 35621



### Europa B autoclave

### class B;

- | microprocessor electronic system;
- | electronic control;
- automatic equalisation of atmospheric pressure;
- | multifunctional graphic display;
- triple door lock (door open only when inside pressure equals 0);
- | 6 automatic sterilisation cycles;
- | with anti-condensation function and automatic shutdown;
- | even heat distribution;
- stainless steel;
- the device always stays dry (thanks to a combination of vacuum and heat);
- clean and waste water containers provide better circulation and prolonged device life.

Cat. No. 35653





Cat. No.	Name	Capacity	Internal dimensions	External dimensions	Weight
35600	Tau 1,6	1,6 l	245 x 115 x 55 mm	310 x 180 x 180 mm	5 kg
35606	Tau 6	61	385 x 270 x 90 mm	495 x 320 x 220 mm	8 kg
35619	Tau 18	18 I	375 x 230 x 200 mm	440 x 310 x 380 mm	14 kg
35621	Dry Gimette 21	20	300 x 250 x 260 mm	440 x 340 x 400 mm	13 kg
35628	Dry Gimette 28	28 I	380 x 250 x 260 mm	550 x 340 x 400 mm	17 kg
32629	Dry Gimette 50	50 l	405 x 340 x 350 mm	605 x 460 x 480 mm	24,5 kg

### NUBYRA autoclave

- class N;
- electronic;
- 2 sterilisation cycles;
- | for sterilising unpacked tools;
- | controlled with a microprocessor;
- | LCD display with temperature and pressure visualisation;
- sterilisation cycle control;
  - easy locking system with a double safety sensor;
  - thermodynamic system;
  - | aluminium chamber;
- | 2 electronic temperature sensors (1 inside, 1 outside);
- | electronic converter for better pressure control;
- | 3 aluminum trays.

Cat. No. 35647

# Hygiene



AUTOCLAVES





### AUTOCLAVES

### Hydra Evo autoclave

class N;

- automatic;2 sterilisation cycles;
- | for sterilising unpacked instruments;
- locking system with a safety axis;
- thermodynamic system;
- preset cycles;
- one water container;
- 2 locks on the door;
- chamber made of stainless steel.
- thermostat;
- | 4 aluminum trays + 1 key for removing them.



Cat. No.	Name	Capacity	Internal dimensions	External dimensions	Weight
35653	Europa B	15 l	ø 245 x 320 mm	510 x 590 x 390 mm	54 kg
35657	Europa B	24	ø 245 x 500 mm	510 x 730 x 390 mm	58 kg
35647	Nubyra	61	ø 17 x 27 cm	450 x 450 x 305 mm	27 kg
35658	Hydra Evo	15 I	ø 24,5 x 32 cm	440 x 380 x 500 mm	35 kg
35659	Hydra Evo	18 I	ø 24,5 x 40 cm	440 x 380 x 500 mm	35 kg

### Autoklave Melag Vacuclav 40B+ (41B+)

Class B +;

- | chamber volume of 18 liters;
- chamber depth 35 cm;
- free-standing version or intended for connection to the water / sewage network;
- it allows to sterilize packed tools in 20 minutes in a short program B and unpackaged in 10 minutes - including pre-vacuum and drying;
- patented intelligent DRYtelligence drying reduces cycle time and helps save resources;
- specially developed algorithms automatically adapt the drying time to each load;
- the intuitive operating concept allows for easy and safe work using the large touch screen;
- software in Polish;
- programs can be selected quickly and intuitively and easily document the result of the cycle;
- integrated batch release function, including tool identification and traceability;
- | equipped with a water-cooled impeller pump

Nr kat. 30555 Vacuclav40B+; 30555 Vacuclav41B+



### Autoklave Melag Vacuclav 43B+(44B+)

Class B +;

- | chamber volume of 23 liters;
- chamber depth 45 cm;
- | free-standing version or intended for connection to the water / sewage network;
- it allows to sterilize packed tools in 20 minutes in a short program B and unpackaged in 10 minutes - including prevacuum and drying;
- patented intelligent DRYtelligence drying reduces cycle time and helps save resources;
- specially developed algorithms automatically adapt the drying time to each load;
- the intuitive operating concept allows for easy and safe work using the large touch screen;
- software in Polish;
- programs can be selected quickly and intuitively and easily document the result of the cycle;
- integrated batch release function, including tool identification and traceability;
- equipped with a water-cooled impeller pump.

### Nr kat. 30555\_Vacuclav43B+; 30555\_Vacuclav44B+

### Autoklave Melag Vacuclav 23B+ (24B+)

- Class B + Profi;
- | chamber volume of 22 liters;
- chamber depth 45 cm;
- free-standing version or intended for connection to the water / sewage network;
- quick sterilization processes in the cabinet: unpacked tools in the Quick Program together with the pre-vacuum phase can be removed after 15 minutes;
- sterilization of packaged tools depending on the size of the load, it only takes 28 minutes (plus drying);
- powerful vacuum pump and powerful steam generator ensure short cycle times, so you can sterilize wrapped hollow instruments in a quick "B" program in just 40 minutes (together with drying);
- integrated quality water treatment and electronic
   monitoring prevent operator errors and ensure high
   sterilization quality and maximum instrument protection;
- patented intelligent DRYtelligence drying reduces cycle time and helps save resources;
- specially developed algorithms automatically adapt the drying time to each load;
- programs can be selected quickly and intuitively and easily document the result of the cycle.

Nr kat. 30555\_Vacuclav23B+; 30555\_Vacuclav24B+

# Hygiene



WELDERS





time to each load; ument the result of the cycle.

### AUTOCLAVES

### Autoklave Melag Vacuclav 30B+(31B+)

Class B + Profi;

- | chamber volume of 17 liters;
- chamber depth 35 cm;
- | free-standing version or intended for connection to the water / sewage network;
- quick sterilization processes in the office: unpacked tools in the Quick Program together with the pre-vacuum phase can be removed after 15 minutes;
- sterilization of packaged tools depending on the size of the load, it only takes 28 minutes (plus drying).
- powerful vacuum pump and powerful steam generator ensure short cycle times, so you can sterilize wrapped hollow instruments in a quick "B" program in just 40 minutes (together with drying);
- Integrated quality water treatment and electronic monitoring prevent operator errors and ensure high sterilization quality and maximum instrument protection.
- patented intelligent DRYtelligence drying reduces cycle time and helps save resources;
- specially developed algorithms automatically adapt the drying time to each load;
- | programs can be selected quickly and intuitively and easily document the result of the cycle.

### Nr kat. 30555\_Vacuclav30B+; 30555\_Vacuclav31B+

### Autoklave Melag Vacuclav 24BL+

### Class B + Profi;

- | chamber volume of 29 liters;
- | chamber depth 60 cm;
- free-standing version or intended for connection to the water / sewage network;
- | quick sterilization processes in the office: unpacked tools in the Quick Program together with the pre-vacuum phase can be removed after 15 minutes;
- sterilization of packaged tools depending on the size of the load, it only takes 28 minutes (plus drying);
- | powerful vacuum pump and powerful steam generator ensure short cycle times, so you can sterilize wrapped hollow instruments in a quick "B" program in just 40 minutes (together with drying);
- Integrated quality water treatment and electronic monitoring prevent operator errors and ensure high sterilization quality and maximum instrument protection;
- patented intelligent DRYtelligence drying reduces cycle time and helps save resources;
- specially developed algorithms automatically adapt the drying time to each load;

| programs can be selected quickly and intuitively and easily document the result of the cycle.

Nr kat. 30555\_Vacuclav24BL+



.

UNION PRIMA

MELAG

### Autoklave Melag Euroklav 29VS+

### Class B;

- | chamber volume 18 l;
- | chamber depth 35 cm;
- | free-standing version;
- offers three "S" class programs for the sterilization of packaged, ordinary hollow instruments (eg scissors, needle holders, clamps, etc.) and porous loads, even when packaged;
- equipped with a powerful vacuum pump, the autoclave effectively removes air from the packaging and interior of the instruments by means of an initial vacuum;
- control and monitoring of parameters necessary for safe sterilization pressure and temperature is carried out through sensors and innovative microprocessor control;
- innovative integrated electronics allow for a comprehensive way of digital documentation - via CF card and MELAtrace, MELAprint42 printer or computer network.

Nr kat. 30555\_Euroklav29VS+

### Autoklave Melag Euroklav 23VS+

### Class B;

- chamber volume 22 l;
- chamber depth 45 cm;
- free-standing version;
- offers three "S" class programs for the sterilization of packaged, ordinary hollow instruments (eg scissors, needle holders, clamps, etc.) and porous loads, even when packaged;
- equipped with a powerful vacuum pump, the autoclave effectively removes air from the packaging and interior of the instruments by means of an initial vacuum;
- control and monitoring of parameters necessary for safe sterilization pressure and temperature is carried out through sensors and innovative microprocessor control;
- innovative integrated electronics allow for a comprehensive way of digital documentation - via CF card and MELAtrace, MELAprint42 printer or computer network.

Nr kat. 30555 Euroklav23VS+

# Hygiene



**AUTOCLAVES** 





### AUTOCLAVES

### Autoklave Melag Melatronic 15 EN+

### Class S;

- | chamber volume 7 l, free-standing;
- compact steam sterilizer with fractional flow and active drying ensures safe sterilization and drying results;
- works with the pulsating-gravity emptying method, meets the requirements of EN 13060;
- equipped with an active drying system;
- | microprocessor control automatically monitors program progress and ensures safe sterilization;
- the liquid crystal display shows the progress of sterilization and displays messages and commands, ensures convenient, easy and safe operation;
- high-quality stainless steel chamber, protective functions and solid parts ensure reliability for many years.

Nr kat. 30555\_Melatronic15EN+



### Class S;

- | 19 liter chamber volume;
- | works by pulsing and gravity emptying;
- equipped with temperature and pressure sensors and microprocessor control;
- | meets the requirements of EN 13060;
- | the whole cycle is automatic;
- | built-in water quality measurement, water measurement in the tank and a two-line alphanumeric display allow for easy and comfortable operation.

### Nr kat. 30555\_Melatronic23EN

### Self-welding autoclave foil

enforced corners to prevent the accumulation of dust;a layer of laminated transparent foil and a layer of moisture

- resistant paper;
- ensures high antibacterial protection;
- external welding of the entire length of the case;certified bags.

Cat. No.	Name	Size
35860	Self-welding autoclave foil	90 x 165 mm
35861	Self-welding autoclave foil	90 x 230 mm
35862	Self-welding autoclave foil	140 x 250 mm
35863	Self-welding autoclave foil	190 x 330 mm
35864	Self-welding autoclave foil	190 x 400 mm
35866	Self-welding autoclave foil	300 x 450 mm





### D-301 welder without accessories

- compact, without accessories;
- | perfect in dentistry as well as medical and veterinary laboratories;
- started manually;
- | vertical design to minimise area;
- removable roll handle;
- | constant welding temperature;
- | constant sealing time.

Cat. No. 35928

### D-351 welder, set

- | with roll handle and cut-off tool;
- perfect in dentistry as well as medical and veterinary laboratories;
- started manually;
- | vertical design to minimise area;
- | removable roll handle;
- constant welding temperature;
- | constant sealing time.

Cat. No. 35931

### Roll of welding sleeve

- | for autoclave sterilisation;
- in transparent polyester-polyethylen foil combined with medical paper;
- vapour and gaseous indicators provide clear changes of colour during sterilisation;
- 200 m rolls.

Cat. No.		
35810	Roll of welding sleeve	200 m x 50 mm
35811	Roll of welding sleeve	200 m x 75 mm
35812	Roll of welding sleeve	200 m x 100 mm
35813	Roll of welding sleeve	200 m x 150 mm
35814	Roll of welding sleeve	200 m x 200 mm
35817	Roll of welding sleeve	200 m x 250 mm
35815	Roll of welding sleeve	200 m x 300 mm
35816	Roll of welding sleeve	200 m x 400 mm


# Hygiene



WELDERS







### Latex gloves

- gloves contain latex;
- disposable;
- | top quality of manufacture;
- elastic;
- tension-resistant;
- | lightly powdered;
- | 100 pcs in package;
- | convenient package.

Cat. No.	Name
260801P	Latex gloves "S" /100 szt.
260802P	Latex gloves "M" /100 szt.
260803P	Latex gloves "L" /100 szt.
260805P	Latex gloves "XL" /100szt
260809P	Latex gloves "XS" /100 szt.



### Nitrile gloves

- contain nitrile;
- | blue colour;
- disposable;
- | top quality of manufacture;
- | slightly elastic;
- tension-resistant;
- | lightly powdered;
- | 100 pcs in package;
- | convenient package.

Cat. No. Name		
		Nitrile gloves "S" /100 szt.
		Nitrile gloves "M" /100 szt.
26	60682-N	Nitrile gloves "L" /100 szt.
26	60683-N	Nitrile gloves "XL" /100 szt.
26	60684-N	Nitrile gloves "XS" /100 szt.

### Mint latex gloves

- | gloves contain latex;
- | with a slight mint scent;
- disposable;
- | top quality of manufacture;
- elastic;
- tension-resistant;
- | lightly powdered;
- | 100 pcs in package;
- | convenient package.

Cat. No.	Name
260801-M	Mint latex gloves "S"/100 szt.
260802-M	Mint latex gloves "M"/100 szt.
260803-M	Mint latex gloves "L"/100 szt.
260805-M	Mint latex gloves "XL"/100 szt.
260809-M	Mint latex gloves "XS"/100 szt.



### Powder-free vinyl gloves

- contain vinyl;
- disposable;
- | top quality of manufacture;
- | slightly elastic;
- tension-resistant;
- powder-free;
- | 100 pcs in package;
- | convenient package.

Cat. No.	
260797P	Powder-free vinyl gloves "XS"/100 szt.
260799P	Powder-free vinyl gloves "S"/100 szt.
260798P	Powder-free vinyl gloves "M"/100 szt.
260795P	Powder-free vinyl gloves "L"/100 szt.
260796P	Powdre-free vinyl gloves "XL"/100 szt.

# Hygiene



GLOVES





### Slightly powdered vinyl gloves

	contain vinyl;				
disposa					
top qua	ality of manufacture;				
slightly	elastic;				
tension	-resistant;				
lightly p	powdered;				
100 pcs	100 pcs in package;				
conven	ient package.				
Cat. No.	Name				
260791P	Slightly powdered vinyl gloves "XS"/100 szt.				
260792P Slightly powdered vinyl gloves "S"/100 szt.					
260793P	Slightly powdered vinyl gloves "M"/100 szt.				
260794P Slightly powdered vinyl gloves "L"/100 szt.					

260795P Slightly powdered vinyl gloves "XL"/100 szt.

### Powdered surgical glove

gloves contain latex;	lightly powdered;
disposable;	50 pairs in package;
sterile;	convenient package.
top quality of manufacture;	
elastic;	
tension-resistant;	

Cat. No.	Clove size
260971	Powdered surgical gloves 6,0
260972	Powdered surgical gloves 6,5
260973	Powdered surgical gloves 7,0
260974	Powderes surgical gloves 7,5
260975	Powdered surgical gloves 8,0
260976	Powdered surgical gloves 8,5
260977	Powdered surgical gloves 9,0.

powder-free;

| 50 pairs in package;

convenient package.

### Powder-free surgical gloves

I	gloves	contain latex;

- disposable;
- sterile;
- | top quality of manufacture;
- elastic;
- tension-resistant;

Cat. No.	Clove size
260971PF	Powder-free surgical gloves 6,0
260972PF	Powder-free surgical gloves 6,5
260973PF	Powder-free surgical gloves 7,0
260974PF	Powder-free surgical gloves 7,5
260975PF	Powder-free surgical gloves 8,0







### Filter fleece coat

filter fleece coat;	colour: eg green;
thickness from 30 g/m <sup>2</sup> ;	finished with elastic cuffs.

Cat. No.	Coat size
260001L	L
260001XL	XL
260001XXL	XXL

### Foil coat

- white coat;
- | made from LDPE foil;
- thickness: 0.03 mm;
- | size: 80 x 160;
- | fixing: carton strip;
- | 100 pcs in package.

### Cat. No. 260007P

### Disposable cap, beret-style

| non-woven polypropylene cap, flat with elastic and on the edge;

- | thickness: 12 g/m<sup>2</sup>;
- size: 53 cm;
- | colour: white, blue and other;
- | box contains 100 pcs.

### Cat. No. 260003WET

### Disposable cap, clip-type

- | disposable folded clip cap;
- thickness: 14 g/m<sup>2</sup>;
- | polypropylene non-woven fabric;
- colour: white, blue and other;
- | box contains 100 pcs.

### Cat. No. 260004P

### Surgical mask

| non-woven fabric mask;

- 3 layers;
- 50 pcs in package;
- | elastic put on ears.

Cat. No. 260002P

# Hygiene



### COATS / MASKS











### WASHING AND DISINFECTION PRODUCTS

### Preparat Sekusept Pulver

| efficient in the washing and disinfection of medical devices and surfaces.

Cat. No. 325001



- | liquid concentrate;
- | dedicated for the washing and disinfection of tools in medical and dental surgeries as well as in hair and beauty salons.

Cat. No. 325002



- | chlorine-based product in the form of fizzy tablets;
- | dedicated for washing and disinfecting surfaces;
- | can also be used to inactivate a dangerous material;
- | efficient disinfection of sanitary facilities and drains.

### Cat. No. **325003**



- | for washing and disinfecting surfaces, with bactericidal and fungicidal action;
- ensures fast and reliable disinfection with perfect surface tolerance;
- | based on glucoprotamin.

Cat. No. 325004

### Laudamonium

- | surface disinfection product;
- | highly tolerable by skin, so its can be used in antifungal prevention.

Cat. No. 325005

| for fast disinfection by spraying

| or wiping small and hard to reach surfaces (eg tables, small equipment, contaminated surfaces).

Cat. No. 325006

# Hygiene









Hygiene 04

#### WASHING AND DISINFECTION PRODUCTS

#### Sani-Cloth Active wipes

| disinfecting and washing non-alcohol wipes, ready to use; | wide spectrum of action (eg noroviruses) with perfect material tolerance.

Cat. No. 325007



#### Incides N wipes

| practical antibacterial wipes soaked in alcohol; | for disinfecting small surfaces.

Cat. No. 325008



| product for hygienic and surgical hand disinfection as well as for the disinfection of patient's body before surgery; | based on chlorhexidine, a bactericidal and fungicidal agent;

a combination of mild, skin-friendly, actively washing

| maintains the natural, slightly acidic skin reaction (pH 5.0);

| inactivates HBV and HIV.

| mild liquid handwash;

| does not contain soap;

pleasant scent.

Cat. No. 325010

synthetic substances (APG),

Cat. No. 325011

- | fast-action disinfecting and washing foam;
- | can be used on surfaces sensitive to alcohols, like plexiglass or acryl.

Cat. No. 325009



- | innovative alcohol-based gel hand disinfectant providing complete antiviral protection;
- | no risk of skin irritation;
- | contains 7 conditioning substances (aloe, panthenol, glycerine, etc), so it is delicate for the skin in providing complete protection;
- recommended for both hygienic and surgical hand disinfection.

Cat. No. 3250012

# Hygiene



#### WASHING AND DISINFECTION PRODUCTS







Hygiene 04

#### WASHING AND DISINFECTION PRODUCTS

#### Skinman Complete

- | disinfectant with substances protecting and lubricating hand skin (pH 5.5);
- | effective against bacteria (also tubercle bacilli), fungi, and viruses, eg HBV and HIV'
- | recommended for sensitive skin;
- tested in accordance with EN 1500 and EN 12054.

Cat. No. 3250013



#### Silonda, Silonda lipid

- | hand moisturiser;
- | regulates the level of skin moisturisation and lubrication; | absorbs perfectly;
- | Silonda Lipid is a skin lubricant with regenerative properties;
- | improves skin elasticity and moisturisation.

Cat. No. 325014



#### Preparat Skinsept Color

- | for disinfecting skin before surgeries, injections, punctures, taking blood samples, and vaccination;
- | complex product including active substances of proven clinical efficacy;
- | disinfects and ungreases the skin, but also visualises the treated area by colouration;
- autosterile.

Cat. No. 325016

#### Dermados

| elbow dispenser for the hygienic dosage of washing, disinfection, and conditioning products, eg Manisoft, Spitaderm, Silonda; | suitable for 0.5 L bottles.

Cat. No. 325017

- | for disinfecting skin before surgeries, injections, punctures, taking blood samples, and vaccination;
- | complex product including active substances
- of proven clinical efficacy;
- effective against bacteria, fungi, and viruses found on skin;
- | disinfects and ungreases the skin;
- autosterile.

Cat. No. 3250015



| a tub for preparing solutions; capacity 2 L.

Cat. No. 325018

#### Hygiene



#### WASHING AND DISINFECTION PRODUCTS







04 Hygiene

#### WASHING AND DISINFECTION PRODUCTS

#### Closed bag dispense

Closed bag dispenser;capacity 800 mL.

Cat. No. 325019



#### Wall-mounted clip

| for 500 ml bottles with washing or conditioning products.

Cat. No. **325020** 




#### NOTE!

Buy a washing and disinfection pack and get a **FREE dispenser**. We adjust washing and disinfection programs to individual needs of Clients and we print out wall-mounted plans of washing and disinfection **FREE OF CHARGE**. For more information call **61 657 67 21**, or go to www.polnet-wet.pl and Facebook.

# Dressing materials

# FUNNYflex end POLflex adhesive bands



# Rozmiar: szerokość 10 cm / długość: 4,5 m Polnoi

#### Adhesive tape

| indicated for sensitive skin;

- | with acrylic glue;
- | highly permeable by air and vapour, so the skin can breathe and does not sweat.

Cat. No.	
161112	Advesive tape 1,25cm x5m
161113	Adhesive tape 2,5cm x 5m
161114	Adhesive tape 5cm x 5m
162001	Adhesive tape 1,25cm x 9,14m
162002	Adhesive tape 2,5cm x 9,14m
162003	Adhesive tape 5cm x 9,14m
161137	Adhesive tape 1,25cm x 9,1m
161138	Adhesive tape 2,5cm x 9,1m
161139	Adhesive tape 5cm x 9,1m
161147	Adhesive tape 1,25cm x 5m
161148	Adhesive tape 2,5cm x 5m
161149	Adhesive tape 5cm x 5m

#### Viscose knit band

| manufactured from quality viscose fibres giving the necessary softness, absorption, and permeability;

| bands knit in a way preventing fraying;

several widths for maximum convenience.

Cat. No.	
161140	5 cm x 4 m viscose knit band
161141	10 cm x 4 m viscose knit band
161142	15 cm x 4 m viscose knit band
161108	6 cm x 5 m viscose knit band
161109	8 cm x 5 m viscose knit band
161110	10 cm x 5 m viscose knit band
161118	14 cm x 5 m viscose knit band

- | for the treatment of venous and lymphatic system insufficiency;
- | adjusts perfectly to the shape of the body and does not slip when used;
- several widths available, with 5 m length.

Cat. No.	
161119	6 cm x 5 m elastic band
161120	8 cm x 5 m elastic band
161121	10 cm x 5 m elastic band
161143	12 cm x 5 m elastic band
161144	15 cm x 5 m elastic band

# Dressing materials

#### ADHESIVE TAPE/KNIT AND ELASTIC BANDS

05







#### ELASTIC BANDS

#### Adhesive elastic band POLflex

- exceptional adhesive properties of the band: individual layers unite, but the band does not stick to fur or hair;
- the stripe facilitates unrolling the band;
- durable material that can be easily put on and fastened but
- also easily removed without using scissors;
- | available in 4 popular colours, 3 widths 4.5 m long.

Cat. No.	Name
164000S	Adhesive elastic band POLflex 2.5 cm (mix)
164002S	Adhesive elastic band POLflex 5 cm (mix)
164004S	Adhesive elastic band POLflex 7,5 cm (mix)
164006S	Adhesive elastic band POLflex 10 cm (mix)
164008S	Adhesive elastic band POLflex 15 cm (mix)



#### FUNNY FLEX adhesive banc

- adhesive band based on a patented technology;
- | no scissors needed to tear off a fragment;
- | waterproof;
- controlled adhesion without compression;
- reusable.

Cat. No.	Name
164000P	FUNNY FLEX adhesive band 2,5 cm mix
164002P	FUNNY FLEX adhesive band 5 cm mix
164004P	FUNNY FLEX adhesive band 7,5 cm mix
164006P	FUNNY FLEX adhesive band 10 cm mix



#### Cotton-viscose wadding

- | made of a mixture of cleaned, degreased cotton and viscose fibers;
- white color;
- | fluffy, soft, delicate;
- highly absorbent;
- non-sterile.

Cat. No.	
161132	Cotton viscose wadding/200g
161133	Cotton viscose wadding/500g

#### Cellulose wadding - rol

for medical, hygienic and economic purposes;not sterile.

Cat. No.		
161136	Cellulose wadding - roll 150 g	

#### Cellulose Wadding

- high absorbency;
- | useful for all medical treatments and local skin
- disinfection;
- non-sterile.

Cat. No.	
161135	Cellulose Wadding 40 x 60 cm / 5 kg

#### Wad-gas tampon

- has an inner layer of wadded absorbent body and an outer layer of gauze;
- sterile.

#### Cat. No.

161101 Wad-gas tampon 9 cm x 9 cm

# Dressing materials



#### WADDING / ROLLS









#### WADDING / GAUZE COMPRESSES

#### Cotton-viscose wadding

| made from a mixture of cotton and viscose;

- | soft, selected for high absorption;
- in foil packaging;
- | folded wadding used in cosmetics and for medical purposes.

Cat. No.	Name
161131	Cotton-viscose wadding 100 g
161132	Cotton-viscose wadding 200 g
161133	Cotton-viscose wadding 500 g
161134	Cotton-viscose wadding 1000 g



#### Non-sterile gauze compress

sheet made of 5 layers;

| 100% cotton gauze, bleached without the use of chloride;

| 17 threads/cm<sup>2</sup> ;

- edges not folded;
- | without X-ray detectable thread.

Cat. No.	Name
160180	Non-sterile gauze compress/17N/8W 5 cm x 5 cm
160175	Non-sterile gauze compress/17N/8W 7,5 cm x 7,5 cm
160170	Non-sterile gauze compress/17N/8W 10 cm x 10 cm



#### Sterile gauze

| sterile gauze sold by the metre is made of 17-thread 100%

- cotton gauze, bleached without chloride (HP);
- sterilised with ethylene oxide (ETO);
- packed in blisters.

Cat. No.	
161128	Sterile gauze, 1/2 m² 17 N
161129	Sterile gauze, 1 m² 17 N
161130	Sterile gauze, 1/4 m² 17 N

#### Plaster bandage

- | bandage used for immobilising extremities after fractures and injuries;
- plaster bandage made of natural gypsum applied on both sides of a cotton base;
- the bandage must be soaked shortly, then initial application takes about 4-8 minutes;
- | manufactured in several widths, 3 m long;
- | individually packed.

Cat. No.	
161122	6 cm x 3 cm plaster bandage
161123	8 cm x 3 cm plaster bandage
161124	10 cm x 3 cm plaster bandage
161125	12 cm x 3 cm plaster bandage
161150	15 cm x 3 cm plaster bandage

# Dressing materials

#### GAUZES / TUPFERS / PLASTER BANDAGES

05





#### PLASTER LINING / BANDS

#### Synthetic plaster lining

| made of non-absorbent polyester, slightly elastic longitudinally and laterally;

- soft, delicate to the touch, 1.5-2 mm thick;
- | intended as a soft filling under a plaster bandage.

Cat. No.	Name
161115	5 cm z 3 m synthetic plaster lining
161116	8 cm z 3 m synthetic plaster lining
161117	10 cm z 3 m synthetic plaster lining
161118	15 cm z 3 m synthetic plaster lining



#### VET-LITE band

| durable, light immobilising band;

- submerged in hot water, it softens and can be formed at will, eg on a limb;
- does not react with urine and faecal matter.

Cat. No.	Name
160610	5 cm z 1.6 m VET-LITE band
160620	7.5 cm z 1.6 m VET-LITE band
160630	10 cm z 1.6 m VET-LITE band
160640	15 cm z 1.6 m VET-LITE band



#### Paw stabiliser

| plastic paw stabiliser;

- | with protective foam exceeding the edges of the stabiliser to protect the animal;
- | with multiple vents that enable air circulation;
- | made of very strong material, with no sharp edges;
- | can be cut to fit the length of paw.

Cat. No.	Name
307800	20 cm paw stabiliser, blue
307805	30 cm paw stabiliser, red
307810	40 cm paw stabiliser, green







# Surgicryl PGA



# PGA







MORE INFORMATION ON PAGE 112

#### Suture



MORE INFORMATION ON PAGE 88

#### PGA SUTURES / NYLON SUTURES

#### Surgicryl<sup>®</sup> PGA

#### Туре

- | Braided and coated multifilament
- Composition Polyglycolic acid
- | Coating Polycaprolactone and calcium stearate (1%)
- | Tissue reaction Minimal
- Absorption The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days.
- | Approximately
- 50% of tensile strength remain after 21 days.
- Presentation
- Box with 12 sutures

#### Cat. No. Name

153319	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 0
153320	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 1
153321	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 2
153322	Surgicryl® PGA 3/8 koła trójkątna 75 cm nr 1/0
153323	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 2/0
153324	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 3/0





#### PDO monofilament

- | Pollydioxanone Monofilament
- | braided and coated multifilament
- | absorbable
- box with 12 sutures.

Cat. No.	Name
153340P	PDO monofilament 0
153341P	PDO monofilament 1
153342P	PDO monofilament 2
153343P	PDO monofilament 1/0
153344P	PDO monofilament 2/0
153345P	PDO monofilament 3/0





#### Surgicryl<sup>®</sup> PGA

#### | Type

- | Braided and coated multifilament
- Composition Polyglycolic acid
- | Coating Polycaprolactone and calcium stearate (1%)
- | Tissue reaction Minimal
- Absorption The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days.
- | Approximately
- | 50% of tensile strength remain after 21 days.
- Presentation
- Box with 12 sutures

Cat. No.	Name
153319	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 0
153320	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 1
153321	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 2
153322	Surgicryl® PGA 3/8 koła trójkątna 75 cm nr 1/0
153323	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 2/0
153324	Surgicryl® PGA 3/8 koła, trójkątna, dł. 75 cm, nr 3/0

#### Daclon Nylon Suture

#### Features:

- | monofilament;
- | non- absorbable;
- | composition: poliamid
- no coating;
- colour: blue;
- approximately 10% a year by rapture of chemical links (hydrolytic action)
- box with 12 sutures;.

#### Advantages:

- excellent tissue support ;
- | monofilament with microscopic uniform diameter;
- | minimal tissue reaction;
- sythetic,
- excellent handling properties, outstanding strenght retention from extended wound support;
- high flexibility;
- | high tensile strength;
- excellent knot security;
- soft passage through the tissues without sawing, tissue drag and trauma;
- | guaranteed seal and product sterility..

#### Indications:

- | general soft tissue closing and/or ligation;
- | Daclon blue: skin closure, general surgery, plastic surgery, gastrointestinal surgery, gynaecology, obstetrics, orthopaedics;
- Daclon black : plastic surgery, ophthalmic surgery.





#### PDO SUTURES / NYLON SUTURES





#### In cassettes:

		Thread thickness	
Cat. No.	EP	USP	— Lenght
9015	1.5	4/O	100 m
9020	2	3/0	100 m
9030	3	2/0	100 m
9035	3.5	0	100 m
9040	4	1	100 m
9050	5	2	75 m
9060	6	3&4	80 m

#### Nici z igłą w blistrach:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
9100116	Round 1/2	16 mm	1	5/0- black	75 cm
9150116	Round 1/2	16 mm	1.5	4/0 - blue	75 cm
9070120	Round 1/2	20 mm	0.7	6/0 - black	75 cm
9100120	Round 1/2	20 mm	1	5/0 - black	75 cm
9150120	Round 1/2	20 mm	1.5	4/0 - blue	75 cm
9200120	Round 1/2	20 mm	2	3/0 - blue	75 cm
9300120	Round 1/2	20 mm	3	2/0 - blue	75 cm
9300122	Round 1/2 Round 1/2	22 mm 26 mm	<u> </u>	2/0- blue 4/0 - blue	75 cm
9150126 9200126	Round 1/2 Round 1/2	26 mm 26 mm	2	3/0 - blue	75 cm 75 cm
9300126	Round 1/2	26 mm	3	2/0 - blue	75 cm
9350126	Round 1/2	26 mm	3.5	0 - blue	75 cm
9400126	Round 1/2	26 mm	4	1 - blue	75 cm
9300130	Round 1/2	30 mm	3	2/0 - blue	75 cm
9350130	Round 1/2	30 mm	3.5	0 - blue	75 cm
9400130	Round 1/2	30 mm	4	1 - blue	75 cm
9300137	Round 1/2	37 mm	3	2/0 - blue	75 cm
9350137	Round 1/2	37 mm	3.5	0 - blue	75 cm
9400137	Round 1/2 Round 1/2	37 mm	<u>4</u> 3	1 - blue 2/0 - blue	75 cm
9300140 9350140	Round 1/2 Round 1/2	40 mm 40 mm	3.5	0 - blue	75 cm 75 cm
9400140	Round 1/2 Round 1/2	40 mm	4	1 - blue	75 cm
9500140	Round 1/2	40 mm	5	2 - blue	75 cm
9400150	Round 1/2	50 mm	4	1 - blue	100 cm
9350930	Round with loop 1/2	30 mm	3.5	0- blue	150 cm
9400930	Round with loop 1/2	30 mm	4	1 - blue	150 cm
9350935	Round with loop 1/2	35 mm	3.5	0 - blue	150 cm
9400935	Round with loop 1/2	35 mm	4	1 - blue	150 cm
9350940	Round with loop 1/2	40 mm	3.5	0 - blue	150 cm
9400940	Round with loop 1/2	40 mm	<u>4</u> 3.5	1 - blue 0 - blue	150 cm
9350950 9400950	Round with loop 1/2 Round with loop 1/2	50 mm 50 mm	<u> </u>	1 - blue	150 cm 150 cm
9300165-2	Double round 1/2	90 cm	3	2/0- blue	90 cm
9051512	Reverse cutting 3/8	12 mm	0.5	7/0 - black	75 cm
9071512	Reverse cutting 3/8	12 mm	0.7	6/0 - black	75 cm
9101512	Reverse cutting 3/8	12 mm	1	5/0 - black	75 cm
9151512	Reverse cutting 3/8	12 mm	1.5	4/0 - blue	75 cm
9071512P	Reverse cutting 3/8	12 mm Easy Pass	0.7	6/0- black	50 cm
9101512P	Reverse cutting 3/8	12 mm Easy Pass	1	5/0 - black	50 cm
9151512P	Reverse cutting 3/8	12 mm Easy Pass	1.5	4/0- blue	50 cm
9071516-45 9071516	Reverse cutting 3/8 Reverse cutting 3/8	<u>16 mm</u> 16 mm	0.7	6/0 -black 6/0 - black	45 cm 75 cm
9101516-45	Reverse cutting 3/8	16 mm	1	5/0 - black	45 cm
91011516	Reverse cutting 3/8	16 mm	1	5/0 - black	75 cm
9151516-45	Reverse cutting 3/8	16 mm	1.5	4/0 - blue	45 cm
9151516	Reverse cutting 3/8	16 mm	1.5	4/0 - blue	75 cm
9201516	Reverse cutting 3/8	16 mm	2	3/0- blue	75 cm
9071516P	Reverse cutting 3/8	16 mm - Easy Pass	0.7	6/0- black	50 cm
9101516P	Reverse cutting 3/8	16 mm- Easy Pass	1	5/0 - black	50 cm
9151516P	Reverse cutting 3/8	16 mm- Easy Pass	1.5	4/0- blue	50 cm
<u>9101519-45</u> 9101519	Reverse cutting 3/8 Reverse cutting 3/8	<u>19 mm</u> 19 mm	1	5/0- black 5/0 - black	45 cm 75 cm
9151519-45	Reverse cutting 3/8	19 mm	1.5	4/0 - blue	45 cm
9151519-45	Reverse cutting 3/8	19 mm	1.5	4/0 - blue	75 cm
9201519-45	Reverse cutting 3/8	19 mm	2	3/0- blue	75 cm
9201519	Reverse cutting 3/8	19 mm	2	3/0 - blue	75 cm
9301519	Reverse cutting 3/8	19 mm	3	2/0 - blue	75 cm
9151519P	Reverse cutting 3/8	19 mm Easy Pass	1.5	4/0 - blue	50 cm
9201519P	Reverse cutting 3/8	19 mm Easy Pass	2	3/0 - blue	50 cm
9301519P	Reverse cutting 3/8	19 mm Easy Pass	3	2/0 - blue	50 cm
9101524	Reverse cutting 3/8	24 mm	1	5/0- black	75 cm
9151524-45	Reverse cutting 3/8	24 mm	<u>1.5</u> 1.5	4/0 - blue 4/0 - blue	45 cm
9151524 9201524-45	Reverse cutting 3/8 Reverse cutting 3/8	24 mm 24 mm	2	3/0- blue	75 cm 45 cm
9201524-45	Reverse cutting 3/8	24 mm	2	3/0 - blue	75 cm
9301524-45	Reverse cutting 3/8	24 mm	3	2/0- blue	45 cm
9301524	Reverse cutting 3/8	24 mm	3	2/0- blue	75 cm
9351524	Reverse cutting 3/8	24 mm	3.5	0 - blue	75 cm
9401524	Reverse cutting 3/8	24 mm	4	1 - blue	75 cm
9501524	Reverse cutting 3/8	24 mm	5	2 - blue	90 cm
)					

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
930152P	Reverse cutting 3/8	24 mm Easy Pass	3	2/0 - blue	50 cm
9201526-45	Reverse cutting 3/8	26 mm	2	3/0 - blue	45 cm
9201526	Reverse cutting 3/8	26 mm	2	3/0 - blue	75 cm
9301526-45	Reverse cutting 3/8	26 mm	3	2/0 - blue	45 cm
9301526	Reverse cutting 3/8	26 mm	3	2/0 - blue	75 cm
9201530	Reverse cutting 3/8	30 mm	2	<u> 3/0 - blue</u>	75 cm
9301530	Reverse cutting 3/8	30 mm	<u> </u>	2/0 - blue	75 cm
9351530 9401530	Reverse cutting 3/8 Reverse cutting 3/8	30 mm 30 mm	4	0 - blue 1 - blue	<u>75 cm</u> 75 cm
9201536	Reverse cutting 3/8	36 mm	2	3/0 - blue	75 cm
9301536	Reverse cutting 3/8	36 mm	3	2/0 - blue	75 cm
9351536	Reverse cutting 3/8	36 mm	3.5	0 - blue	75 cm
9401536	Reverse cutting 3/8	36 mm	4	1 - blue	90 cm
9301540	Reverse cutting 3/8	40 mm	3	2/0 - blue	75 cm
9351540	Reverse cutting 3/8	40 mm	3.5	0 - blue	75 cm
9401540	Reverse cutting 3/8	40 mm	4	1 - blue	75 cm
9501540	Reverse cutting 3/8	40 mm	5	2 - blue	90 cm
9301545	Reverse cutting 3/8	45 mm	3	2/0 - blue	100 cm
9401545	Reverse cutting 3/8	45 mm	4	1 - blue	75 cm
9301560	Reverse cutting 3/8	60 mm	<u> </u>	2/0 - blue	75 cm
9351560 9301590	Reverse cutting 3/8 Reverse cutting 3/8	60 mm 90 mm	<u> </u>	0 - blue 2/0 - blue	75 cm 100 cm
9351590	Reverse cutting 3/8	90 mm	3.5	0 - blue	100 cm
9401590	Reverse cutting 3/8	90 mm	4	1 - blue	100 cm
9151417	Reverse cutting 1/2	17 mm	1.5	4/0- blue	75 cm
9201417	Reverse cutting 1/2	17 mm	2	3/0 - blue	75 cm
9201422	Reverse cutting 1/2	22 mm	2	3/0 - blue	75 cm
9201426	Reverse cutting 1/2	26 mm	2	3/0 - blue	75 cm
9301426	Reverse cutting 1/2	26 mm	3	2/0 - blue	75 cm
9351430	Reverse cutting 1/2	30 mm	3.5	0 - blue	75 cm
9351436	Reverse cutting 1/2	36 mm	3.5	0 -blue	75 cm
9401436	Reverse cutting 1/2	36 mm	4	1 - blue	75 cm
9501436 9351440	Reverse cutting 1/2 Reverse cutting 1/2	<u>36 mm</u> 40 mm	<u>5</u> 3.5	2 -blue 0 - blue	<u>75 cm</u> 75 cm
9401440	Reverse cutting 1/2	40 mm	4	1 - blue	75 cm
9401448	Reverse cutting 1/2	48 mm	4	1 - blue	100 cm
9202160	Straight cutting	60 mm	2	3/0 - blue	75 cm
9302160	Straight cutting	60 mm	3	2/0 - blue	75 cm
9201726	Cutting 3/8	26 mm	2	3/0 - blue	75 cm
9301730	Cutting 3/8	30 mm	3	2/0 - blue	75 cm
9301740	Cutting 3/8	40 mm	3	2/0 - blue	75 cm
MN102	Round micro 3/8	3.8 mm, 50µm	0.2	10/0 - black	10 cm
MN124	Round micro 3/8	5.1 mm, 70µm	0.01	12/0 - black	10 cm
MN114 MN104	Round micro 3/8 Round micro 3/8	<u>5.1 mm, 70µm</u> 5.1 mm, 70µm	0.1	<u>11/0 - black</u> 10/0 - black	<u>10 cm</u> 10 cm
MN104 MN113	Round micro 3/8	3.8 mm, 70µm	0.1	11/0- black	10 cm
MN105	Round micro 3/8	3.8 mm, 100µm	0.2	10/0- black	15 cm
MN106	Round micro 3/8	5.1 mm, 100µm	0.2	10/0 - black	15 cm
MN096	Round micro 3/8	5.1 mm, 100µm	0.3	9/0 - black	15 cm
MN087	Round micro 3/8	6.4 mm, 140µm	0.4	8/0 - black	15 cm
MN078	Round micro 3/8	9.7 mm, 200µm	0.5	7/0 - black	45 cm
ON81	Double round 1/4	2 x 8 mm, 220 µm	0.4	8/0- black	45 cm
ON102	Double armed spatula 1/2	2 x 6.2 mm, 150 µm	0.2	10/0 - black	30 cm
ON106	Double armed spatula 1/2	2 x 6.4 mm, 220 µm	0.2	10/0 - black	<u>30 cm</u>
ON62	Double armed spatula 1/2 Double armed spatula 3/8	2 x 8 mm, 300 µm 2 x 6 mm, 220 µm	0.7	<u>6/0 - black</u> 10/0 - black	<u>45 cm</u> 30 cm
ON103 ON111	Double armed spatula 3/8	2 x 6.2 mm, 220 µm	0.2	11/0- black	20 cm
- · · · · · · ·	2 Subic annea Spatala 5/0	2 x 6.2 mm, 150 µm	0.2	10/0 - black	30 cm
	Double armed spatula 3/8		0.3	9/0 - black	30 cm
ON101 ON95	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm	0.5	5/6 black	
ON101			0.3	9/0 - black	45 cm
ON101 ON95	Double armed spatula 3/8	2 x 6.2 mm, 150 µm			
ON101 ON95 ON95-45	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 μm 2 x 6.2 mm, 150 μm	0.3 0.4 0.4	9/0 - black	45 cm
ON101 ON95 ON95-45 ON85	Double armed spatula 3/8 Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 μm 2 x 6.2 mm, 150 μm 2 x 6.2 mm, 150 μm	0.3	9/0 - black 8/0 - black	45 cm 30 cm
ON101 ON95 ON95-45 ON85 ON85-45 ON75 ON105	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm 2 x 6.2 mm, 200 µm 2 x 6.4 mm, 220 µm	0.3 0.4 0.4 0.5 0.2	9/0 - black 8/0 - black 8/0 - black 7/0- black 10/0- black	45 cm 30 cm 45 cm 30 cm 30 cm
ON101 ON95 ON95-45 ON85 ON85-45 ON75 ON105 ON105-38	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm 2 x 6.2 mm, 200 µm 2 x 6.4 mm, 220 µm 2 x 6.4 mm, 220 µm	0.3 0.4 0.4 0.5 0.2 0.2	9/0 - black 8/0 - black 8/0 - black 7/0- black 10/0- black 10/0- black 10/0- black	45 cm 30 cm 45 cm 30 cm 30 cm 38 cm
ON101           ON95           ON95-45           ON85           ON85-45           ON75           ON105           ON105-38           ON105-45	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm 2 x 6.2 mm, 200 µm 2 x 6.4 mm, 220 µm 2 x 6.4 mm, 220 µm 2 x 6.4 mm, 220 µm	0.3 0.4 0.5 0.2 0.2 0.2 0.2	9/0 - black 8/0 - black 8/0 - black 7/0- black 10/0- black 10/0- black 10/0- black 10/0 - black	45 cm 30 cm 45 cm 30 cm 30 cm 38 cm 45 cm
ON101           ON95           ON95-45           ON85           ON85-45           ON75           ON105           ON105-38           ON105-45           ON61	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm 2 x 6.2 mm, 200 µm 2 x 6.4 mm, 220 µm	0.3 0.4 0.5 0.2 0.2 0.2 0.2 0.2 0.7	9/0 - black 8/0 - black 8/0 - black 7/0- black 10/0- black 10/0- black 10/0 - black 6/0- black	45 cm 30 cm 45 cm 30 cm 30 cm 38 cm 45 cm 45 cm
ON101           ON95           ON95-45           ON85           ON85-45           ON75           ON105           ON105-38           ON105-45	Double armed spatula 3/8 Double armed spatula 3/8	2 x 6.2 mm, 150 µm 2 x 6.2 mm, 200 µm 2 x 6.4 mm, 220 µm 2 x 6.4 mm, 220 µm 2 x 6.4 mm, 220 µm	0.3 0.4 0.5 0.2 0.2 0.2 0.2	9/0 - black 8/0 - black 8/0 - black 7/0- black 10/0- black 10/0- black 10/0- black 10/0 - black	45 cm 30 cm 45 cm 30 cm 30 cm 38 cm 45 cm

# Suture



NYLON SUTURES

#### Catgut Plain Suture

#### Features:

- twisted multifilament with a monofilament appearance;absorbable;
- colour: dark brown;
- strands of purified collagen taken from the serosal layer of selected bovines, origin from BSE-free classified countries;
- absorption by phagocytosis. Complete mass absorption is slowed down by chromizing - approximately 90 days.
   Approximately 50% of tensile strength remain after 14
- days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix,
- vagina) it is more quickly absorber; | box with 12 sutures.

#### Advantages:

- uniform absorption;
- excellent handling properties;
- high flexibility;
- high tensile strength;
- excellent knot security;
- oft passage through the tissues without sawing, tissue drag and trauma;
- excellent knot tie down;
- | guaranteed seal and product sterility

#### Indications:

- general soft tissue closing and/or ligation especially in general surgery;
- | gastrointestinal surgery;
- gynaecology;
- obstetrics;
- urology;
- ophthalmic surgery.

#### **Contraindications:**

- cardiovascular surgery
- neurosurgery being absorbable, not to use when an extended wound support is required;
- caution in patients with delayed wound healing (e.g. compromised immune system), in the region of infected wounds;
- known sensitivities or allergies to collagen.



#### Sutures without needle:

		ead thickness	I an alat
Cat. No.	EP	USP	Lenght
2200015	2	4/0	150 cm
2300015	3	3/0	150 cm
2350015	3.5	2/0	150 cm
2400015	4	0	150 cm
2500015	5	1	150 cm
2600015	6	2	150 cm
2700015	7	3	150 cm
2800015	8	4	150 cm
21000015	10	6	150 cm
2050	5	1	75 m
20505	5	1	50 m
2060	6	2	50 m
2070	7	3	50 m
2080	8	4	30 m
20802	8	4	25 m
2090	9	5	30 m

#### In cassettes:

Cat. No.         EP         USP           2015         1.5         5/0           2020         2         4/0	Lenght . 100 m 100 m
2020 2 4/0	100 m
2020 Z 4/0	
20205 2 4/0	100 m
<b>20305</b> 3 3/0	50 m
<b>2035</b> 3.5 2/0	100 m
<b>20355</b> 3.5 2/0	100 m
<b>2040</b> 4 0	50 m
<b>20407</b> 4 0	75 m
<b>20405</b> 4 0	50 m
<b>2050</b> 5 1	75 m
<b>20505</b> 5 1	50 m
<b>2060</b> 6 2	50 m
<b>2070</b> 7 3	50 m
2080 8 4	30 m
<b>20802</b> 8 4	25 m
<b>2090</b> 9 5	30 m
<b>2100</b> 10 6	30 m

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
2100115	Round ½	15 mm	1	6/0	75 cm
150115	Round ½	15 mm	1.5	5/0	75 cm
200115	Round ½	15 mm	2	4/0	75 cm
150120	Round ½	20 mm	1.5	5/0	75 cm
200120	Round ½	20 mm	2	4/0	75 cm
2300120	Round ½	20 mm	3	3/0	75 cm
2350120	Round ½	20 mm	3.5	2/0	75 cm
2400120	Round ½	20 mm	4	0	75 cm
2150126	Round ½	26 mm	1.5	5/0	75 cm
2200126	Round ½	26 mm	2	4/0	75 cm
2300126	Round ½	26 mm	3	3/0	75 cm
2350126	Round ½	26 mm	3.5	2/0	75 cm
2400126	Round ½	26 mm	4	0	75 cm
2500126	Round ½	26 mm	5	1	75 cm
2200130	Round ½	30 mm	2	4/0	75 cm
300130	Round ½	30 mm	3	3/0	75 cm
2350130	Round ½	30 mm	3.5	2/0	75 cm
400130	Round ½	30 mm	4	0	75 cm
2500130	Round ½	30 mm	4 4	1	75 cm
2200130	Round ½	37 mm	2	4/0	75 cm
2300137	Round ½	37 mm	3	3/0	75 cm
2350137	Round ½	37 mm	3.5	2/0	75 cm
	Round ½	37 mm	4	0	75 cm
2400137					
2500137	Round ½	37 mm	<u> </u>	1	75 cm
2350140	Round ½	40 cm		2/0	75 cm
400140	Round ½	40 cm	4	0	<u>75 cm</u>
2500140	Round ½	40 cm	5	1	75 cm
2500140-90		40 cm	5	1	<u>90 cm</u>
2600140	Round ½	40 cm	6	2	75 cm
2400145	Round ½	45 cm	4	0	75 cm
2600145	Round ½	45 cm	6	2	75 cm
2500148	Round ½	48 cm	5	1	75 cm
2600150	Round ½	50 cm	6	2	75 cm
2600150-90		50 cm	6	2	90 cm
2101512	Reverse cutting 3/8	12 mm	1	6/0	45 cm
2151512	Reverse cutting 3/8	12 mm	1.5	5/0	45 cm
2151516	Reverse cutting 3/8	16 mm	1.5	5/0	75 cm
201516	Reverse cutting 3/8	16 mm	2	4/0	75 cm
301516	Reverse cutting 3/8	16 mm	3	3/0	75 cm
2151519	Reverse cutting 3/8	19 mm	1.5	5/0	75 cm
2201519	Reverse cutting 3/8	19 mm	2	4/0	75 cm
2301519	Reverse cutting 3/8	19 mm	3	3/0	75 cm
2351519	Reverse cutting 3/8	19 mm	3.5	2/0	75 cm
2201524	Reverse cutting 3/8	24 mm	2	4/0	75 cm
301524	Reverse cutting 3/8	24 mm	3	3/0	75 cm
351524	Reverse cutting 3/8	24 mm	3.5	2/0	75 cm
401524	Reverse cutting 3/8	24 mm	2	0	75 cm
2501524	Reverse cutting 3/8	24 mm	3	1	45 cm
301530	Reverse cutting 3/8	30 mm	3.5	3/0	45 cm
351530	Reverse cutting 3/8	30 mm	4	2/0	75 cm
401530	Reverse cutting 3/8	30 mm	5	0	75 cm
501530	Reverse cutting 3/8	30 mm	3	1	75 cm
351537	Reverse cutting 3/8	37 mm	3.5	2/0	75 cm
401537	Reverse cutting 3/8	37 mm	4	0	75 cm
2501537	Reverse cutting 3/8	37 mm	5	1	75 cm
2351540	Reverse cutting 3/8	40 mm	3.5	2/0	75 cm
2601540	Reverse cutting 3/8	40 mm	6	2	75 cm
2401437	Reverse cutting 1/2	37 mm	4	0	75 cm
	Reverse cutting I/Z	37 mm	5	1	75 cm

#### Suture



#### CATGUT Plain Suture

#### Features:

- | twisted multifilament with a monofilament appearance; absorbable;
- colour: ivory;
- | strands of purified collagen taken from the serosal layer of selected bovines, origin from BSE-free classified countries;
- | absorption by phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorber;
- box with 12 sutures

#### Advantages:

excellent tissue support; uniform absorption; excellent handling properties; high flexibility; high tensile strength; excellent knot security; | soft passage through the tissues without sawing, tissue drag and trauma; excellent knot tie down; | guaranteed seal and product sterility

#### Indications:

- | general soft tissue closing and/or ligation especially in general surgery;
- gastrointestinal surgery;
- gynaecology;
- obstetrics;
- urology;
- ophthalmic surgery.

#### Contraindications:

- cardiovascular surgery
- | neurosurgery being absorbable, not to use when an extended wound support is required;
- caution in patients with delayed wound healing (e.g.
- compromised immune system), in the region of infected wounds; In cassettes:
- known sensitivities or allergies to collagen.

#### Sutures without needle:

		Thread thickness	
Cat. No.	EP	USP	Lenght
1200015	2	4/O	150 cm
1300015	3	3/0	150 cm
1250015	3.5	2/0	150 cm
1400015	4	0	150 cm
1500015	5	1	150 cm
1600015	6	2	150 cm
1700015	7	3	150 cm



#### Nici z igłą w blistrach:

		Needle size	Metric	USP	
Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
1150116	Round ½	16 mm	1.5	5/0	75 cm
1200116	Round ½	16 mm	2	4/0	75 cm
1200120	Round ½	20 mm	2	4/0	75 cm
1300120	Round ½	20 mm	3	3/0	75 cm
1350120	Round ½	20 mm	3.5	2/0	75 cm
1200126	Round ½	26 mm	2	4/0	75 cm
1300126	Round ½	26 mm	3	3/0	75 cm
1350126	Round ½	26 mm	3.5	2/0	75 cm
1400126	Round ½	26 mm	4	0	75 cm
1300130	Round ½	30 mm	3	3/0	75 cm
1300037	Round ½	37 mm	3	3/0	75 cm
1350137	Round ½	37 mm	3.5	2/0	75 cm
1400137	Round ½	37 mm	4	0	75 cm
1500137	Round ½	37 mm	5	1	75 cm
1500140-90	Round ½	40 mm	5	1	90 cm
1400145	Round ½	45 mm	4	0	75 cm
1500145	Round ½	45 mm	5	1	75 cm
1600145	Round ½	45 mm	6	2	75 cm
1151512	Reverse cutting 3/8	12 mm	1.5	5/0	45 cm
1101516	Reverse cutting 3/8	16 mm	1	6/0	45 cm
1151516	Reverse cutting 3/8	16 mm	1.5	5/0	75 cm
1201516	Reverse cutting 3/8	16 mm	2	4/0	75 cm
1301516	Reverse cutting 3/8	16 mm	3	3/0	75 cm
1201519	Reverse cutting 3/8	19 mm	2	4/0	75 cm
1301519	Reverse cutting 3/8	19 mm	3	3/0	75 cm
1351519	Reverse cutting 3/8	19 mm	3.5	2/0	75 cm
1201524	Reverse cutting 3/8	24 mm	2	4/0	75 cm
1301524	Reverse cutting 3/8	24 mm	3	3/0	75 cm
1351524	Reverse cutting 3/8	24 mm	3.5	2/0	75 cm
1401524	Reverse cutting 3/8	24 mm	4	0	75 cm
1351537	Reverse cutting 3/8	37 mm	3.5	2/0	75 cm

		Thread thickness	lopabt
Cat. No.	EP	USP	Lenght
1020	2	4/0	100 m
1030	3	3/0	100 m
1035	3.5	2/0	100 m
1040	4	0	100 m
1050	5	1	75 m
1060	6	2	50 m
1070	7	3	50 m
1080	8	4	30 m
1090	9	5	30 m
1100	10	6	30 m

#### Suture



#### Steel Suture - Dermafil, Vitafil

#### Features:

| braided multifilament; non-absorbable; composition: poliester | teflon / silicone coating; | colour: green- Dermafil/ white- Vitafil; box with 12 sutures;

#### Advantages:

excellent tissue support; coated multifilament | minimal tissue reaction; sythetic, excellent handling properties, outstanding strenght retention from extended wound support; high flexibility; high tensile strength; excellent knot security; | soft passage through the tissues without sawing, tissue drag and trauma; guaranteed seal and product sterility.

#### Indications:

general soft tissue closing and/or ligation; | cardiovascular surgery; ophthalmic surgery; skin closure; | gastrointestinal surgery; obstetrics orthopaedics.

#### Sutures without needle:

Cat. No.		Thread thickness	
Cat. NO.	EP	USP	- Lenght
7150015	1.5	4/O	150 cm
7200015	2	3/0	150 cm
7300015	3	2/0	150 cm
7350015	3.5	0	150 cm
7400015	4	1	150 cm
7500015	5	2	150 cm
7600015	6	3&4	150 cm
7150045	1.5	4/O	6 x 45 cm
7200045	2	3/0	6 x 45 cm
7300045	3	2/0	6 x 45 cm
7350045	3.5	0	6 x 45 cm
7400045	4	1	6 x 45 cm
7500045	5	2	6 x 45 cm
7600045	6	3&4	6 x 45 cm
7150025	1.5	4/O	250 cm ( Ligapack)
7200025	2	3/0	250 cm ( Ligapack)
7300025	3	2/0	250 cm ( Ligapack)
7350025	3.5	0	250 cm ( Ligapack)
7400025	4	1	250 cm ( Ligapack)
7500025	5	2	250 cm ( Ligapack)
7600025	6	3&4	250 cm ( Ligapack)



#### In cassettes:

		Thread thickness	
Cat. No.	EP	USP	- Lenght
7015	1.5	4/O	100 m
7020	2	3/0	100 m
7025	2.5	2/0	100 m
7030	3	2/0	100 m
7035	3.5	0	100 m
7040	4	1	100 m
7050	5	2	75 m
7060	6	3&4	50 m
7070	7	5	40 m
7080	8	6	40 m
6030	3	2/0	100 m
6035	3.5	0	100 m
6040	4	1	100 m
6050	5	2	100 m
6060	6	3&4	80m
6070	7	5	60 m
6080	8	6	50 m
6090	9	7	40 m
6100	10	8	30 m

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
7100116	Round 1/2	16 mm	1	5/0	75 cm
7150116	Round 1/2	16 mm	1.5	4/0	75 cm
7200116	Round 1/2	16 mm	2	3/0	75 cm
7150120	Round 1/2	20 mm	1.5	4/0	75 cm
7200120	Round 1/2	20 mm	2	3/0	75 cm
7200126	Round 1/2	26 mm	2	3/0	75 cm
7300126	Round 1/2	26 mm	3	2/0	75 cm
7350126	Round 1/2	26 mm	3.5	0	75 cm
7400126	Round 1/2	26 mm	4	1	75 cm
7200130	Round 1/2	30 mm	2	3/0	75 cm
7300130	Round 1/2	30 mm	3	2/0	75 cm
7350130	Round 1/2	30 mm	3.5	0	75 cm
7400130	Round 1/2	30 mm	4	1	75 cm
7300136	Round 1/2	36 mm	3	2/0	75 cm
7350136	Round 1/2	36 mm	3.5	0	75 cm
7400136	Round 1/2	36 mm	4	1	75 cm
7350140	Round 1/2	40 mm	3.5	0	75 cm
7400140	Round 1/2	40 mm	4	1	75 cm
7500140	Round 1/2	40 mm	5	2	75 cm
7700140	Round 1/2	40 mm	7	5	75 cm
7400148	Round 1/2	48 mm	4	1	75 cm
7500148	Round 1/2	48 mm	5	2	75 cm
7300117-2	Round ½ double	2 x 17 mm	3	2/0	90 cm
7300120-2	Round ½ double	2 x 20 mm	3	2/0	90 cm
7300126-2	Round ½ double	2 x 26 mm	3	2/0	90 cm
7200116-36	Round 1/2 - podwójna	16 , 36 mm	2	3/0	60 cm
7152019	Straight roung bodied taper point 1/2	19 mm	1.5	4/0	45 cm
7202019	Straight roung bodied taper point 1/2	19 mm	2	3/0	50 cm
7150517-2	Straight roung bodied taper point ½- double	2 x 17 mm	1.5	4/0	90 cm
7200517-2	Straight roung bodied taper point ½- double	2 x 17 mm	2	3/0	90 cm
7300517-2	Straight roung bodied taper point ½- double	2 x 17 mm	3	2/0	90 cm
7300517-2-10	Straight roung bodied taper point ½- double	2 x 17 mm	3	2/0	90 cm
7200525-2	Straight roung bodied taper point ½- double	2 x 25 mm	2	3/0	90 cm
7300525-2	Straight roung bodied taper point ½- double	2 x 25 mm	3	2/0	90 cm
7300525-2-10	Straight roung bodied taper point ½- double	2 x 25 mm	3	2/0	90 cm
7150518	Round ½ cutting	18 mm	1.5	4/0	75 cm
7200518	Round ½ cutting	18 mm	2	3/0	75 cm
7500540	Round ½ cutting	40 mm	5	2	75 cm
			7	5	

96

#### Suture



#### Polypropylene Suture

#### Features:

| monofilament;
| non- absorbable;
| composition: a polimer of propylene
| no coating;
| colour: blue;
| box with 12 sutures;

#### Advantages:

| excellent tissue support ;
| monofilament with microscopic uniform diameter;
| minimal tissue reaction;
| sythetic,
| xcellent handling properties, outstanding strenght retention from extended wound support;
| high flexibility;
| high tensile strength;
| excellent knot security;
| soft passage through the tissues without sawing, tissue drag and trauma;
| guaranteed seal and product sterility.

#### Indications:

| general soft tissue closing and/or ligation;
| cardiovascular and neurosurgery surgery;
| ophthalmic surgery, microsurgery
| skin closure, plastic surgery;
| gastrointestinal surgery;
| bstetrics, orthopaedics, gynaecology and ostetrics.

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
5050110	Round 1/2	10 mm	0.5	7/0	45 cm
5070110	Round 1/2	10 mm	0.7	6/0	75 cm
5070115	Round 1/2	15 mm	0.7	6/0	75 cm
5100115	Round 1/2	15 mm	1	5/0	75 cm
5150115	Round 1/2	15 mm	1.5	4/0	75 cm
5100117	Round 1/2	17 mm	1	5/0	75 cm
5200117	Round 1/2	17 mm	2	3/0	75 cm
5100120	Round 1/2	20 mm	1	5/0	75 cm
5150120	Round 1/2	20 mm	1.5	4/0	75 cm
5200120	Round 1/2	20 mm	2	3/0	75 cm
5100125	Round 1/2	25 mm	1	5/0	75 cm
5150125	Round 1/2	25 mm	1.5	4/0	75 cm
5200125	Round 1/2	25 mm	2	3/0	75 cm
5300125	Round 1/2	25 mm	3	2/0	75 cm
5350125	Round 1/2	25 mm	3.5	0	75 cm
5400125	Round 1/2	25 mm	4	1	75 cm
5150216	Round 3/8	16 mm	1.5	4/0	75 cm
5040108-2	Round 1/2	2 x 8 mm	0.4	8/0	60 cm
5050108-2	Round 1/2	2 x 8 mm	0.5	7/0	60 cm
5070108-2	Round 1/2	2 x 8 mm	0.7	6/0	60 cm
5050110-2	Round 1/2	2 x 10 mm	0.5	7/0	60 cm
5070110-2	Round 1/2	2 x 10 mm	0.7	6/0	60 cm
5100110-2	Round 1/2	2 x 10mm	1	5/0	60 cm
5050113-2	Round 1/2	2 x 13 mm	0.5	7/0	60 cm
5070113-2	Round 1/2	2 x 13 mm	0.7	6/0	60 cm



Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
5100113-2	Round 1/2	2 x 13 mm	1	5/0	75 cm
5100116-2	Round 1/2	2 x 16 mm	1	5/0	70 0111
5150116-2	Round 1/2	2 x 16 mm	1.5	4/0	90 cm
5200116-2	Round 1/2	2 x 16 mm	2	3/0	90 cm
5100120-2	Round 1/2	2 x 20 mm	1	5/0	90 cm
5150120-2	Round 1/2	2 x 20 mm	1.5	4/0	90 cm
5200120-2	Round 1/2	2 x 20 mm	2	3/0	90 cm
5150125-2	Round 1/2	2 x 25 mm	1.5	4/0	90 cm
5200125-2	Round 1/2	2 x 25 mm	2	3/0	90 cm
5300125-1	Round 1/2	2 x 25 mm	2	2/0	90 cm
5200130-2 5300130-2	Round 1/2 Round 1/2	2 x 30 mm 2 x 30 mm	3	<u> </u>	<u> </u>
5300130-2	Round 1/2	2 x 35 mm	3	2/0	90 cm
5040208-2	Round double 3/8	2 x 8 mm	0.4	8/0	90 cm
5050208-2	Round double 3/8	2 x 8 mm	0.5	7/0	45 cm
5070210-2	Round double 3/8	2 x 10 mm	0.7	6/0	60 cm
5100213-2	Round double 3/8	2 x 13 mm	1	5/0	60 cm
507021302	Round double 3/8	2 x 13mm	0.7	6/0	60 cm
5100218-2	Round double 3/8	2 x 18 mm	1	5/0	60 cm
5150218-2	Round double 3/8	2 x 18 mm	1.5	4/0	90 cm
5150222-2	Round double 3/8	2 x 22mm	1.5	4/0	90 cm
5071510	Reverse cutting 3/8	10 mm	0.7	6/0	45 cm
5051512	Reverse cutting 3/8	12 mm	0.5	7/0	45 cm
5071512	Reverse cutting 3/8	12 mm	0.7	6/0	45 cm
5101512	Reverse cutting 3/8	12 mm	1	5/0	45 cm
5151512	Reverse cutting 3/8	12 mm	1.5	4/0	75 cm
5071512P	Reverse cutting 3/8	12 mm Easy Pass	0.7	<u>6/0</u>	50 cm
5101512P	Reverse cutting 3/8	12 mm Easy Pass	0.7	5/0	50 cm
5071516 5101516	Reverse cutting 3/8 Reverse cutting 3/8	<u>16 mm</u> 16 mm	1	<u> </u>	<u>45 cm</u> 45 cm
5151516	Reverse cutting 3/8	16 mm	1.5	4/0	75 cm
5201516	Reverse cutting 3/8	16 mm	2	3/0	75 cm
5071516P	Reverse cutting 3/8	16 mm Easy Pass	0.7	6/0	50 cm
5101516P	Reverse cutting 3/8	16 mm Easy Pass	1	5/0	50 cm
5151516P	Reverse cutting 3/8	16 mm Easy Pass	1.5	4/0	50 cm
5101516VB	Reverse cutting 3/8	16 mm Black	1	5/0	45 cm
5071519	Reverse cutting 3/8	19 mm	0.7	6/0	75 cm
5101519	Reverse cutting 3/8	19 mm	1	5/0	75 cm
5151519-45	Reverse cutting 3/8	19 mm	1.5	4/0	45 cm
5151519	Reverse cutting 3/8	19 mm	1.5	4/0	75 cm
5201519-45	Reverse cutting 3/8	19 mm	2	3/0	45 cm
5201519	Reverse cutting 3/8	19 mm	2	3/0	75 cm
5301519	Reverse cutting 3/8	19 mm	3	2/0	75 cm
5351519	Reverse cutting 3/8	19 mm	3.5	0	75 cm
5101519P	Reverse cutting 3/8	19 mm Easy Pass	1	5/0	50 cm
5151519P	Reverse cutting 3/8	19 mm Easy Pass	1.5	4/0	50 cm
5201519P	Reverse cutting 3/8	19 mm Easy Pass	2	3/0	50 cm
5101519VB 5151519VB	Reverse cutting 3/8 Reverse cutting 3/8	19 mm black 19 mm black	1.5	<u> </u>	<u>75 cm</u> 75 cm
515151976	Reverse cutting 3/8	24 mm	1.5	4/0	75 cm
5201524-45	Reverse cutting 3/8	24 mm	2	3/0	45 cm
5201524-45	Reverse cutting 3/8	24 mm	2	3/0	75 cm
5301524-45	Reverse cutting 3/8	24 mm	3	2/0	45 cm
5301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
5351524	Reverse cutting 3/8	24 mm	3.5	0	75 cm
5401524	Reverse cutting 3/8	24 mm	4	1	75 cm
5201524P	Reverse cutting 3/8	24 mm Easy Pass	2	3/0	50 cm
5301524P	Reverse cutting 3/8	24 mm Easy Pass	3	2/0	50 cm
5201526-45	Reverse cutting 3/8	26 mm	2	3/0	45 cm
5301526-45	Reverse cutting 3/8	26 mm	3	2/0	45 cm
5301526	Reverse cutting 3/8	26 mm	3	2/0	75 cm
5201530	Reverse cutting 3/8	30 mm	2	3/0	75 cm
5301530	Reverse cutting 3/8	30 mm	3	2/0	75 cm
5351530	Reverse cutting 3/8	30 mm	3.5	0	75 cm
5401530	Reverse cutting 3/8	<u>30 mm</u>	4	1	75 cm
5201536	Reverse cutting 3/8	36 mm	2	3/0	75 cm
5301536	Reverse cutting 3/8	36 mm	3	2/0	75 cm
5351536	Reverse cutting 3/8	36 mm	3.5	0	75 cm
5401536	Reverse cutting 3/8	36 mm	4	1	75 cm

#### Suture



06 Suture

#### Polypropylene Suture

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
5301540	Reverse cutting 3/8	40 mm	3	2/0	75 cm
5351540	Reverse cutting 3/8	40 mm	3.5	0	75 cm
5401540	Reverse cutting 3/8	40 mm	4	1	90 cm
5501540	Reverse cutting 3/8	40 mm	5	2	90 cm
5401545	Reverse cutting 3/8	45 mm	4	1	90 cm
5401590	Reverse cutting 3/8	90 mm	4	1	90 cm
5201422	Reverse cutting 1/2	22 mm	2	3/0	45 cm
5301242	Reverse cutting 1/2	24 mm	3	2/0	75 cm
5202160	Reverse cutting 1/2	60 mm	2	3/0	75 cm
5302160	Reverse cutting 1/2	60 mm	3	2/0	75 cm
5302175	Reverse cutting 1/2	75 mm	3	2/0	75 cm
5100516-2	Round bodied taper cutting point 1/2 double	2 x 16 mm	1	5/0	75 cm
5150522-2	Round bodied taper cutting point 1/2 double	2 x 22 mm	1.5	4/0	90 cm
5150525-2	Round bodied taper cutting point 1/2 double	2 x 25 mm	1.5	4/0	90 cm
5200525-2	Round bodied taper cutting point 1/2 double	2 x 25 mm	2	3/0	90 cm
5300525-2	Round bodied taper cutting point 1/2 double	2 x 25 mm	3	2/0	90 cm
5300530-2	Round bodied taper cutting point 1/2 double	2 x 30 mm	3	2/0	90 cm
5070510	Round bodied taper cutting point 1/2	10 mm	0.7	6/0	75 cm
5300525	Round bodied taper cutting point 1/2	25 mm	3	2/0	75 cm
MP106	Micro taper point 3/8	5,1 mm, 100 µm	0.2	10/0	15 cm
MP097	Micro taper point 3/8	6,4 mm, 150 µm	0.3	9/0	15 cm
OPP105	Side cutting spatula double 3/0	2 x 6.4 mm, 200 µm	0.2	10/0	30 cm
OPP95	Side cutting spatula double 3/0	2 x 6.4 mm, 200 µm	0.3	9/0	30 cm
OPP103	Spatula 1/4	16 mm	0.2	10/0	10 cm
OPP103L	Loop spatula 1/4	16 mm, 150µm	0.2	10/0	15 cm
OPP101	Double cutting spatula	2 x 16,15 mm, 150µm	0.2	10/0	20 cm
OPP91	Double cutting spatula	2 x 16,15 mm, 150µm	0.3	9/0	20 cm
OPP110	Straight loop spatula	16,15 mm, 150 µm	0.2	10/0	15 cm

#### Features:

- | braided multifilament; non- absorbable; composition: silk silicone coating; colour: black; | Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues. box with 12 sutures; Advantages: | excellent tissue support for up to 3 months; coated multifilament minimal tissue reaction;
- sythetic,

100-

- excellent handling properties, outstanding strenght retention from extended wound support;
- high flexibility;
- high tensile strength;
- excellent knot security;
- | oft passage through the tissues without sawing, tissue drag and trauma;
- guaranteed seal and product sterility.



#### Indications:

general soft tissue closing and/or ligation;
cardiovascular surgery;
ophthalmic surgery;
skin closure;
gastrointestinal surgery;
obstetrics
orthopaedics.

#### In cassettes:

		Thread thickness	
Cat. No.	EP	USP	Lenght
8015	1.5	4/0	100 m
8020	2	3/0	100 m
8030	3	2/0	100 m
8035	3.5	0	100 m
8040	4	1	100 m
8050	5	2	75 m
8060	6	3&4	80 m
8070	7	5	60 m
8080	8	6	50 m
8090	9	7	40 m
8100	10	8	30 m

#### Suture



#### **Contraindications:**

- | urinary tract;
- biliary tract;
- | known allergies or sensitivities to silk.

#### Sutures without needle:

Cat No		Thread thickness	Lenght
Cat. No.	EP	USP	
8200015	2	3/0	150 cm
8300015	3	2/0	150 cm
8350015	3.5	0	150 cm
8400015	4	1	150 cm
8500015	5	2	150 cm
8600015	6	3&4	150 cm
8800015	8	6	150 cm
8200044	2	3/0	6 x 45 cm
8300044	3	2/0	6 x 45 cm
8350044	3.5	0	6 x 45 cm
8400044	4	1	6 x 45 cm
8500044	5	2	6 x 45 cm
8150046	1.5	4/0	17 x 45 cm
8200046	2	3/0	17 x 45 cm
8300046	3	2/0	17 x 45 cm
8350046	3.5	0	17 x 45 cm
8400046	4	1	17 x 45 cm
8500046	5	2	17 x 45 cm
8150060	1.5	4/0	13 x 60 cm
8200060	2	3/0	13 x 60 cm
8300060	3	2/0	13 x 60 cm
8350060	3.5	0	13 x 60 cm
8400060	4	1	13 x 60 cm
8200075	2	3/0	12 x 75 cm
8300075	3	2/0	12 x 75 cm
8350075	3.5	0	12 x 75 cm
8400075	4	1	12 x 75 cm
8500075	5	2	12 x 75 cm
8200025	2	3/0	250 cm ( Ligapack)
8300025	3	2/0	250 cm (Ligapack)
8350025	3.5	0	250 cm ( Ligapack)
8400025	4	1	250 cm ( Ligapack)
8500025	5	2	250 cm ( Ligapack)
8600025	6	3&4	250 cm ( Ligapack)

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
8070112	Round 1/2	12 mm	0.7	6/0	45 cm
8100112	Round 1/2	12 mm	1	5/0	45 cm
8150112	Round 1/2	12 mm	1.5	4/0	75 cm
8070112	Round 1/2	16 mm	0.7	6/0	45 cm
8100116	Round 1/2	16 mm	1	5/0	75 cm
8150116	Round 1/2	16 mm	1.5	4/0	75 cm
8200116	Round 1/2	16 mm	2	3/0	75 cm
8100120	Round 1/2	20 mm	1	5/0	75 cm
8150120	Round 1/2	20 mm	1.5	4/0	75 cm
8200120 8300120	Round 1/2	20 mm	2 3	3/0 2/0	75 cm 75 cm
8350120	Round 1/2 Round 1/2	20 mm 20 mm	3.5	0	75 cm
8200130	Round 1/2	30 mm	2	3/0	75 cm
8300130	Round 1/2	30 mm	3	2/0	75 cm
8350130	Round 1/2	30 mm	3.5	0	75 cm
8400130	Round 1/2	30 mm	4	1	75 cm
8500130	Round 1/2	30 mm	5	2	75 cm
8200137	Round 1/2	37 mm	2	3/0	75 cm
8300137	Round 1/2	37 mm	3	2/0	75 cm
8350137	Round 1/2	37 mm	3.5	0	75 cm
8400137	Round 1/2	37 mm	4	1	75 cm
8500137	Round 1/2	37 mm	5	2	75 cm
8350140	Round 1/2	40 mm	3.5	0	75 cm
8400140	Round 1/2	40 mm	4	1	75 cm
8500140	Round 1/2	40 mm	5	2	75 cm
8350145	Round 1/2	45 mm	3.5	0	75 cm
8400145	Round 1/2	45 mm	4	1	75 cm
8500140	Round 1/2	45 mm	5	2	75 cm
8150216	Round 3/8	16 mm	1.5	4/0	75 cm
8200216	Round 3/8	16 mm	2	3/0	75 cm
8200225	Round 3/8	25 mm	2	3/0	75 cm
8300225	Round 3/8	25 mm	3	2/0	75 cm
8300230	Round 3/8	30 mm	3	2/0	75 cm
8100518	Round ½ cutting	18 mm	1	5/0	75 cm
8150518	Round ½ cutting	18 mm	1.5	4/0	75 cm
8200518	Round ½ cutting	18 mm	2	3/0	75 cm
8300518	Round ½ cutting	18 mm	3	2/0	75 cm
8200527	Round ½ cutting	27 mm	2	3/0	75 cm
8300527 8150418	Round ½ cutting Round 3/8 tnaca	27 mm 18 mm	<u>3</u> 1.5	<u>2/0</u> 4/0	75 cm 75 cm
8200418	Round 3/8 thaca	18 mm	2	3/0	75 cm
8302040	Round bodied taper point	40 mm	3	2/0	75 cm
8051512-2	Double reverse cutting 3/8	2 x 12 mm	0.5	7/0	45 cm
8071512-2	Double reverse cutting 3/8	2 x 12 mm	0.7	6/0	45 cm
8051512	Reverse cutting 3/8	12 mm	0.5	7/0	45 cm
8071512	Reverse cutting 3/8	12 mm	0.7	6/0	45 cm
8101512	Reverse cutting 3/8	12 mm	1	5/0	45 cm
8151512	Reverse cutting 3/8	12 mm	1.5	4/0	75 cm
8201512	Reverse cutting 3/8	12 mm	2	3/0	75 cm
8071512P	Reverse cutting 3/8	12 mm Easy Pass	0.7	6/0	50 cm
8071512P	Reverse cutting 3/8	14 mm	0.7	6/0	45 cm
8071516	Reverse eatering 5,6				
8101516	Reverse cutting 3/8	16 mm	0.7	6/0	75 cm
8151516		16 mm 16 mm	1	6/0 5/0	75 cm 75 cm
0201516	Reverse cutting 3/8				
8201516	Reverse cutting 3/8 Reverse cutting 3/8	16 mm	1 1.5 2	5/0 4/0 3/0	75 cm
8301516	Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm	1 1.5 2 3	5/0 4/0 3/0 2/0	75 cm 75 cm 75 cm 75 cm 75 cm
8301516 8101516P	Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm 16 mm Easy Pass	1 1.5 2 3 1	5/0 4/0 3/0 2/0 5/0	75 cm 75 cm 75 cm 75 cm 75 cm 50 cm
8301516 8101516P 8151516P	Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass	1 1.5 2 3 1 1.5	5/0 4/0 3/0 2/0 5/0 4/0	75 cm 75 cm 75 cm 75 cm 75 cm 50 cm 50 cm
8301516 8101516P 8151516P 8071519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass 19 mm	1 1.5 2 3 1 1.5 0.7	5/0 4/0 3/0 2/0 5/0 4/0 6/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm
8301516 8101516P 8151516P 8071519 8101519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass 19 mm 19 mm	1 1.5 2 3 1 1.5 0.7 1	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm
8301516 8101516P 8151516P 8071519 8101519 8151519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass 19 mm 19 mm 19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm 75 cm
8301516 8101516P 8151516P 8071519 8101519 8151519 8201519-45	Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass 19 mm 19 mm 19 mm 19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm 75 cm 45 cm
8301516 8101516P 8151516P 8071519 8101519 8151519 8201519-45 8201519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm 16 mm 16 mm 16 mm Easy Pass 16 mm Easy Pass 19 mm 19 mm 19 mm 19 mm 19 mm 19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm 45 cm 75 cm
8301516 8101516P 8151516P 8071519 8101519 8151519 8201519-45 8201519 8301519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm         16 mm         16 mm         16 mm Easy Pass         16 mm Easy Pass         19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2 3	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0 3/0 2/0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm 45 cm 75 cm 75 cm 75 cm
8301516 8101516P 8151516P 8071519 8101519 8151519 8201519-45 8201519 8301519 8351519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm         16 mm         16 mm         16 mm Easy Pass         16 mm Easy Pass         19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2 3 3.5	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0 3/0 2/0 0	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm 75 cm 75 cm 45 cm 75 cm 75 cm 75 cm 75 cm 75 cm 75 cm
8301516 8101516P 8071519 8101519 8151519 8201519-45 8201519 8301519 8351519 8401519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm         16 mm         16 mm         16 mm Easy Pass         16 mm Easy Pass         19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2 3 3.5 4	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0 3/0 2/0 0 1	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm
8301516 8101516P 8071519 8101519 8151519 8201519-45 8201519 8301519 8351519 8401519 8151519P	Reverse cutting 3/8 Reverse cutting 3/8	16 mm         16 mm         16 mm         16 mm Easy Pass         16 mm Easy Pass         19 mm         19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2 3 3.5 4 1.5	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0 3/0 2/0 0 1 4/0	75 cm         75 cm         75 cm         50 cm         50 cm         75 cm
8301516 8101516P 8151516P 8071519 8101519 8151519 8201519-455 8201519 8301519 8351519 8401519	Reverse cutting 3/8 Reverse cutting 3/8	16 mm         16 mm         16 mm         16 mm Easy Pass         16 mm Easy Pass         19 mm	1 1.5 2 3 1 1.5 0.7 1 1.5 2 2 2 3 3.5 4	5/0 4/0 3/0 2/0 5/0 4/0 6/0 5/0 4/0 3/0 3/0 3/0 2/0 0 1	75 cm 75 cm 75 cm 75 cm 50 cm 50 cm 75 cm

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
8151519VB	Reverse cutting 3/8	19 mm black	1.5	4/0	75 cm
8201519VB	Reverse cutting 3/8	19 mm black	2	3/0	75 cm
8101524	Reverse cutting 3/8	24 mm	1	5/0	75 cm
8151524	Reverse cutting 3/8	24 mm	1.5	4/0	75 cm
8201524	Reverse cutting 3/8	24 mm	2	3/0	75 cm
8301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
8351524	Reverse cutting 3/8	24 mm	3.5	0	75 cm
8401524	Reverse cutting 3/8	24 mm	4	1	75 cm
8201526	Reverse cutting 3/8	26 mm	2	3/0	75 cm
8301526	Reverse cutting 3/8	26 mm	3	2/0	75 cm
8151530	Reverse cutting 3/8	30 mm	1.5	4/0	75 cm
8201530	Reverse cutting 3/8	30 mm	2	3/0	75 cm
8301530	Reverse cutting 3/8	30 mm	3	2/0	75 cm
8351530	Reverse cutting 3/8	30 mm	3.5	0	75 cm
8401530	Reverse cutting 3/8	30 mm	4	1	75 cm
8501530	Reverse cutting 3/8	30 mm	5	2	75 cm
8201537	Reverse cutting 3/8	37 mm	2	3/0	75 cm
8301537	Reverse cutting 3/8	37 mm	3	2/0	75 cm
8351537	Reverse cutting 3/8	37 mm	3.5	0	75 cm
8401537	Reverse cutting 3/8	37 mm	4	1	75 cm
8501537	Reverse cutting 3/8	37 mm	5	2	75 cm
8601537	Reverse cutting 3/8	37 mm	6	3&4	75 cm
8301545	Reverse cutting 3/8	45 mm	3	2/0	75 cm
8351545	Reverse cutting 3/8	45 mm	3.5	0	75 cm
8401545	Reverse cutting 3/8	45 mm	4	1	75 cm
8501545	Reverse cutting 3/8	45 mm	5	2	75 cm
8401560	Reverse cutting 3/8	60 mm	4	1	75 cm
8351590	Reverse cutting 3/8	90 mm	3.5	0	100 cm
8071412	Reverse cutting 1/2	12 mm	0.7	6/0	75 cm
8101417	Reverse cutting 1/2	17 mm	1	5/0	75 cm
8151417	Reverse cutting 1/2	<u>17 mm</u>	1.5	4/0	75 cm
8201417	Reverse cutting 1/2	17 mm	2	3/0	75 cm
8151420	Reverse cutting 1/2	20 mm	1.5	4/0	75 cm
8201420	Reverse cutting 1/2	20 mm	2	3/0	75 cm
8101422	Reverse cutting 1/2	22 mm	1	5/0	75 cm
8151422	Reverse cutting 1/2	22 mm	1.5	4/0	75 cm
8201422	Reverse cutting 1/2	22 mm	2	3/0	75 cm
8301422	Reverse cutting 1/2	22 mm	3	2/0	75 cm
8101426	Reverse cutting 1/2	26 mm		5/0	75 cm
8151426 8201426	Reverse cutting 1/2 Reverse cutting 1/2	26 mm 26 mm	1.5 2	<u>4/0</u> 3/0	75 cm 75 cm
8301426	Reverse cutting 1/2	26 mm	3	2/0	75 cm
8351426	Reverse cutting 1/2	26 mm	<u> </u>	0	75 cm
8351436	Reverse cutting 1/2	36 mm	3.5	0	75 cm
8401436	Reverse cutting 1/2	36 mm	4	1	75 cm
8302140	Straight cutting	40 mm	3	2/0	75 cm
8202140	Straight cutting	60 mm	2	3/0	75 cm
8302160	Straight cutting	60 mm	3	2/0	75 cm
8352160	Straight cutting	60 mm	3.5	0	75 cm
8402160	Straight cutting	60 mm	4	1	75 cm
O\$81B	Double spatula 3/8	2 x 6,5 mm, 200µm	0.4	8/0	45 cm
OS71	Double spatula 3/8	2 x 6,5 mm, 240µm	0.5	7/0	45 cm
OS61	Double spatula 3/8	2 x 6,5 mm, 240µm	0.7	6/0	45 cm
OS52	Double spatula 3/6 Double spatula 1/4	2 x 8 mm, 350µm	1	5/0	30 cm
OS62	Double spatula 1/4	2 x 8 mm, 250 µm	0.7	6/0	45 cm
OS51	Double spatula 1/4	2 x 10 mm, 350 µm	1	5/0	45 cm
5551	Double spatula I/4	2 Λ 10 ΠΠΠ, JJU μΠΤ	1	5/0	

#### Suture



#### Steel Suture

#### Features:

monofilament; absorbable; composition: iron, nickel and chromium alloy no coating; colour: metal; box with 12 sutures;

#### Advantages:

excellent tissue support; | Monofilament with microscopic uniform diameter minimal tissue reaction; sythetic, | excellent handling properties, outstanding strenght retention from extended wound support; high flexibility; high tensile strength; excellent knot security; | soft passage through the tissues without sawing, tissue drag and trauma; guaranteed seal and product sterility.

#### Indications:

closure of sternum; orthopaedics ; abdominal wall closure and hernia repair

#### Contraindications:

| known sensitivities or allergies to stainless or metals such as chromium and nickel

#### In cassettes:

Cot No	Thread	d thickness	
Cat. No.	EP	USP	-Lenght
10010	1	5/0	50 m
10015	1.5	4/0	50 m
10020	2	3/0	50 m
10030	3	2/0	50 m
10035	3.5	0	50 m
10040	4	1	50 m

# SUTURES

#### Sutures with needle:

	Thread thickness					
Cat. No.	Needle shape	Needle size	Metric	USP	-Lenght	
ST45	Round ½ cutting	40 mm	6	3&4	4 x 45 cm	
ST55	Round ½ cutting	48 mm	7	5	4 x 45 cm	
ST65	Round ½ cutting	48 mm	8	6	4 x 45 cm	
ST54	Round ½ cutting	54 mm	7	5	4 x 45 cm	
ST56	Cutting 1/2	48 mm	7	5	4 x 45 cm	
ST63	Cutting 1/2	48 mm	8	6	4 x 45 cm	
ST72	Cutting 1/2	48 mm	9	7	4 x 45 cm	
ST01	Reverse cutting 1/2	25 mm	3.5	0	4 x 45 cm	
STII	Reverse cutting 1/2	40 mm	4	1	4 x 45 cm	
ST21	Reverse cutting 1/2	40 mm	5	2	4 x 45 cm	
ST41	Reverse cutting 1/2	48 mm	6	3&4	4 x 45 cm	
ST71	Reverse cutting 1/2	48 mm	9	7	4 x 45 cm	
ST53	Reverse cutting 1/2	50 mm	7	5	4 x 45 cm	
ST51	Reverse cutting 1/2	50 mm	7	5	4 x 45 cm	
ST61	Reverse cutting 1/2	50 mm	8	6	4 x 45 cm	
ST52	Reverse cutting 1/2	58 mm	7	5	4 x 45 cm	
ST62	Reverse cutting 1/2	58 mm	8	6	4 x 45 cm	

#### Supramid black Suture

#### Features:

twisted and coated multifilament; 5/0 i 6/0 - monofilament non-absorbable; | composition: polyamide silicone coating; colour: black; box with 12 sutures;

#### Advantages:

excellent tissue support; coated multifilament | minimal tissue reaction; sythetic, excellent handling properties, outstanding strenght retention from extended wound support; high flexibility; high tensile strength; excellent knot security; soft passage through the tissues without sawing, tissue drag and trauma; guaranteed seal and product sterility.

#### Indications:

general soft tissue closing and/or ligation; skin closure.

#### In cassettes:

	Thread thickness		- Lowelat
Cat. No.	EP	USP	-Lenght
4015	1.5	4/0	100 m
4020	2	3/0	100 m
4030	3	2/0	100 m
4035	3.5	0	100 m
4040	4	1	100 m
4050	5	2	75 m
4060	6	3&4	50 m
4070	7	5	40 m
4080	8	6	40 m
4090	9	7	30 m







#### Sutures with needle:

ISOB         Bound 1/2         16 mm         15         4/0         75 cm           420016         Bound 1/2         20 mm         2         3/0         45 cm           4200120         Bound 1/2         20 mm         2         3/0         45 cm           6300125         Bound 1/2         26 mm         3         2/0         75 cm           415018         Bound 1/2         26 mm         15         4/0         50 cm           415018         Round 1/2 cutting         18 mm         15         4/0         75 cm           420018         Round 3/8 cutting         18 mm         2         3/0         75 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           415156         Reverse cutting 3/8         16 mm <t< th=""><th>Cat. No.</th><th>Needle shape</th><th>Needle size</th><th>Metric</th><th>USP</th><th>Lenght</th></t<>	Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
1200120         Round 1/2         20 mm         2         3/0         45 cm           6300126         Round 1/2         26 mm         3         2/0         75 cm           6300126         Round 1/2 cutting         18 mm         15         4/0         50 cm           4150818         Round 1/2 cutting         18 mm         15         4/0         75 cm           4200518         Round 3/2 cutting         18 mm         2         3/0         75 cm           4200418         Round 3/2 cutting         18 mm         2         3/0         75 cm           4200418         Round 3/2 cutting         18 mm         2         3/0         75 cm           4071512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4071512         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4071516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4051516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4201516         Reverse cutting 3/8         16 mm         1.5         4/0         75 cm           4201516         Reverse cuttin	4150116	Round 1/2	16 mm	1.5	4/0	75 cm
4300126         Round 1/2         26 mm         3         2/0         75 cm           4152022         Stright round bodied taper point 22 mm         15         4/0         75 cm           41200318         Round 1/2 cutting         18 mm         15         4/0         75 cm           4120418         Round 3/2 cutting         18 mm         15         4/0         75 cm           4120418         Round 3/2 cutting         18 mm         15         4/0         75 cm           4200418         Round 3/2 cutting 3/8         12 mm         0.7         6/0         75 cm           4101512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4101512         Reverse cutting 3/8         12 mm         1         5/0         75 cm           4101516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4101516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         45 cm           4201516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4201516         R	4200116	Round 1/2	16 mm	2	3/0	75 cm
I452022       Straight round bodie taper point 22 mm       15       4/0       50 cm         I450518       Round ½ cutting       18 mm       1.5       4/0       75 cm         I420518       Round ½ cutting       18 mm       2       3/0       75 cm         I430418       Round ¾ cutting       18 mm       2       3/0       75 cm         I400418       Round ¾ cutting       18 mm       2       3/0       75 cm         I400418       Round ¾ cutting ¾       18 mm       2       3/0       75 cm         I401512-5       Reverse cutting ¾       12 mm       1       5/0       45 cm         I01512       Reverse cutting ¾       12 mm       1       5/0       45 cm         I01516       Reverse cutting ¾       16 mm       1       5/0       75 cm         I01516       Reverse cutting ¾       16 mm       1       5/0       75 cm         I01516       Reverse cutting ¾       16 mm       15       4/0       45 cm         I01516       Reverse cutting ¾       16 mm       15       4/0       75 cm         I01516       Reverse cutting ¾       16 mm       2       3/0       75 cm         I01516       Reverse cut	4200120	Round 1/2	20 mm	2	3/0	45 cm
ISDB         Round k cutting         IB mm         15         4/0         75 cm           4200518         Round 3/8 cutting         IB mm         2         3/0         75 cm           4200418         Round 3/8 cutting         IB mm         2         3/0         75 cm           4200418         Round 3/8 cutting         IB mm         2         3/0         75 cm           4200418         Reverse cutting 3/8         12 mm         0         75 cm         15           401512-45         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         75 cm           401516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           401516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           401516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse	4300126	Round 1/2	26 mm	3	2/0	75 cm
4200518         Pound ½ cutting         18 mm         2         3/0         75 cm           4150418         Round 3/8 cutting         18 mm         15         4/0         75 cm           4200418         Round 3/8 cutting         18 mm         2         3/0         75 cm           4071512         Reverse cutting 3/8         12 mm         0.7         6/0         75 cm           401512.5         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4101512.7         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4071516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4001516.45         Reverse cutting 3/8         16 mm         1         5/0         75 cm           401516         Reverse cutting 3/8         16 mm         15         4/0         45 cm           4151516         Reverse cutting 3/8         16 mm         15         4/0         75 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4201516	4152022	Straight round bodied taper point	: 22 mm	1.5	4/0	50 cm
HSONB         Round 3/8 cutting         18 mm         15         4/0         75 cm           4200418         Round 3/8 cutting         18 mm         2         3/0         75 cm           4200418         Round 3/8 cutting         18 mm         2         3/0         75 cm           4001512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4001512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401512         Reverse cutting 3/8         12 mm         1         5/0         45 cm           401516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           401516         Reverse cutting 3/8         16 mm         1         5/0         45 cm           415156         Reverse cutting 3/8         16 mm         15         4/0         45 cm           420156         Reverse cutting 3/8         16 mm         2         3/0         75 cm           420156         Reverse cutting 3/8         16 mm         2         3/0         75 cm           420156         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4201597         Rever	4150518	Round ½ cutting	18 mm	1.5	4/0	75 cm
420018         Round 3/8 cutting         18 mm         2         3/0         75 cm           4071512         Reverse cutting 3/8         12 mm         0.7         6/0         75 cm           4010152.45         Reverse cutting 3/8         12 mm         1         5/0         46 cm           4101512         Reverse cutting 3/8         12 mm         1         5/0         75 cm           4101516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4001516         Reverse cutting 3/8         16 mm         1         5/0         45 cm           401516         Reverse cutting 3/8         16 mm         15         4/0         45 cm           4151516         Reverse cutting 3/8         16 mm         15         4/0         75 cm           4201516.75         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516.75         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4201516.75         Reverse cutting 3/8         19 mm         0.7         6/0         75 cm           4001519         Reverse cutting 3/8         19 mm         1         5/0         75 cm           410	4200518	Round ½ cutting	18 mm	2	3/0	75 cm
4071912         Reverse cutting 3/8         12 mm         0.7         6/0         75 cm           4101912         Reverse cutting 3/8         12 mm         1         5/0         45 cm           4101912         Reverse cutting 3/8         12 mm         1         5/0         75 cm           4101912         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4101916         Reverse cutting 3/8         16 mm         1         5/0         45 cm           4101916         Reverse cutting 3/8         16 mm         1         5/0         45 cm           4101916         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4101916         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4201916         Reverse cutting 3/8         16 mm         2         3/0         45 cm           4201916         Reverse cutting 3/8         16 mm black         1         5/0         75 cm           4201916         Reverse cutting 3/8         19 mm         0.7         6/0         75 cm           4101919         Reverse cutting 3/8         19 mm         1.5         4/0         75 cm           42019	4150418	Round 3/8 cutting	18 mm	1.5	4/O	75 cm
4101512-45       Reverse cutting 3/8       12 mm       1       5/0       45 cm         4101512       Reverse cutting 3/8       12 mm       1       5/0       75 cm         4101516       Reverse cutting 3/8       16 mm       0.7       6/0       75 cm         4071516       Reverse cutting 3/8       16 mm       0.7       6/0       75 cm         401516-45       Reverse cutting 3/8       16 mm       1       5/0       45 cm         4101516-45       Reverse cutting 3/8       16 mm       1.5       4/0       75 cm         4151516       Reverse cutting 3/8       16 mm       1.5       4/0       75 cm         4201516       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4001519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         401519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         4151519       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4301519       Reverse cutting 3/8       19 mm black       1       5/0       75 cm	4200418	Round 3/8 cutting	18 mm	2	3/0	75 cm
4101512         Reverse cutting 3/8         12 mm         1         5/0         75 cm           4151512         Reverse cutting 3/8         12 mm         1.5         4/0         45 cm           4071516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4101516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4101516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4151516         Reverse cutting 3/8         16 mm         1.5         4/0         75 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         45 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4201516         Reverse cutting 3/8         19 mm         1.5         4/0         45 cm           4101519         Reverse cutting 3/8         19 mm         1.5         4/0         75 cm           4201519         Reverse cutting 3/8         19 mm back         1         5/0         75 cm           42	4071512	Reverse cutting 3/8	12 mm	0.7	6/0	75 cm
HSISI2         Reverse cutting 3/8         12mm         15         4/0         45 cm           4071516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4101516-45         Reverse cutting 3/8         16 mm         1         5/0         45 cm           4101516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4151516         Reverse cutting 3/8         16 mm         1.5         4/0         75 cm           4201516-75         Reverse cutting 3/8         16 mm         2         3/0         45 cm           4201516-75         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516-75         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516-75         Reverse cutting 3/8         19 mm         0.7         6/0         75 cm           4201519         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4071519         Reverse cutting 3/8         19 mm         1.5         4/0         75 cm           4301519         Reverse cutting 3/8         19 mm Easy Pass         1.5         4/0         75 cm <t< td=""><td>4101512-45</td><td>Reverse cutting 3/8</td><td>12 mm</td><td>1</td><td>5/0</td><td>45 cm</td></t<>	4101512-45	Reverse cutting 3/8	12 mm	1	5/0	45 cm
4071516         Reverse cutting 3/8         16 mm         0.7         6/0         75 cm           4101516-45         Reverse cutting 3/8         16 mm         1         5/0         45 cm           4101516         Reverse cutting 3/8         16 mm         1         5/0         75 cm           4151516         Reverse cutting 3/8         16 mm         1.5         4/0         45 cm           4151516         Reverse cutting 3/8         16 mm         2         3/0         45 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         16 mm         2         3/0         75 cm           4201516         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4101519         Reverse cutting 3/8         19 mm         1         5/0         75 cm           4101519         Reverse cutting 3/8         19 mm         1.5         4/0         45 cm           4151519         Reverse cutting 3/8         19 mm         3         2/0         75 cm           4201519         Reverse cutting 3/8         19 mm black         1         5/0         75 cm           4301	4101512	Reverse cutting 3/8	12 mm	1	5/0	75 cm
4101516-45       Reverse cutting 3/8       16 mm       1       5/0       45 cm         4101516       Reverse cutting 3/8       16 mm       1       5/0       75 cm         4101516       Reverse cutting 3/8       16 mm       1.5       4/0       45 cm         415156       Reverse cutting 3/8       16 mm       1.5       4/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       45 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4151519       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm back       1.5       4/0       75 cm         4301519       Reverse cutting 3/8       19 mm black       1.5       4/0       7	4151512	Reverse cutting 3/8	12mm	1.5	4/O	45 cm
H01516       Reverse cutting 3/8       16 mm       1       5/0       75 cm         H31516       Reverse cutting 3/8       16 mm       1.5       4/0       45 cm         H31516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm black       1       5/0       75 cm         401519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         401519       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4151919       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4151919       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4151919-75       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       <	4071516	Reverse cutting 3/8	16 mm	0.7	6/0	75 cm
H01516       Reverse cutting 3/8       16 mm       1       5/0       75 cm         H31516       Reverse cutting 3/8       16 mm       1.5       4/0       45 cm         H31516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4201516-75       Reverse cutting 3/8       16 mm black       1       5/0       75 cm         401519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         401519       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4151919       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4151919       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4151919-75       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       <	4101516-45	Reverse cutting 3/8	16 mm	1	5/0	45 cm
4151516-75       Reverse cutting 3/8       16 mm       1.5       4/0       75 cm         4201516       Reverse cutting 3/8       16 mm       2       3/0       45 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4101516VB       Reverse cutting 3/8       16 mm black       1       5/0       75 cm         4071519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         4151919       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4301519       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       24 mm       3       2/0		Reverse cutting 3/8	16 mm	1	5/0	75 cm
4151516-75       Reverse cutting 3/8       16 mm       1.5       4/0       75 cm         4201516       Reverse cutting 3/8       16 mm       2       3/0       45 cm         4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4101516VB       Reverse cutting 3/8       16 mm black       1       5/0       75 cm         4071519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         41519       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         41519       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4301519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       24 mm       3       2/0	4151516	Reverse cutting 3/8	16 mm	1.5	4/0	
4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4101516VB       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4071519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         4151519.75       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4201519       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm bacy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4201524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301536       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3.5       0       75 c	4151516-75	Reverse cutting 3/8	16 mm	1.5	4/0	
4201516-75       Reverse cutting 3/8       16 mm       2       3/0       75 cm         4101516VB       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4071519       Reverse cutting 3/8       19 mm       0.7       6/0       75 cm         4101519       Reverse cutting 3/8       19 mm       1       5/0       75 cm         4151519.75       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4201519       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm bacy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4201524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301536       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3.5       0       75 c	4201516	Reverse cutting 3/8	16 mm	2	3/0	45 cm
4071519Reverse cutting 3/819 mm0.76/075 cm4101519Reverse cutting 3/819 mm15/075 cm4151519Reverse cutting 3/819 mm1.54/045 cm4151519.75Reverse cutting 3/819 mm1.54/075 cm4201519Reverse cutting 3/819 mm23/075 cm4301519Reverse cutting 3/819 mm32/075 cm4301519Reverse cutting 3/819 mm Easy Pass1.54/075 cm4101519VBReverse cutting 3/819 mm black15/075 cm4101519VBReverse cutting 3/819 mm black1.54/075 cm4101519VBReverse cutting 3/824 mm23/075 cm4201524Reverse cutting 3/824 mm3.5075 cm4301524Reverse cutting 3/836 mm3.5075 cm4351526Reverse cutting 3/836 mm3.5075 cm4351536Reverse cutting 3/836 mm4175 cm440147Reverse cutting 1/217 mm1.54/075 cm440147Reverse cutting 1/220 mm23/075 cm4401420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/	4201516-75	Reverse cutting 3/8	16 mm	2		75 cm
4071519Reverse cutting 3/819 mm0.76/075 cm4101519Reverse cutting 3/819 mm15/075 cm4151519Reverse cutting 3/819 mm1.54/045 cm4151519.75Reverse cutting 3/819 mm1.54/075 cm4201519Reverse cutting 3/819 mm23/075 cm4301519Reverse cutting 3/819 mm32/075 cm4301519Reverse cutting 3/819 mm Easy Pass1.54/075 cm4101519VBReverse cutting 3/819 mm black15/075 cm4101519VBReverse cutting 3/819 mm black1.54/075 cm4101519VBReverse cutting 3/824 mm23/075 cm4201524Reverse cutting 3/824 mm3.5075 cm4301524Reverse cutting 3/836 mm3.5075 cm4351526Reverse cutting 3/836 mm3.5075 cm4351536Reverse cutting 3/836 mm4175 cm440147Reverse cutting 1/217 mm1.54/075 cm440147Reverse cutting 1/220 mm23/075 cm4401420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/220 mm23/075 cm4101420Reverse cutting 1/	4101516VB	Reverse cutting 3/8	16 mm black	1	5/0	75 cm
4151519       Reverse cutting 3/8       19 mm       1.5       4/0       45 cm         4151519-75       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4151519-75       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301524       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       4       1       75 cm         4301536       Reverse cutting 3/8       36 mm       4       1       75 cm         4151417       Reverse cutting 1/2       17 mm       1.5       4/0       75	4071519	Reverse cutting 3/8		0.7		75 cm
4151519-75       Reverse cutting 3/8       19 mm       1.5       4/0       75 cm         4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         41515197-75       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4351524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4351536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4351536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         440137       Reverse cutting 1/2       17 mm       1.5       4/0	4101519	Reverse cutting 3/8	19 mm	1	5/0	75 cm
4201519       Reverse cutting 3/8       19 mm       2       3/0       75 cm         4301519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4151519P-75       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4201524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301536       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3       2/0       75 cm         4301536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401540       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4101420       Reverse cutting 1/2       20 mm       1.5       4/0       75 c	4151519	Reverse cutting 3/8	19 mm	1.5		45 cm
4301519       Reverse cutting 3/8       19 mm       3       2/0       75 cm         4151519P-75       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4351519VB       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4351524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4351536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4351536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4401547       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4151417       Reverse cutting 1/2       20 mm       2       5/0       75 c	4151519-75	Reverse cutting 3/8	19 mm	1.5	4/0	75 cm
4151519P-75       Reverse cutting 3/8       19 mm Easy Pass       1.5       4/0       75 cm         4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4401536       Reverse cutting 3/8       36 mm       4.1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4401547       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4401420       Reverse cutting 1/2       20 mm       2       5/0       75 cm         4101420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm <td>4201519</td> <td>Reverse cutting 3/8</td> <td>19 mm</td> <td>2</td> <td>3/0</td> <td>75 cm</td>	4201519	Reverse cutting 3/8	19 mm	2	3/0	75 cm
4101519VB       Reverse cutting 3/8       19 mm black       1       5/0       75 cm         4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3       2/0       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4151417       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4101420       Reverse cutting 1/2       20 mm       2       5/0       75 cm         4101420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm	4301519	Reverse cutting 3/8	19 mm	3	2/0	
4151519VB       Reverse cutting 3/8       19 mm black       1.5       4/0       75 cm         4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4351524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3       2/0       75 cm         4351526       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4401536       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4151417       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4101420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm         4151420       Reverse cutting 1/2       20 mm       2       3/0       75 cm <tr< td=""><td>4151519P-75</td><td>Reverse cutting 3/8</td><td>19 mm Easy Pass</td><td>1.5</td><td>4/O</td><td>75 cm</td></tr<>	4151519P-75	Reverse cutting 3/8	19 mm Easy Pass	1.5	4/O	75 cm
4201524       Reverse cutting 3/8       24 mm       2       3/0       75 cm         4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4351524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3       2/0       75 cm         4301536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4351536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4401536       Reverse cutting 1/2       17 mm       2       3/0       75 cm         44201417       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4101420       Reverse cutting 1/2       20 mm       2       5/0       75 cm         4151420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm         4151420       Reverse cutting 1/2       20 mm       2       3/0       75 cm         4151422       Reverse cutting 1/2       20 mm       2       3/0       75 cm         <	4101519VB	Reverse cutting 3/8	19 mm black	1	5/0	75 cm
4301524       Reverse cutting 3/8       24 mm       3       2/0       75 cm         4351524       Reverse cutting 3/8       24 mm       3.5       0       75 cm         4301536       Reverse cutting 3/8       36 mm       3       2/0       75 cm         4351536       Reverse cutting 3/8       36 mm       3.5       0       75 cm         4401536       Reverse cutting 3/8       36 mm       4       1       75 cm         4401536       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4151417       Reverse cutting 1/2       17 mm       1.5       4/0       75 cm         4201417       Reverse cutting 1/2       17 mm       2       3/0       75 cm         4101420       Reverse cutting 1/2       20 mm       2       5/0       75 cm         4101420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm         4151420       Reverse cutting 1/2       20 mm       1.5       4/0       75 cm         4201420       Reverse cutting 1/2       20 mm       2       3/0       75 cm         4201420       Reverse cutting 1/2       26 mm       2       3/0       75 cm	4151519VB	Reverse cutting 3/8	19 mm black	1.5	4/0	75 cm
4351524Reverse cutting 3/824 mm3.5075 cm4301536Reverse cutting 3/836 mm32/075 cm4351536Reverse cutting 3/836 mm3.5075 cm4401536Reverse cutting 3/836 mm4175 cm4401536Reverse cutting 1/217 mm1.54/075 cm4151417Reverse cutting 1/217 mm23/075 cm4201417Reverse cutting 1/217 mm25/075 cm4101420Reverse cutting 1/220 mm25/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4201524	Reverse cutting 3/8	24 mm	2	3/0	75 cm
4301536Reverse cutting 3/836 mm32/075 cm4351536Reverse cutting 3/836 mm3.5075 cm4401536Reverse cutting 3/836 mm4175 cm4151417Reverse cutting 1/217 mm1.54/075 cm4201417Reverse cutting 1/217 mm23/075 cm4101420Reverse cutting 1/220 mm25/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
4351536Reverse cutting 3/836 mm3.5075 cm4401536Reverse cutting 3/836 mm4175 cm4151417Reverse cutting 1/217 mm1.54/075 cm4201417Reverse cutting 1/217 mm23/075 cm4101420Reverse cutting 1/220 mm25/075 cm4101420Reverse cutting 1/220 mm1.54/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm					0	
4351536Reverse cutting 3/836 mm3.5075 cm4401536Reverse cutting 3/836 mm4175 cm4151417Reverse cutting 1/217 mm1.54/075 cm4201417Reverse cutting 1/217 mm23/075 cm4101420Reverse cutting 1/220 mm25/075 cm4101420Reverse cutting 1/220 mm1.54/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201420Reverse cutting 1/220 mm23/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4301536	Reverse cutting 3/8	36 mm		2/0	75 cm
4401536Reverse cutting 3/836 mm4175 cm4151417Reverse cutting 1/217 mm1.54/075 cm4201417Reverse cutting 1/217 mm23/075 cm4101420Reverse cutting 1/220 mm25/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4151422Reverse cutting 1/222 mm1.54/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4351536	Reverse cutting 3/8	36 mm	3.5	0	75 cm
4201417Reverse cutting 1/217 mm23/075 cm4101420Reverse cutting 1/220 mm25/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4151422Reverse cutting 1/220 mm1.54/075 cm4151422Reverse cutting 1/222 mm1.54/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4401536	Reverse cutting 3/8	36 mm		1	75 cm
4101420Reverse cutting 1/220 mm25/075 cm4151420Reverse cutting 1/220 mm1.54/075 cm4201420Reverse cutting 1/220 mm23/075 cm4151422Reverse cutting 1/222 mm1.54/075 cm4201426Reverse cutting 1/226 mm23/075 cm4301426Reverse cutting 1/226 mm32/075 cm4152120Straight cutting20 mm1.54/075 cm	4151417	Reverse cutting 1/2	17 mm	1.5	4/0	75 cm
4151420         Reverse cutting 1/2         20 mm         1.5         4/0         75 cm           4201420         Reverse cutting 1/2         20 mm         2         3/0         75 cm           4151422         Reverse cutting 1/2         22 mm         1.5         4/0         75 cm           4151422         Reverse cutting 1/2         22 mm         1.5         4/0         75 cm           4201426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm						
4201420         Reverse cutting 1/2         20 mm         2         3/0         75 cm           4151422         Reverse cutting 1/2         22 mm         1.5         4/0         75 cm           4201426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm	4101420	Reverse cutting 1/2	20 mm	2	5/0	75 cm
4151422         Reverse cutting 1/2         22 mm         1.5         4/0         75 cm           4201426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm	4151420	Reverse cutting 1/2	20 mm	1.5	4/0	75 cm
4151422         Reverse cutting 1/2         22 mm         1.5         4/0         75 cm           4201426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm	4201420	Reverse cutting 1/2	20 mm	2	3/0	75 cm
4201426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm	4151422	Reverse cutting 1/2	22 mm			75 cm
4301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           4152120         Straight cutting         20 mm         1.5         4/0         75 cm	4201426	Reverse cutting 1/2	26 mm		3/0	75 cm
	4301426			3	2/0	
	4152120	Straight cutting	20 mm	1.5	4/0	75 cm
	4152125		25 mm		4/0	

#### Features:

| twisted and coated multifilament; non-absorbable; composition: polyamide silicone coating; colour: white; box with 12 sutures;

#### Advantages:

excellent tissue support; coated multifilament | minimal tissue reaction; sythetic, excellent handling properties, outstanding strenght retention from extended wound support; high pliability and softness; high tensile strength; excellent knot security; | soft passage through the tissues without sawing, tissue drag and trauma; | guaranteed seal and product sterility.

#### Indications:

| general soft tissue closing and/or ligation; skin closure

#### In cassettes:

	Thread thickness		
Cat. No.	EP	USP	- Lenght
3015	1.5	4/0	100 m
3020	2	3/0	100 m
3025	2.5	2/0	100 m
3030	3	2/0	100 m
3035	3.5	0	100 m
3040	4	1	100 m
3050	5	2	75 m
3060	6	3&4	50 m
3070	7	5	40 m
3080	8	6	40 m
3090	9	7	30 m
3100	10	8	30 m







#### Surgicryl Monofast Suture

#### Features:

| monofilament;

- absorbable;
- | composition: Polyglecaprone 25, a copolymer made of 75% glycolide and 25%  $\epsilon$ -caprolactone
- no coating;
- | colour: violet/beige/ undyed;
- box with 12 sutures;
- The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.

#### Advantages:

excellent tissue support;
Monofilament with microscopic uniform diameter
minimal tissue reaction;
sythetic, coated multifilament;
breaks down by hydrolysis
excellent handling properties;
high flexibility;
high tensile strength;
excellent knot security;
soft passage through the tissues without sawing, tissue drag and trauma;
excellent knot tie down;
guaranteed seal and product sterility.

#### Indications:

superficial soft tissue closing and/or ligation especially in general surgery;gastrointestinal surgery;

- gynaecology;
- urology;
- | plastic surgery;

#### Contraindications:

| cardiovascular surgery, ophthalmic surgery

- | neurosurgery being absorbable, not to use when an extended wound support is required;
- caution in patients with delayed wound healing (e.g. compromised immune system), in the region of infected wounds

#### Sutures without needle:

Cat. No.	Thread thickness		- Lowerlat	
Cal. NO.	EP	USP	<sup>—</sup> Lenght	
16300015	3	2/0	150 cm	
14350015	3.5	0	150 cm	
16400015	4	1	150 cm	
16500012	5	2	120 cm	
16500015	5	2	150 cm	
16500025	5	2	250 cm ( Ligapack)	



#### In cassettes:

Cat. No.		hread thickness	
Cat. NO.	EP	USP	Lengin
16151	1.5	4/0	15 m
16201	2	3/0	15 m
16202	2	3/0	25 m
16301	3	2/0	15 m
16302	3	2/0	25 m
16351	3.5	0	15 m
16352	3.5	0	25 m
16401	4	1	15 m
16402	4	1	25 m
16405	4	1	50 m
16501	5	2	15 m
16502	5	2	25 m
16505	5	2	50 m
165010	5	2	100 m
16601	6	3&4	15 m
16602	6	3&4	25 m
16605	6	3&4	50 m

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
16070113	Round 1/2	13 mm	0.7	6/0	45 cm
6100113	Round 1/2	13 mm	1	5/0	45 cm
6150113	Round 1/2	13 mm	1.5	4/0	75 cm
6100117	Round 1/2	17 mm	1	5/0	45 cm
6150117	Round 1/2	17 mm	1.5	4/0	75 cm
6200117	Round 1/2	17 mm	2	3/0	75 cm
6150122	Round 1/2	22 mm	1.5	4/0	75 cm
6200122	Round 1/2	22 mm	2	3/0	75 cm
6300122	Round 1/2	22 mm	3	2/0	75 cm
6150126	Round 1/2	26 mm	1.5	4/0	75 cm
6200126	Round 1/2	26 mm	2	3/0	75 cm
6300126	Round 1/2	26 mm	3	2/0	75 cm
6350126	Round 1/2	26 mm	3.5	0	75 cm
6200130	Round 1/2	30 mm	2	3/0	75 cm
16300130	Round 1/2	30 mm	3	2/0	75 cm
6350130	Round 1/2	30 mm	3.5	0	75 cm
6400130	Round 1/2	30 mm	4	1	75 cm
6300136	Round 1/2	36 mm	3	2/0	75 cm
	Round 1/2		3.5	0	90 cm
6350136		36 mm		1	
6400136	Round 1/2	36 mm	4		<u>90 cm</u>
6300140	Round 1/2	40 mm	3	2/0	75 cm
6350140	Round 1/2	40 mm	3.5	0	90 cm
16400140	Round 1/2	40 mm	4	1	90 cm
6350148	Round 1/2	48 mm	3.5	0	<u>90 cm</u>
16400148	Round 1/2	48 mm	4	1	90 cm
6300165	Round 1/2	65 mm	3	2/0	90 cm
16071512	Reverse cutting 3/8	12 mm- niebarwiona	0.7	6/0	
16071512-75	Reverse cutting 3/8	12 mm	0.7	6/0	45 cm
16071512V	Reverse cutting 3/8	12 mm Easy Pass	0.7	6/0	75 cm
16101512	Reverse cutting 3/8	12 mm Easy Pass	1	5/0	75 cm
16101512V	Reverse cutting 3/8	12 mm Easy Pass	1	5/0	45 cm
6071512UP	Reverse cutting 3/8	16 mm	0.7	6/0	45 cm
16101512UP	Reverse cutting 3/8	16 mm	1	5/0	50 cm
6151512UP	Reverse cutting 3/8	16 mm	1.5	4/0	50 cm
6101516	Reverse cutting 3/8	16 mm Easy Pass	1	5/0	50 cm
6151516	Reverse cutting 3/8	16 mm Easy Pass	1.5	4/0	75 cm
6201516	Reverse cutting 3/8	16 mm Black	2	3/0	75 cm
6101516P	Reverse cutting 3/8	19 mm	1	5/0	75 cm
I615156P	Reverse cutting 3/8	19 mm- niebarwiona	1.5	4/0	50 cm
6101516VB	Reverse cutting 3/8	19 mm	1	5/0	50 cm
6101519V	Reverse cutting 3/8	19 mm- niebarwiona	1	5/0	75 cm
6101519	Reverse cutting 3/8	19 mm	1	5/0	75 cm
6151519V	Reverse cutting 3/8	19 mm	1.5	4/0	75 cm
6151519	Reverse cutting 3/8	19 mm - niebarwiona	1.5	4/0	75 cm
6201519V	Reverse cutting 3/8	19 mm Easy Pass- niebarwiona	2	3/0	75 cm
6201519	Reverse cutting 3/8	19 mm Easy Pass- niebarwiona	2	3/0	75 cm
6301519	Reverse cutting 3/8	19 mm czarna	3	2/0	75 cm
6151519P	Reverse cutting 3/8	19 mm czarna	1.5	4/0	75 cm
6251519P	Reverse cutting 3/8	24 mm	2	3/0	50 cm
6201519P	Reverse cutting 3/8	24 mm	1	5/0	<u>50 cm</u>
6151519VB	Reverse cutting 3/8	24 mm Easy Pass	1.5	4/0	75 cm
6201524	Reverse cutting 3/8	26 mm	2	3/0	75 cm
6301524	Reverse cutting 3/8	26 mm	3	2/0	75 cm
6201524P-75	Reverse cutting 3/8	26 mm	2	3/0	75 cm
6151526	Reverse cutting 3/8	26 mm	1.5	4/0	75 cm
6201526	Reverse cutting 3/8	26 mm	2	3/0	75 cm
6201526V	Reverse cutting 3/8	30 mm	2	3/0	75 cm
6301526	Reverse cutting 3/8	40 mm	3	2/0	75 cm
6301526V	Reverse cutting 3/8	40 mm	3	2/0	75 cm
6351530V	Reverse cutting 1/2	36 mm	3.5	0	75 cm
6351540	Reverse cutting 1/2	36 mm	3.5	0	75 cm
6401540	Reverse cutting 1/2	40 mm	4	1	75 cm
6202160	Straight cutting	60 mm	2	3/0	75 cm
6302160	Straight cutting	60 mm	3	2/0	75 cm
6351840	Blunt 1/2	40 mm	3.5	0	90 cm
	Blunt 1/2		4	1	90 cm

#### Suture



#### Surgicryl Monofilament Suture

#### Features:

#### | monofilament;

absorbable;

composition: polydioxanone, a polymer made from polyester poly (p-dioxanone)

- no coating;
- colour: violet;
- box with 12 sutures;
- | The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.

#### Advantages:

excellent tissue support;
Monofilament with microscopic uniform diameter
minimal tissue reaction;
sythetic, coated multifilament;
breaks down by hydrolysis
excellent handling properties, outstanding strenght retention from extended wound support;
high flexibility;
high tensile strength;
excellent knot security;
soft passage through the tissues without sawing, tissue drag and trauma;
excellent knot tie down;
guaranteed seal and product sterility.

#### Indications:

- general soft tissue closing and/or ligation especially in general surgery;
- | gastrointestinal surgery;
- gynaecology;
- urology;
- | plastic surgery;

#### **Contraindications**:

| cardiovascular surgery, prosthetic devices;

- | neurosurgery being absorbable, not to use when an extended wound support is required;
- caution in patients with delayed wound healing (e.g. compromised immune system), in the region of infected wounds

#### Sutures without needle:

		Thread thickness		
Cat. No.	EP	USP	Lenght	
13500015	5	2	150 cm	
13600015	6	3&4	150 cm	



#### In cassettes:

Cat. No.		hread thickness	
Cal. NO.	EP	USP	Lengin
13151	1.5	4/0	15 m
13201	2	3/0	15 m
13202	2	3/0	25 m
13301	3	2/0	15 m
13302	3	2/0	25 m
13351	3.5	0	15 m
13352	3.5	0	25 m
13401	4	1	15 m
13402	4	1	25 m
13405	4	1	50 m
13501	5	2	15 m
13502	5	2	25 m
13505	5	2	50 m
135010	5	2	100 m
13601	6	3&4	15 m
13602	6	3&4	25 m
13605	6	3&4	50 m

#### Sutures with needle:

Cat. No.	Needle shape	Needle size
13070113	Round 1/2	13 mm
13100113	Round 1/2	13 mm
13070117	Round 1/2	17 mm
3100117	Round 1/2	17 mm
3150117	Round 1/2	17 mm
3200117	Round 1/2	17 mm
3150122	Round 1/2	22 mm
3200122	Round 1/2	22 mm
13150126	Round 1/2	26 mm
13200126	Round 1/2	26 mm
13300126	Round 1/2	26 mm
3350126	Round 1/2	26 mm
13200130	Round 1/2	30 mm
13300130	Round 1/2	30 mm
13350130	Round 1/2	30 mm
13300136	Round 1/2	36 mm
13350136	Round 1/2	36 mm
13400136	Round 1/2	36 mm
3300140	Round 1/2	40 mm
3350140	Round 1/2	40 mm
13400140	Round 1/2	40 mm
13500140	Round 1/2	40 mm
13400148	Round 1/2	48 mm
13500148	Round 1/2	48 mm
13100113-2	Round ½- double	2 x 13 mm
13100117-2	Round ½- double	2 x 17 mm
13150117-2	Round ½- double	2 x 17 mm
13350940	Round ½ loop	40 mm
13400940	Round ½ loop	40 mm
13350950	Round ½ loop	50 mm
13400950	Round ½ loop	50 mm
13500950	Round ½ loop	50 mm
13400965	Round ½ loop	65 mm
13071512	Reverse cutting 3/8	12 mm
13101512	Reverse cutting 3/8	12 mm
13151512	Reverse cutting 3/8	12 mm
13101516	Reverse cutting 3/8	16 mm
13151516	Reverse cutting 3/8	16 mm
13101516VB	Reverse cutting 3/8	16 mm Black
13101519	Reverse cutting 3/8	19 mm
13101519-75	Reverse cutting 3/8	
13151519	Reverse cutting 3/8	19 mm
13201519		<u>19 mm</u> 19 mm
	Reverse cutting 3/8	
13151519VB	Reverse cutting 3/8	19 mm - black
13101516P	Reverse cutting 3/8	19 mm Easy Pass
13151524	Reverse cutting 1/2	24 mm
13201524	Reverse cutting 1/2	24 mm
13301524	Reverse cutting 1/2	24 mm
13351524	Reverse cutting 1/2	24 mm
13201526	Reverse cutting 1/2	26 mm
13301526	Reverse cutting 1/2	26 mm
13301530	Reverse cutting 1/2	30 mm
13351530	Reverse cutting 1/2	30 mm
13401540	Reverse cutting 1/2	40 mm
13501540	Reverse cutting 1/2	40 mm
13301436	Reverse cutting 1/2	
13351436	Reverse cutting 1/2	36 mm
13401436	Reverse cutting 1/2	36 mm
13401436	Reverse cutting 1/2	36 mm

#### Suture



Metric	USP	Lenght
0.7	6/0	75 cm
1	5/0	75 cm
0.7	6/0	75 cm
1	5/0	75 cm
1.5	4/0	75 cm
2	3/0	75 cm
1.5	4/0	75 cm
2	3/0	75 cm
1.5	4/0	75 cm
2	3/0	75 cm
3	2/0	75 cm
3.5	0	75 cm
2	3/0	75 cm
3	2/0	75 cm
3.5	0	75 cm
3	2/0	75 cm
3.5	0	90 cm
4	1	90 cm
3	2/0	75 cm
3.5	0	90 cm
4	1	90 cm
5	2	90 cm
4	1	90 cm
5	2	90 cm
1	5/0	75 cm
1	5/0	90 cm
1.5	4/0	75 cm
3.5	0	150 cm
4	1	150 cm
3.5	0	150 cm
4	1	150 cm
5	2	150 cm
4	1	150 cm
0.7	6/0	45 cm
1	5/0	45 cm
1.5	4/0	75 cm
1	5/0	75 cm
1.5	4/0	75 cm
1	5/0	75 cm
1	5/0	45 cm
1	5/0	75 cm
1.5	4/0	75 cm
2	3/0	75 cm
1.5	4/0	75 cm
1	5/0	50 cm
1.5	4/0	75 cm
2 3	3/0	75 cm
	2/0	75 cm
3.5	0 3/0	75 cm
2 3	2/0	75 cm
3	2/0	75 cm 75 cm
<u> </u>	0	75 cm
5	2	90 cm 90 cm
5	<u> </u>	75 cm
3	2/0	75 cm
3.5	0	75 cm
4	1	75 cm
4	1	150 cm
		100 0111

#### Surgicryl PGA Suture

#### Features:

Braided and coated multifilament;
absorbable;
composition: polyglycolic acid
coating: polycaprolactone and calcium stearate (1%)
colour: violet/beige/ undyed;
box with 12 sutures;
he hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days.
Approximately 50% of tensile strength remain after 21 days.

#### Advantages:

- excellent tissue support;
- minimal tissue reaction;
- sythetic, coated multifilament;
- breaks down by hydrolysis
- excellent handling properties;
- high flexibility;
- high tensile strength;
- excellent knot security;
- soft passage through the tissues without sawing, tissue drag and trauma;
- | excellent knot tie down;
- | guaranteed seal and product sterility.

#### Indications:

- | general soft tissue closing and/or ligation especially in
- general surgery;
- gastrointestinal surgery;
- gynaecology;
- obstetrics;
- urology;
- ophthalmic surgery;
- | plastic surgery;
- orthopaedics.

#### Contraindications:

- | cardiovascular surgery
- neurosurgery being absorbable, not to use when an extended wound support is required;
- caution in patients with delayed wound healing (e.g. compromised immune system), in the region of infected wounds.



#### In cassettes:

Cat. No.	Thread thickness		-longht
Cdl. NO.	EP	USP	-Lenght
11151	1.5	4/0	15 m
11201	2	3/0	15 m
11202	2	3/0	25 m
11301	3	2/0	15 m
11302	3	2/0	25 m
11305	3	2/0	50 m
11351	3.5	0	15 m
11352	3.5	0	25 m
11355	3.5	0	50 m
11401	4	1	15 m
11402	4	1	25 m
11405	4	1	50 m
11501	5	2	15 m
11502	5	2	25 m
11505	5	2	50 m
115010	5	2	100 m
17601	6	3&4	15 m
17602	6	3&4	25 m
17605	6	3&4	50 m
116010	6	3&4	100m
11701	7	5	15 m
11702	7	5	25 m
11705	7	5	50 m
11801	8	6	15 m
11802	8	6	25 m
11805	8	6	50 m

#### Suture



#### ABSORBABLE SUTURES IN CASETTES

#### Sutures without needle:

		Thread thickness	
Cat. No.	EP	USP	<sup>—</sup> Lenght
11200015	2	3/0	150 cm
11300015	3	2/0	150 cm
11350015	3.5	0	150 cm
11400015	4	1	150 cm
11500015	5	2	150 cm
11600015	6	3&4	150 cm
11700015	7	5	150 cm
11800015	8	6	150 cm
11550050	1.5	4/O	3 x 50 cm
11200050	2	3/0	3 x 50 cm
11300050	3	2/0	3 x 50 cm
11350050	3.5	0	3 x 50 cm
11400050	4	1	3 x 50 cm
11500050	5	2	3 x 50 cm
11200045	2	3/0	12 x 45 cm
11300045	3	2/0	12 x 45 cm
11350045	3.5	0	12 x 45 cm
11400045	4	1	12 x 45 cm
11500045	5	2	12 x 45 cm
1150045-6	1.5	4/O	6 x 45 cm
11200045-6	2	3/0	6 x 45 cm
11300045-6	3	2/0	6 x 45 cm
11350045-6	3.5	0	6 x 45 cm
11400045-6	4	1	6 x 45 cm
11500045-6	5	2	6 x 45 cm
11200060	2	3/0	10 x 60 cm
11300060	3	2/0	10 x 60 cm
11350060	3.5	0	10 x 60 cm
11200070	2	3/0	5 x 70 cm
11300070	3	2/0	5 x 70 cm
11350070	3.5	0	5 x 70 cm
11400070	4	1	5 x 70 cm
11500070	5	2	5 x 70 cm
11150025	1.5	4/O	250 cm
11200025	2	3/0	250 cm
11300025	3	2/0	250 cm
11350025	3.5	0	250 cm
11400025	4	1	250 cm
11500025	5	2	250 cm
11600025	6	3&4	250 cm

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
11070113	Round 1/2	13 mm	0.7	6/0	75 cm
11100113	Round 1/2	13 mm	1	5/0	75 cm
11150113	Round 1/2	13 mm	1.5	4/0	75 cm
1070117	Round 1/2	17 mm	0.7	6/0	75 cm
1150117	Round 1/2	17 mm	1.5	4/0	75 cm
1200117	Round 1/2	17 mm	2	3/0	75 cm
1300117 1150122	Round 1/2 Round 1/2	17 mm	<u>3</u> 1.5	2/0 4/0	75 cm 75 cm
1200122	Round 1/2	22 mm 22 mm	2	3/0	75 cm
11300122	Round 1/2	22 mm	3	2/0	75 cm
11150126	Round 1/2	26 mm	1.5	4/0	75 cm
1200126	Round 1/2	26 mm	2	3/0	75 cm
1300126	Round 1/2	26 mm	3	2/0	75 cm
1350126	Round 1/2	26 mm	3.5	0	75 cm
1400126	Round 1/2	26 mm	4	2	75 cm
15500126 11200130	Round 1/2 Round 1/2	26 mm 30 mm	5	3/0	75 cm 75 cm
1300130	Round 1/2	30 mm	3	2/0	75 cm
1350130	Round 1/2	30 mm	3.5	0	75 cm
1400130	Round 1/2	30 mm	4	1	75 cm
1200136	Round 1/2	36 mm	2	3/0	75 cm
1300136	Round 1/2	36 mm	3	2/0	75 cm
1350136	Round 1/2	36 mm	3.5	0	90 cm
1400136 1300140	Round 1/2 Round 1/2	36 mm 40 mm	4 3	1 2/0	90 cm 75 cm
1300140	Round 1/2 Round 1/2	40 mm 40 mm	3.5	0	75 cm 90 cm
1400140-75		40 mm	4	1	75 cm
1400140	Round 1/2	40 mm	4	1	90 cm
1500140-75		40 mm	5	2	75 cm
1500140	Round 1/2	40 mm	5	2	90 cm
1300148	Round 1/2	48 mm	3	2/0	75 cm
11350148	Round 1/2	48 mm	3.5	0	90 cm
11400148	Round 1/2	48 mm	4	2	90 cm
11500148 11600148	Round 1/2 Round 1/2	48 mm 48 mm	5	3&4	90 cm 90 cm
11400950	Round ½ loop	50 mm	4	]	150 cm
11300327	Round 5/8	27 mm	3	2/0	75 cm
11200217	Round 5/8	17 mm	2	3/0	75 cm
11150217	Round 5/8	17 mm	1.5	4/0	75 cm
11071512	Reverse cutting 3/8	12 mm	0.7	6/0	45 cm
11101512	Reverse cutting 3/8	12 mm	1	5/0	45 cm
11151512	Reverse cutting 3/8	12 mm	1.5 0.7	4/0 6/0	45 cm 50 cm
11071512UP 11101512UP	Reverse cutting 3/8 Reverse cutting 3/8	12 mm Easy Pass 12 mm Easy Pass	1	5/0	50 cm
11151512UP	Reverse cutting 3/8	12 mm Easy Pass	1.5	4/0	50 cm
11071516	Reverse cutting 3/8	16 mm	0.7	6/0	75 cm
11101516	Reverse cutting 3/8	16 mm	1	5/0	75 cm
11151516	Reverse cutting 3/8	16 mm	1.5	4/0	75 cm
11201516	Reverse cutting 3/8	16 mm	2	3/0	75 cm
11151516UP	Reverse cutting 3/8	16 mm Easy Pass	2	3/0	50 cm
11101519	Reverse cutting 3/8	19 mm	1.5	5/0 4/0	75 cm
1151519 11201519	Reverse cutting 3/8 Reverse cutting 3/8	19 mm 19 mm	2	3/0	75 cm 75 cm
11201319 11151519UP	Reverse cutting 3/8	19 mm Easy Pass	1.5	4/0	50 cm
5151519VB	Reverse cutting 3/8	19 mm czarna	1.5	4/0	75 cm
1201524	Reverse cutting 3/8	24 mm	2	3/0	75 cm
1301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
1351524	Reverse cutting 3/8	24 mm	3.5	0	75 cm
1201524UP	Reverse cutting 3/8	24 mm Easy Pass	2	3/0	50 cm
1151530	Reverse cutting 3/8	30 mm	1.5	4/0	75 cm
1201530	Reverse cutting 3/8 Reverse cutting 3/8	30 mm 30 mm	2 3	3/0 2/0	75 cm 75 cm
1351530	Reverse cutting 3/8	30 mm	3.5	0	75 cm
		30 mm	4	1	75 cm
1401530	Reverse cutting 3/8				
	Reverse cutting 3/8 Reverse cutting 3/8	30 mm	5	2	75 cm
1501530			5 3	2 2/0	75 cm
11501530 11301536 11351536	Reverse cutting 3/8 Reverse cutting 3/8 Reverse cutting 3/8	30 mm 36 mm 36 mm	3 3.5	2/0 0	75 cm 75 cm
11401530 11501530 11301536 11351536 11401536 11401536	Reverse cutting 3/8 Reverse cutting 3/8	30 mm 36 mm	3	2/0	75 cm

1401540         Reverse cutting 3/8         40 mm         4         1         75 cm           11501540         Reverse cutting 1/2         22 mm         2         3/0         75 cm           11201422         Reverse cutting 1/2         22 mm         3         2/0         75 cm           11301426         Reverse cutting 1/2         26 mm         1.5         4/0         75 cm           11301426         Reverse cutting 1/2         26 mm         2         3/0         75 cm           11301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           11301426         Reverse cutting 1/2         30 mm         5         2         90 cm           11301436         Reverse cutting 1/2         36 mm         3         2/0         90 cm           11301436         Reverse cutting 1/2         36 mm         4         1         90 cm           11301436         Reverse cutting 1/2         36 mm         5         2         90 cm           11301440         Reverse cutting 1/2         40 mm         3         2/0         90 cm           11301440         Reverse cutting 1/2         40 mm         5         2         90 cm           11301440	Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
11201422       Reverse cutting 1/2       22 mm       3/0       75 cm         11301422       Reverse cutting 1/2       22 mm       3       2/0       75 cm         11301426       Reverse cutting 1/2       26 mm       1.5       4/0       75 cm         11301426       Reverse cutting 1/2       26 mm       2       3/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3       2/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3       2/0       75 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301436       Reverse cutting 1/2       40 mm       4       1       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       4       1       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301448 <td>11401540</td> <td>Reverse cutting 3/8</td> <td>40 mm</td> <td>4</td> <td>1</td> <td>75 cm</td>	11401540	Reverse cutting 3/8	40 mm	4	1	75 cm
11301422       Reverse cutting 1/2       22 mm       3       2/0       75 cm         11151426       Reverse cutting 1/2       26 mm       1.5       4/0       75 cm         11301426       Reverse cutting 1/2       26 mm       2       3/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3       2/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3.5       0       75 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301436       Reverse cutting 1/2       36 mm       3.5       0       90 cm         11301436       Reverse cutting 1/2       36 mm       3.5       0       90 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       5.5       0       90 cm         11401440       Reverse cutting 1/2       40 mm       5.5       0       90 cm         11401440       Reverse cutting 1/2       48 mm       1       90 cm         11501444	11501540	Reverse cutting 3/8	40 mm	5	2	75 cm
11151426       Reverse cutting 1/2       26 mm       1.5       4/0       75 cm         11201426       Reverse cutting 1/2       26 mm       2       3/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3       2/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3.5       0       75 cm         11301430       Reverse cutting 1/2       30 mm       5       2       90 cm         11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11301436       Reverse cutting 1/2       36 mm       4       1       90 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       48 mm       35       0       90 cm	11201422	Reverse cutting 1/2	22 mm	2	3/0	75 cm
11201426       Reverse cutting 1/2       26 mm       2       3/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3       2/0       75 cm         11301426       Reverse cutting 1/2       26 mm       3.5       0       75 cm         11501430       Reverse cutting 1/2       30 mm       5       2       90 cm         11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11301436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3.5       0       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301448       Reverse cutting 1/2       48 mm       4       1       90 cm         11301448       Reverse cutting 1/2       48 mm       5       2       90 cm         1	11301422	Reverse cutting 1/2	22 mm	3	2/0	75 cm
11301426         Reverse cutting 1/2         26 mm         3         2/0         75 cm           11351426         Reverse cutting 1/2         26 mm         35         0         75 cm           11501430         Reverse cutting 1/2         30 mm         5         2         90 cm           11301436         Reverse cutting 1/2         36 mm         3         2/0         90 cm           11301436         Reverse cutting 1/2         36 mm         4         1         90 cm           11301436         Reverse cutting 1/2         36 mm         4         1         90 cm           11301436         Reverse cutting 1/2         36 mm         5         2         90 cm           11301440         Reverse cutting 1/2         40 mm         3         2/0         90 cm           11301440         Reverse cutting 1/2         40 mm         5         2         90 cm           11301440         Reverse cutting 1/2         40 mm         5         2         90 cm           11301440         Reverse cutting 1/2         48 mm         35         0         90 cm           11301440         Reverse cutting 1/2         48 mm         5         2         90 cm           11301448         <	11151426	Reverse cutting 1/2	26 mm	1.5	4/O	75 cm
11351426       Reverse cutting 1/2       26 mm       3.5       0       75 cm         11501430       Reverse cutting 1/2       30 mm       5       2       90 cm         11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11351436       Reverse cutting 1/2       36 mm       3.5       0       90 cm         11401436       Reverse cutting 1/2       36 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3.5       0       90 cm         11401440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm	11201426	Reverse cutting 1/2	26 mm	2	3/0	75 cm
11501430         Reverse cutting 1/2         30 mm         5         2         90 cm           11301436         Reverse cutting 1/2         36 mm         3.5         0         90 cm           11351436         Reverse cutting 1/2         36 mm         3.5         0         90 cm           11401436         Reverse cutting 1/2         36 mm         4         1         90 cm           11501436         Reverse cutting 1/2         36 mm         5         2         90 cm           11301440         Reverse cutting 1/2         40 mm         3         2/0         90 cm           11301440         Reverse cutting 1/2         40 mm         3.5         0         90 cm           11401440         Reverse cutting 1/2         40 mm         4         1         90 cm           11501440         Reverse cutting 1/2         40 mm         5         2         90 cm           11501448         Reverse cutting 1/2         48 mm         3.5         0         90 cm           11501448         Reverse cutting 1/2         48 mm         5         2         90 cm           11501448         Reverse cutting 1/2         48 mm         5         7/0         90 cm           11501448	11301426	Reverse cutting 1/2	26 mm	3	2/0	75 cm
11301436       Reverse cutting 1/2       36 mm       3       2/0       90 cm         11351436       Reverse cutting 1/2       36 mm       4       1       90 cm         11401436       Reverse cutting 1/2       36 mm       5       2       90 cm         11501435       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501446       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11401443       Reverse cutting 1/2       48 mm       5       2       90 cm         11601448       Reverse cutting 1/2       48 mm       6       384       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         0SU	11351426	Reverse cutting 1/2	26 mm	3.5	0	75 cm
11351436       Reverse cutting 1/2       36 mm       3.5       0       90 cm         11401436       Reverse cutting 1/2       36 mm       4       1       90 cm         11501436       Reverse cutting 1/2       36 mm       5       2       90 cm         11501440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11351440       Reverse cutting 1/2       40 mm       3.5       0       90 cm         11351440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       4       1       90 cm         11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11501448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         0SU07	11501430	Reverse cutting 1/2	30 mm	5	2	90 cm
11401436       Reverse cutting 1/2       36 mm       4       1       90 cm         11501436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       35       0       90 cm         11401440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       38/4       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       38/4       90 cm         11202160       Straight cutting       60 mm       3       2/0       75 cm         0SU08	11301436	Reverse cutting 1/2	36 mm	3	2/0	90 cm
11501436       Reverse cutting 1/2       36 mm       5       2       90 cm         11301440       Reverse cutting 1/2       40 mm       3       2/0       90 cm         11301440       Reverse cutting 1/2       40 mm       3.5       0       90 cm         11401440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11501440       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm         11601448       Reverse cutting 1/2       48 mm       6       384       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         0SU07       Double spatula 3/8       2 x 6 mm.200 µm       0.5       7/0       45 cm         0SU08       Double spatula 3/8       2 x 6 cm.150 µm       0.2       10/0       30 cm	11351436	Reverse cutting 1/2	36 mm	3.5	0	90 cm
11301440         Reverse cutting 1/2         40 mm         3         2/0         90 cm           11351440         Reverse cutting 1/2         40 mm         3.5         0         90 cm           11401440         Reverse cutting 1/2         40 mm         4         1         90 cm           11501440         Reverse cutting 1/2         40 mm         5         2         90 cm           11501448         Reverse cutting 1/2         48 mm         35         0         90 cm           11501448         Reverse cutting 1/2         48 mm         4         1         90 cm           11501448         Reverse cutting 1/2         48 mm         6         3&4         90 cm           11501448         Reverse cutting 1/2         48 mm         6         3&4         90 cm           11202160         Straight cutting         60 mm         2         3/0         90 cm           11302160         Straight cutting         60 mm         3         2/0         75 cm           OSU07         Double spatula 3/8         2 x 6 mm.200 µm         0.5         7/0         45 cm           OSU08         Double spatula 3/8         2 x 6 mm.200 µm         0.4         8/0         30 cm           OSU0	11401436	Reverse cutting 1/2	36 mm	4	1	90 cm
11351440       Reverse cutting 1/2       40 mm       3.5       0       90 cm         11401440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11351448       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11401448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm         11501448       Reverse cutting 1/2       48 mm       6       384       90 cm         11202160       Straight cutting       60 mm       3       2/0       75 cm         0SU07       Double spatula 1/2       2 x 8 mm, 200 µm       0.5       7/0       45 cm         0SU08       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       30 cm         0SU08       Double spatula 3/8       2 x 6 mm, 200 µm       0.2       10/0       30 cm <td>11501436</td> <td>Reverse cutting 1/2</td> <td>36 mm</td> <td>5</td> <td>2</td> <td>90 cm</td>	11501436	Reverse cutting 1/2	36 mm	5	2	90 cm
11401440       Reverse cutting 1/2       40 mm       4       1       90 cm         11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11501440       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11501448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 µm       0.5       7/0       45 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.3       9/0       30 cm         OSU071       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.5       7/0 <td>11301440</td> <td>Reverse cutting 1/2</td> <td>40 mm</td> <td>3</td> <td>2/0</td> <td>90 cm</td>	11301440	Reverse cutting 1/2	40 mm	3	2/0	90 cm
11501440       Reverse cutting 1/2       40 mm       5       2       90 cm         11351448       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11401448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 µm       0.5       7/0       45 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6.2 mm, 150 µm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.5       7/0<	11351440	Reverse cutting 1/2	40 mm	3.5	0	90 cm
11351448       Reverse cutting 1/2       48 mm       3.5       0       90 cm         11401448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11501448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 µm       0.5       7/0       45 cm         OSU08-30       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6,2 mm, 150 µm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6,4 mm, 200 µm       0.5       7/0       45 cm         OSU091       Double spatula 3/8       2 x 6,4 mm, 200 µm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 6,4 mm, 200 µm       0.5 <td>11401440</td> <td>Reverse cutting 1/2</td> <td>40 mm</td> <td>4</td> <td>1</td> <td>90 cm</td>	11401440	Reverse cutting 1/2	40 mm	4	1	90 cm
11401448       Reverse cutting 1/2       48 mm       4       1       90 cm         11501448       Reverse cutting 1/2       48 mm       5       2       90 cm         11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 µm       0.5       7/0       45 cm         OSU08-30       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 µm       0.4       8/0       45 cm         OSU08       Double spatula 3/8       2 x 6.2 mm, 150 µm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.5       7/0       45 cm         OSU072       Double spatula 3/8       2 x 6.4 mm, 200 µm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 8 mm, 200 µm       0.5       7/0       20 cm         OSU071       Double spatula 1/4       2 x 8 mm, 300,350 µm<	11501440	Reverse cutting 1/2	40 mm	5	2	90 cm
11501448         Reverse cutting 1/2         48 mm         5         2         90 cm           11601448         Reverse cutting 1/2         48 mm         6         3&4         90 cm           11202160         Straight cutting         60 mm         2         3/0         90 cm           11202160         Straight cutting         60 mm         2         3/0         90 cm           11302160         Straight cutting         60 mm         3         2/0         75 cm           OSU07         Double spatula 1/2         2 x 8 mm, 200 μm         0.5         7/0         45 cm           OSU08-30         Double spatula 3/8         2 x 6 mm, 200 μm         0.4         8/0         30 cm           OSU08         Double spatula 3/8         2 x 6 mm, 200 μm         0.4         8/0         45 cm           OSU08         Double spatula 3/8         2 x 6,2 mm, 150 μm         0.2         10/0         30 cm           OSU09         Double spatula 3/8         2 x 6,4 mm, 200 μm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         20 cm           OSU06         Double spatula 3/8         2 x 8 mm, 300,350 μm         0.7         6/0<	11351448	Reverse cutting 1/2	48 mm	3.5	0	90 cm
11601448       Reverse cutting 1/2       48 mm       6       3&4       90 cm         11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 μm       0.5       7/0       45 cm         OSU08-30       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       45 cm         OSU08       Double spatula 3/8       2 x 6,2 mm, 150 μm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.3       9/0       30 cm         OSU072       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 8 mm, 200 μm       0.5       7/0       20 cm         OSU06       Double spatula 1/4       2 x 8 mm, 300,350 μm       0.7       6/0       45 cm         OSU05       Double spatula 1/4       2 x 8 mm,300,350 μm       0.7       6/0       45 cm         OSU05       Double spatula 1/4	11401448	Reverse cutting 1/2	48 mm	4	1	90 cm
11202160       Straight cutting       60 mm       2       3/0       90 cm         11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 μm       0.5       7/0       45 cm         OSU08-30       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       45 cm         OSU08       Double spatula 3/8       2 x 6,2 mm, 150 μm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.3       9/0       30 cm         OSU072       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 8 mm, 200 μm       0.5       7/0       20 cm         OSU06       Double spatula 1/4       2 x 8 mm, 300,350 μm       0.7       6/0       45 cm         OSU05       Double spatula 1/4       2 x 8 mm,300,350 μm       0.7       6/0       45 cm         OSU05       Double spatula 3/8       6.2 mm,150 μm       0.2       10/0       10 cm	11501448	Reverse cutting 1/2	48 mm	5	2	90 cm
11302160       Straight cutting       60 mm       3       2/0       75 cm         OSU07       Double spatula 1/2       2 x 8 mm, 200 μm       0.5       7/0       45 cm         OSU08-30       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       30 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       45 cm         OSU08       Double spatula 3/8       2 x 6 mm, 200 μm       0.4       8/0       45 cm         OSU09       Double spatula 3/8       2 x 6.2 mm, 150 μm       0.2       10/0       30 cm         OSU09       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.3       9/0       30 cm         OSU072       Double spatula 3/8       2 x 6.4 mm, 200 μm       0.5       7/0       45 cm         OSU071       Double spatula 3/8       2 x 8 mm, 200 μm       0.5       7/0       20 cm         OSU06       Double spatula 1/4       2 x 8 mm, 300,350 μm       0.7       6/0       45 cm         OSU05       Double spatula 1/4       2 x 8 mm,300,350 μm       1       5/0       45 cm         OSU05       Double spatula 3/8       6.2 mm,150 μm       0.2       10/0       10 cm	11601448	Reverse cutting 1/2	48 mm	6	3&4	90 cm
OSU07         Double spatula 1/2         2 x 8 mm, 200 µm         0.5         7/0         45 cm           OSU08-30         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         30 cm           OSU08         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         45 cm           OSU08         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         45 cm           OSU09         Double spatula 3/8         2 x 6,2 mm, 150 µm         0.2         10/0         30 cm           OSU09         Double spatula 3/8         2 x 6.4 mm, 200 µm         0.3         9/0         30 cm           OSU072         Double spatula 3/8         2 x 6.4 mm, 200 µm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 µm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 µm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300,350 µm         1         5/0         45 cm           OSU101         Spatula 3/8         6.2 mm,150 µm         0.2         10/0         10 cm	11202160	Straight cutting	60 mm	2	3/0	90 cm
OSU08-30         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         30 cm           OSU08         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         45 cm           OSU08         Double spatula 3/8         2 x 6 mm, 200 µm         0.4         8/0         45 cm           OSU09         Double spatula 3/8         2 x 6.2 mm, 150 µm         0.2         10/0         30 cm           OSU09         Double spatula 3/8         2 x 6.4 mm, 200 µm         0.3         9/0         30 cm           OSU072         Double spatula 3/8         2 x 6.4 mm, 200 µm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 µm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 µm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300, 350 µm         1         5/0         45 cm           OSU101         Spatula 3/8         6.2 mm,150 µm         0.2         10/0         10 cm	11302160	Straight cutting	60 mm	3	2/0	75 cm
OSU08         Double spatula 3/8         2 x 6 mm, 200 μm         0.4         8/0         45 cm           OSU10         Double spatula 3/8         2 x 6,2 mm, 150 μm         0.2         10/0         30 cm           OSU09         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.3         9/0         30 cm           OSU072         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU07	Double spatula 1/2	2 x 8 mm, 200 µm	0.5	7/0	45 cm
OSU10         Double spatula 3/8         2 x 6,2 mm, 150 μm         0.2         10/0         30 cm           OSU09         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.3         9/0         30 cm           OSU072         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU08-30	Double spatula 3/8	2 x 6 mm, 200 µm	0.4	8/0	30 cm
OSU09         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.3         9/0         30 cm           OSU072         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         45 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm, 300, 350 μm         1         5/0         45 cm           OSU101         Spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU08	Double spatula 3/8	2 x 6 mm, 200 µm	0.4	8/0	45 cm
OSU072         Double spatula 3/8         2 x 6.4 mm, 200 μm         0.5         7/0         45 cm           OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm, 300,350 μm         1         5/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300,350 μm         1         5/0         45 cm           OSU101         Spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU10	Double spatula 3/8	2 x 6,2 mm, 150 µm	0.2	10/0	30 cm
OSU071         Double spatula 3/8         2 x 8 mm, 200 μm         0.5         7/0         20 cm           OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm, 300,350 μm         1         5/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm,300,350 μm         0.2         10/0         10 cm	OSU09	Double spatula 3/8	2 x 6.4 mm, 200 µm	0.3	9/0	30 cm
OSU06         Double spatula 1/4         2 x 8 mm, 300,350 μm         0.7         6/0         45 cm           OSU05         Double spatula 1/4         2 x 8 mm, 300, 350 μm         1         5/0         45 cm           OSU01         Spatula 3/8         6.2 mm, 150 μm         0.2         10/0         10 cm	OSU072	Double spatula 3/8	2 x 6.4 mm, 200 µm	0.5	7/0	45 cm
OSU05         Double spatula 1/4         2 x 8 mm,300, 350 μm         1         5/0         45 cm           OSU101         Spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU071	Double spatula 3/8	2 x 8 mm, 200 µm	0.5	7/0	20 cm
OSUI01         Spatula 3/8         6.2 mm,150 μm         0.2         10/0         10 cm	OSU06	Double spatula 1/4	2 x 8 mm, 300,350 µm	0.7	6/0	45 cm
	OSU05	Double spatula 1/4	2 x 8 mm,300, 350 µm	1	5/0	45 cm
<b>OSU09-1</b> Spatula 3/8 6.4 mm,200 μm 0.3 9/0 30 cm	OSU101	Spatula 3/8	6.2 mm,150 µm	0.2	10/0	10 cm
	OSU09-1	Spatula 3/8	6.4 mm,200 µm	0.3	9/0	30 cm

# Suture



#### Surgicryl Rapid Suture

#### Features:

- Braided and coated multifilament;absorbable;
- | composition: polyglycolic acid
- | coating: polycaprolactone and calcium stearate (1%)
- colour: beige/ undyed, violet;
- box with 12 sutures;
- The hydrolytic action by which the material is broken down results in total absorption in approximately 42 days.
   Approximately 50% of tensile strength remain after 7 days.

#### Advantages:

- excellent tissue support;
- | minimal tissue reaction;
- | sythetic, coated multifilament;
- | breaks down by hydrolysis
- excellent handling properties;
- high flexibility;
- high tensile strength;
- excellent knot security;
- soft passage through the tissues without sawing, tissue drag and trauma;
- excellent knot tie down;
- | guaranteed seal and product sterility.

#### Indications:

- | superficial soft tissue closing of skin and mucosa;
- | plastic surgery,obstetrics , gynaecology, urology, ophthalmic surgery, pediatric surgery

#### **Contraindications**:

- | cardiovascular surgery, neurosurgery, ligation
- | when wound support for more than 7 days is required.
- caution in patients with delayed wound healing
- (e.g. compromised immune system).

#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
14100113	Round 1/2	13 mm	1	5/0	75 cm
14100117	Round 1/2	17 mm	1	5/0	75 cm
14150117	Round 1/2	17 mm	1.5	4/0	75 cm
14200117	Round 1/2	17 mm	2	3/0	75 cm
14150122	Round 1/2	22 mm	1.5	4/0	75 cm
14200122	Round 1/2	22 mm	2	3/0	75 cm
14300122	Round 1/2	22 mm	3	2/0	75 cm
14200126	Round 1/2	26 mm	2	3/0	75 cm
14300126	Round 1/2	26 mm	3	2/0	75 cm
14350126	Round 1/2	26 mm	3.5	0	75 cm
14200130	Round 1/2	30 mm	2	3/0	75 cm
14300130	Round 1/2	30 mm	3	2/0	75 cm
14350130	Round 1/2	30 mm	3.5	0	75 cm
14200136	Round 1/2	36 mm	2	3/0	75 cm
14300136	Round 1/2	36 mm	3	2/0	75 cm
14350136	Round 1/2	36 mm	3.5	0	90 cm



Sutures without needle:	
-------------------------	--

	1	Thread thickness	I associate
Cat. No.	EP	USP	Lenght
14200015	2	3/0	150 cm
14300015	3	2/0	150 cm
14350015	3.5	0	150 cm
14400015	4	1	150 cm
14500015	5	2	150 cm

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
14400136	Round 1/2	36 mm	4	1	90 cm
4300140	Round 1/2	40 mm	3	2/0	75 cm
4350140	Round 1/2	40 mm	3.5	0	90 cm
4400140	Round 1/2	40 mm	4	1	90 cm
4500140	Round 1/2	40 mm	5	2	90 cm
4300148	Round 1/2	48 mm	3	2/0	75 cm
4350148	Round 1/2	48 mm	3.5	0	90 cm
4400148	Round 1/2	48 mm	4	1	90 cm
4500148	Round 1/2	48 mm	5	2	90 cm
4200517	Round cutting 1/2	17 mm	2	3/0	75 cm
4300526	Round cutting 1/2	26 mm	3	2/0	75 cm
43000536	Round cutting 1/2	26 mm	3	2/0	90 cm
4350536	Round cutting 1/2	36 mm	3.5	0	90 cm
44000536	Round cutting 1/2	36 mm	4	1	90 cm
43000540	Round cutting 1/2	40 mm	3	2/0	90 cm
43505540	Round cutting 1/2	40 mm	3.5	0	90 cm
44000540	Round cutting 1/2	40 mm	4	1	90 cm
4100416	Round cutting 3/8	16 mm	1	5/0	75 cm
4150416	Round cutting 3/8	16 mm	1.5	4/0	75 cm
4200416	Round cutting 3/8	16 mm	2	3/0	75 cm
4150419	Round cutting 3/8	19 mm	1.5	4/0	75 cm
4200419	Round cutting 3/8	19 mm	2	3/0	75 cm
4150317	Round 5/8	17 mm	1.5	4/0	75 cm
4071512	Reverse cutting 3/8	12 mm	0.7	6/0	75 cm
4101512	Reverse cutting 3/8	12 mm	1	5/0	75 cm
4151512	Reverse cutting 3/8	12 mm	1.5	4/0	75 cm
4071512P	Reverse cutting 3/8	12 mm Easy Pass	0.7	6/0	50 cm
4101512P	Reverse cutting 3/8	12 mm Easy Pass	1	5/0	50 cm
4151512P	Reverse cutting 3/8	12 mm Easy Pass	1.5	4/0	50 cm
4101516	Reverse cutting 3/8	16 mm	1	5/0	75 cm
4151516	Reverse cutting 3/8	16 mm	1.5	4/0	75 cm
4201516	Reverse cutting 3/8	16 mm	2	3/0	45 cm
4101516VB	Reverse cutting 3/8	16 mm Black	1	5/0	75 cm
4101519	Reverse cutting 3/8	19 mm	1	5/0	75 cm
4151519	Reverse cutting 3/8	19 mm	1.5	4/0	75 cm
4201519	Reverse cutting 3/8	19 mm	2	3/0	75 cm
4301519	Reverse cutting 3/8	19 mm	3	4/0	75 cm
4151519VB	Reverse cutting 3/8	19 mm czarna	1.5	4/0	75 cm
4201519VB	Reverse cutting 3/8	19 mm czarna	2	3/0	75 cm
4201524	Reverse cutting 3/8	24 mm	2	3/0	75 cm
4301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
4201524P	Reverse cutting 3/8	24 mm Easy Pass	2	3/0	50 cm
4201526	Reverse cutting 3/8	26 mm	2	3/0	75 cm
4301526	Reverse cutting 3/8	26 mm	3	2/0	75 cm
4351526	Reverse cutting 3/8	26 mm	3.5	0	75 cm
4351530	Reverse cutting 3/8	30 mm	3.5	0	75 cm
4301536	Reverse cutting 3/8	36 mm	3	2/0	75 cm
4351536	Reverse cutting 3/8	36 mm	3.5	0	75 cm
4151417	Reverse cutting 1/2	17 mm	1.5	4/0	45 cm
4201417	Reverse cutting 1/2	17 mm	2	3/0	45 cm
4301436	Reverse cutting 1/2	36 mm	3	2/0	75 cm
4351436	Reverse cutting 1/2	36 mm	3.5	0	75 cm
4351440	Reverse cutting 1/2	40 mm	3.5	0	90 cm
4401440	Reverse cutting 1/2	40 mm	4	1	90 cm
14501440	Reverse cutting 1/2	40 mm	5	2	90 cm
	Straight cutting	60 mm	2	3/0	75 cm

#### Suture



#### Surgicryl Suture

#### Features:

- | Braided and coated multifilament;
- absorbable;
- colour: violet;
- Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide; coating Poly-glycolide-co-L-lactide and calcium stearate (<1%)</li>
   box with 12 sutures
- box with 12 sutures.

#### Advantages:

- excellent tissue support;
- | minimal tissue reaction;
- sythetic, coated multifilament;
- breaks down by hydrolysis
- | excellent handling properties;
- high flexibility;
- high tensile strength;
- | excellent knot security;
- soft passage through the tissues without sawing, tissue drag and trauma:
- excellent knot tie down;
- | guaranteed seal and product sterility.

#### In cassettes:

	Thread	thickness	
Cat. No.	EP	USP	- Lenght
17151	1.5	4/0	15 m
17201	2	3/0	15 m
17202	2	3/0	25 m
17301	3	2/0	15 m
17302	3	2/0	25 m
17305	3	2/0	50 m
17351	3.5	0	15 m
17352	3.5	0	25 m
17355	3.5	0	50 m
17401	4	1	15 m
17402	4	1	25 m
17405	4	1	50 m
17501	5	2	15 m
17502	5	2	25 m
17505	5	2	50 m
175010	5	2	100 m
17601	6	3&4	15 m
17602	6	3&4	25 m
17605	6	3&4	50 m
17701	7	5	15 m
17702	7	5	25 m
17705	7	5	50 m
17801	8	6	15 m
17802	8	6	25 m
17805	8	6	50 m

#### Indications:

general soft tissue closing and/or ligation especially in general surgery;
gastrointestinal surgery;
gynaecology;
obstetrics;
urology;
ophthalmic surgery;
plastic surgery;
orthopaedics.

#### Contraindications:

- cardiovascular surgery
  neurosurgery being absorbable, not to use when an extended wound support is required;
  caution in patients with delayed wound healing (e.g.
- compromised immune system), in the region of infected wounds.

#### Sutures without needle:

		Thread thickness	
Cat. No.	EP	USP	— Lenght
15150015	1.5	4/0	150 cm
15200015	2	3/0	150 cm
15300015	3	2/0	150 cm
15350015	3.5	0	150 cm
15400015	4	1	150 cm
15500015	5	2	150 cm
15300045-3	3	2/0	3 x 45 cm
1350045-3	3.5	0	3 x 45 cm
15150045-6	1.5	4/O	6 x 45 cm
15150045-6	2	3/0	6 x 45 cm
15200045-6	3	2/0	6 x 45 cm
15300045-6	3.5	0	6 x 45 cm
15150045	1.5	4/O	12 x 45 cm
15200045	2	3/0	12 x 45 cm
15300045	3	2/0	12 x 45 cm
15350045	3.5	0	12 x 45 cm
15400045	4	1	12 x 45 cm
15500045	5	2	12 x 45 cm
15200070-2	2	3/0	2 x 70 cm
15300070-2	3	2/0	2 x 70 cm
15350070-2	3.5	0	2 x 70 cm
15400070-2	4	1	2 x 70 cm
15500070-2	5	2	2 x 70 cm
15200025	2	3/0	250 cm
15300025	3	2/0	250 cm
15350025	3.5	0	250 cm
15400025	4	1	250 cm
15500025	5	2	250 cm



#### Sutures with needle:

Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
5071412	Reverse cutting 1/2	12 mm	0.7	6/0	45 cm
5071412	Reverse cutting 1/2	17 mm	1	5/0	75 cm
5151417	Reverse cutting 1/2	17 mm	1.5	4/0	75 cm
5201422	Reverse cutting 1/2	22 mm	2	3/0	75 cm
5501423	Reverse cutting 1/2	23 mm	5	2	90 cm
5201426	Reverse cutting 1/2	26 mm	2	3/0	75 cm
5301426	Reverse cutting 1/2	26 mm	3	2/0	75 cm
5351426	Reverse cutting 1/2	26 mm	3.5	0	75 cm
5201430	Reverse cutting 1/2	30 mm	2	3/0	75 cm
5301436	Reverse cutting 1/2	36 mm	3	2/0	75 cm
5351436	Reverse cutting 1/2	36 mm	3.5	0	90 cm
5401436	Reverse cutting 1/2	36 mm	4	1	90 cm
5501436	Reverse cutting 1/2	36 mm	5	2	90 cm
5201440	Reverse cutting 1/2	40 mm	2	3/0	90 cm
5301440	Reverse cutting 1/2	40 mm	3	2/0	90 cm
5351440	Reverse cutting 1/2	40 mm	3.5	0	90 cm
5401440	Reverse cutting 1/2	40 mm	4	1	90 cm
5501440	Reverse cutting 1/2	40 mm	5	2	90 cm
5351448	Reverse cutting 1/2	48 mm	3.5	0	90 cm
5401448	Reverse cutting 1/2	48 mm	4	1	90 cm
5501448	Reverse cutting 1/2	48 mm	5	2	90 cm
5071512	Reverse cutting 3/8	12 mm	0.7	6/0	45 cm
5101512	Reverse cutting 3/8	12 mm	1	5/0	45 cm
5151512	Reverse cutting 3/8	12 mm	1.5	4/0	45 cm
5071516	Reverse cutting 3/8	16 mm	0.7	6/0	45 cm
5101516	Reverse cutting 3/8	16 mm	1	5/0	75 cm
5151516	Reverse cutting 3/8	16 mm	1.5	4/0	75 cm
5201516-45		16 mm	2	3/0	45 cm
5201516UP		16 mm Easy Pass	2	3/0	50 cm
5071516VB	Reverse cutting 3/8	16 mm czarma	0.7	6/0	45 cm
5101519	Reverse cutting 3/8	19 mm	1	5/0	75 cm
5151519	Reverse cutting 3/8	19 mm	1.5	4/0	75 cm
5201519	Reverse cutting 3/8	19 mm	2	3/0	75 cm
5151519UP	Reverse cutting 3/8	19 mm Easy Pass	1.5	4/0	50 cm
51015190P	Reverse cutting 3/8	19 mm czarna	1.5	5/0	75 cm
5151519VB	Reverse cutting 3/8	19 mm czarna	1.5	4/0	75 cm
5151519VB			1.5	4/0	
	Reverse cutting 3/8	24 mm			75 cm
5201524	Reverse cutting 3/8	24 mm	2	3/0	75 cm
5301524	Reverse cutting 3/8	24 mm	3	2/0	75 cm
5351524	Reverse cutting 3/8	24 mm	3.5	0	75 cm
5201524UP		24 mm Easy Pass	2	3/0	50 cm
5201526	Reverse cutting 3/8	26 mm	2	3/0	75 cm

#### Suture



Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
15301526	Reverse cutting 3/8	26 mm	3	2/0	75 cm
15201526UP	Reverse cutting 3/8	26 mm Easy Pass	2	3/0	75 cm
15201530	Reverse cutting 3/8	30 mm	2	3/0	75 cm
15301530	Reverse cutting 3/8	30 mm	3	2/0	75 cm
15351530	Reverse cutting 3/8	30 mm	3.5	0	75 cm
15301536	Reverse cutting 3/8	<u>36 mm</u>	3	2/0	75 cm
15351536	Reverse cutting 3/8	36 mm	3.5	0	75 cm
15401536	Reverse cutting 3/8	<u>36 mm</u>	4	1	75 cm
15351540 15401540	Reverse cutting 3/8 Reverse cutting 3/8	40 mm 40 mm	3.5	0	<u>75 cm</u> 75 cm
15301545	Reverse cutting 3/8	45 mm	3	2/0	75 cm
15701580	Reverse cutting 3/8	80 mm	7	5	150 cm
15150513	Reverse cutting 1/2	13 mm	1.5	4/0	75 cm
15150517	Reverse cutting 1/2	17 mm	1.5	4/0	75 cm
15150517	Reverse cutting 1/2	17 mm	2	3/0	75 cm
15150526	Reverse cutting 1/2	26 mm	1.5	4/0	75 cm
15200526	Reverse cutting 1/2	26 mm	2	3/0	75 cm
15300526	Reverse cutting 1/2	26 mm	3	2/0	75 cm
15300530	Reverse cutting 1/2	30 mm	3	2/0	75 cm
15350530	Reverse cutting 1/2	30 mm	3.5	0	75 cm
15400530	Reverse cutting 1/2	30 mm	4	1	75 cm
15300536	Reverse cutting 1/2	36 mm	3	2/0	90 cm
15350536	Reverse cutting 1/2	36 mm	3.5	0	90 cm
15400536	Reverse cutting 1/2	36 mm	4	1	90 cm
15350540	Reverse cutting 1/2	40 mm	3.5	0	90 cm
15400540	Reverse cutting 1/2	40 mm	4	1	90 cm
15500540	Reverse cutting 1/2	40 mm	5	2	90 cm
15800548	Reverse cutting 1/2	48 mm	8	6	150 cm
15070113	Round 1/2	13 mm	0.7	6/0	75 cm
15100113	Round 1/2	13 mm	1	5/0	75 cm
15150113	Round 1/2	<u>13 mm</u> 17 mm	1.5 1.5	<u>4/0</u> 5/0	75 cm 75 cm
15100117 15150117	Round 1/2 Round 1/2		2	4/0	75 cm
15200117	Round 1/2	<u>17 mm</u> 17 mm	2	3/0	75 cm
15300117	Round 1/2	17 mm	3	2/0	75 cm
15150122	Round 1/2	22 mm	1.5	4/0	75 cm
15200122	Round 1/2	22 mm	2	3/0	75 cm
15300122	Round 1/2	22 mm	3	2/0	75 cm
15400122	Round 1/2	22 mm	4	1	75 cm
15150126	Round 1/2	26 mm	1.5	4/0	75 cm
15200126	Round 1/2	26 mm	2	3/0	75 cm
15300126	Round 1/2	26 mm	3	2/0	75 cm
15350126	Round 1/2	26 mm	3.5	0	75 cm
15400126	Round 1/2	26 mm	4	1	75 cm
15500126	Round 1/2	26 mm	5	2	75 cm
15200130	Round 1/2	30 mm	2	3/0	75 cm
15300130	Round 1/2	30 mm	3	2/0	75 cm
15350130	Round 1/2	<u>30 mm</u>	3.5	0	75 cm
15400130	Round 1/2	30 mm	4	1	75 cm
15500136	Round 1/2	30 mm	5	2	90 cm
15200136	Round 1/2	36 mm	2	<u>3/0</u> 2/0	75 cm
15300136 15350136	Round 1/2 Round 1/2	36 mm 36 mm	3 3.5	0	75 cm 90 cm
15400136	Round 1/2	36 mm	4	1	90 cm
15500136	Round 1/2	36 mm	5	2	90 cm
15200140	Round 1/2	40 mm	2	3/0	75 cm
15300140	Round 1/2	40 mm	3	2/0	75 cm
15350140-75		40 mm	3.5	0	75 cm
15350140	Round 1/2	40 mm	3.5	0	90 cm
15400140-75	-	40 mm	4	1	75 cm
15400140	Round 1/2	40 mm	4	1	90 cm
		40 mm	5	2	75 cm
15500140-75	Round I/2	40 11111			
15500140-75 15500140	Round 1/2 Round 1/2	40 mm	5	2	90 cm
			5 5	2 2	90 cm 90 cm
15500140	Round 1/2	40 mm			
15500140 15500145	Round 1/2 Round 1/2	40 mm 45 mm	5	2	90 cm
15500140 15500145 15300148 15350148 1540148	Round 1/2 Round 1/2 Round 1/2 Round 1/2 Round 1/2	40 mm 45 mm 48 mm 48 mm 48 mm	5 3 3.5 4	2 2/0 0 1	90 cm 90 cm 90 cm 90 cm
15500140 15500145 15300148 15350148	Round 1/2 Round 1/2 Round 1/2 Round 1/2	40 mm 45 mm 48 mm 48 mm	5 3 3.5	2 2/0 0	90 cm 90 cm 90 cm


Cat. No.	Needle shape	Needle size	Metric	USP	Lenght
15400160	Round 1/2	60 mm	4	1	75 cm
15300165	Round 1/2	65 mm	3	2/0	90 cm
15350165	Round 1/2	65 mm	3.5	0	90 cm
15400180	Round 1/2	80 mm	4	1	100 cm
15200126-2	Round double 1/2	2 x 26 mm	2	3/0	75 cm
15250945	Round ½ loop	45 mm	3.5	0	150 cm
15300327	Round 5/8	27 mm	3	2/0	75 cm
15350327	Round 5/8	27 mm	3.5	0	75 cm
15150217	Round 3/8	17 mm	1.5	4/0	45 cm
15200222	Round 3/8	22 mm	2	3/0	75 cm
15400230	Round 3/8	30 mm	4	1	75 cm
15202160	Cutting straight	60 mm	2	3/0	75 cm
15302160	Cutting straight	60 mm	3	2/0	75 cm
15301830	Blunt 1/2	30 mm	3	2/0	75 cm
15351836	Blunt 1/2	36 mm	3.5	0	90 cm
15401836	Blunt 1/2	36 mm	4	1	90 cm
15501845	Blunt 1/2	45 mm	5	2	90 cm
15401963	Blunt 3/8	63 mm	4	1	100 cm
15501985	Blunt 3/8	85 mm	5	2	75 cm
M91087	Micro round		0.4	8/0	20 cm
0910081	Spatula double 1/2	6,4 mm 2 x 8 mm, 200 μm	0.5	7/0	45 cm
0910081	Spatula double 1/2		0.7	6/0	45 cm
091007	Spatula double 1/2 Spatula double 1/2	2 x 8 mm, 200, 300 μm 2 x 8 mm, 200, 300 μm	0.7	6/0	45 cm
O910061U			0.7	6/0	45 cm
O9100810	Spatula double 1/2	2 x 8 mm, 300 µm	0,7	8/0 8/0	45 cm
O9100B	Spatula double 3/8	2 x 6 mm, 200,220 µm	0.5	7/0	
O910072 O910051	Spatula double 3/8	2 x 6 mm, 200,220 µm	1	5/0	45 cm
0910051	Spatula double 3/8	2 x 11 mm, 350 µm	0.7	5/0 6/0	45 cm
	Spatula double 3/8	2 x 8 mm, 300 µm	0.7	6/0 6/0	45 cm
	Spatula double 1/4	2 x 8 mm, 300 µm	0.7	6/0 6/0	20 cm
091006U	Spatula double 1/4	2 x 8 mm, 300 µm	0.7	5/0	45 cm
091005	Spatula double 1/4	2 x 8 mm, 350 µm	1		45 cm
O91005U	Spatula double 1/4	2 x 8 mm, 350 µm	1	5/0	45 cm

# Suture



#### HIGH-QUALITY PACKAGING – FULL PROFESSIONALISM IN SURGERY

SAFE AND HANDY CASETTE



#### | made of the highest quality stainless steel;

| 1/2 wheel, triangular, with one eyelet.

Suture needle spring Post-MORTEM

Cat. No.	Name
147692	Section needle 1/2 circle, triangular, 115 mm 12 pcs.
147690	Section needle 1/2 triangular circles, 71 mm 12 pcs.
147691	Section needle 1/2 circle, triangular, 97 mm, 12 pcs.

#### Standard suture needles

- | made of the highest quality stainless steel;
- extremely sharp and durable;
- a smoothly polished surface minimizes the risk of inflammation.

Cat. No.	Name	Shape
157708	Surgical needle 3/8 circle, 103 mm 12 pcs.	Round
157746	Surgical needle 3/8 circle, 13 mm 12 pcs.	Round
157745	Surgical needle 3/8 circle, 15 mm 12 pcs.	Round
157744	Surgical needle 3/8 circle, 24 mm 12 pcs.	Round
157743	Surgical needle 3/8 circle, 32 mm 12 pcs.	Round
157713	Surgical needle 3/8 circle, 32 mm 12 pcs.	Round
157742	Surgical needle 3/8 circle, 36 mm 12 pcs.	Round
157712	Surgical needle 3/8 circle, 36 mm, one eye 12 pcs.	Round
157741	Surgical needle 3/8 circle, 54 mm 12 pcs.	Round
157711	Surgical needle 3/8 circle, 54 mm, one eye 12 pcs.	Round
157740	Surgical needle 3/8 circle, 70 mm 12 pcs.	Round
157710	Surgical needle 3/8 circle, 70 mm one mesh 12 pcs.	Round
157698	Surgical needle 3/8 circle, 100 mm one eye 12 pcs.	Round
157727	Surgical needle 3/8 circle, 13 mm 12 pcs.	Triangular
157726	Surgical needle 3/8 circle, 15 mm 12 pcs.	Triangular
157725	Surgical needle 3/8 circle, 24 mm 12 pcs.	Triangular
157724	Surgical needle 3/8 circle, 32 mm 12 pcs.	Triangular
157704	Surgical needle 3/8 circle, 35 mm one eye 12 pcs.	Triangular
157723	Surgical needle 3/8 circle, 36 mm 12 pcs.	Triangular
157722	Surgical needle 3/8 circle, 54 mm 12 pcs.	Triangular
157701	Surgical needle 3/8 circle, 54 mm, one eye 12 pcs.	Triangular
157728	Surgical needle 3/8 circle, 63 mm 12 pcs.	Triangular
157702	Surgical needle 3/8 circle, 70 mm one eye 12 pcs.	Triangular
157703	Surgical needle 3/8 circle, 72 mm, one eye 12 pcs.	Triangular
157721	Surgical needle 3/8 circle, 73 mm 12 pcs.	Triangular
157720	Surgical needle 3/8 circle, 80 mm 12 pcs.	Triangular
157700	Surgical needle 3/8 circle, 80 mm one mesh 12 pcs.	Triangular



In order to ensure the highest quality, thread production is carried out in accordance with the requirements of ISO 9001: 2008; ISO 13485: 2003 and in accordance with European (EDQM) and American Pharmacopoeia standards.





#### SUTURE NEEDLE





#### STAPLER / SKIN STAPLE REMOVER

#### POLSTAP - stapler disposable

Simple and easy to operate.

- High sewing up and stapling speed help shorten the time of sewing and avoid cross infection of contagious diseases.
- Reduced tissue reaction, thinner scar after operation and faster wound healing.Painless and convenient staples removal after the incision has healed up.

Cat. No. **340755** 



#### Disposable skin staple remover

Original forceps for removing staples.Sterile.

Cat. No. 340756





# DIATERMO MB 122 electrocauter

# Mathieu Olsen Premium needle forceps





MORE INFORMATION ON PAGE 140



MORE INFORMATION ON PAGE 136

# SCISSORS

07

#### Standard straight scissors

| made of stainless steel;

- standard hardness;
- straight model.

Cat. No.			
Type sharp/sharp	Type sharp/blunt	Type blunt/blunt	Name/Size
130445S	130446S	130447S	Standard straight scissors/14 cm
130455S	130456S	130457S	Standard straight scissors/16 cm
130465S	130466S	130467S	Standard straight scissors/18 cm



#### Premium straight scissors

made of hardened steel;reinforced cutting parts;precisely made..

Cat. No.			
Type sharp/sharp	Type sharp/blunt	Type blunt/blunt	Name/Size
130445P	130446P	130447P	Premium straight scissors/14 cm
130455P	130456P	130457P	Premium straight scissors/16 cm
130465P	130466P	130467P	Premium straight scissors/18 cm


#### Standard bent scissors

| made of stainless steel;

- standard hardness;
- | bent model (Cooper).

Cat. No.			
Type sharp/sharp	Type sharp/blunt	Type blunt/blunt	Name/Size
130115S	130116S	130117S	Premium bent scissors/14 cm
130215S	130216S	130217S	Premium bent scissors/16 cm
130315S	130316S	130317S	Premium bent scissors/18 cm





#### Premium bent scissors

- | made of hardened steel;
- | reinforced cutting parts;
- precisely made;
- | bent model (Cooper).

Cat. No.			
Type sharp/sharp	Type sharp/blunt	Type blunt/blunt	Name/Size
130115P	130116P	130117P	Premium bent scissors/14 cm
130215P	130216P	130217P	Premium bent scissors/16 cm
130315P	130316P	130317P	Premium bent scissors/18 cm

#### Scissors for cutting bandage Lister Premium type

- | made of hardened steel;
- | reinforced cutting parts;
- intended for cutting various kinds of dressings, bandages, and bandaids.

Cat. No.	Name/Size
34125	Scissors for cutting bandage Lister type/11 cm
34129	Scissors for cutting bandage Lister type/15 cm
34130	Scissors for cutting bandage Lister type/18 cm

#### Scissors for cutting bandage Lister Premium type

- | made of hardened steel;
- | cutting part reinforced with wolfram;
- highly reliable for hardness and corrosion resistance;intended for cutting various kinds of dressings, bandages,
- and bandaids.

Name/Size
Scissors for cutting bandage Lister Premium type/14 cm
Scissors for cutting bandage Lister Premium type/18 cm

#### Surgery



SCISSORS







**07** S

# Surgery

#### SCISSORS

#### Stitch removal scissors

13 cm long;Spencer model.

Cat. No. 26737



#### Standard Mayo scissors

made of stainless steel;both ends rounded.

Cat. No.		– Name/rozmiar
Proste	Zagięte	Name/IO2IIIIai
26846	26849	Standard Mayo scissors/14,5 cm
26847	26850	Standard Mayo scissors/18 cm
26848	26851	Standard Mayo scissors/20 cm



#### Mayo Premium scissors

| made of hardened steel;

- | cutting part reinforced with wolfram;
- | both ends rounded.

Cat. No.		
Proste	Zagięte	— Name/rozmiar
26520	26522	Mayo Premium scissors/14,5 cm
26521	26523	Mayo Premium scissors/18 cm



#### Iris scissors Premium

| made of hardened steel; | cutting part reinforced with wolfram; | highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26526	Iris scissors Premium- straight/11.5 cm
26527	Iris scissors Premium- curved/11.5 cm

#### Goldman Fox scissors Premium

| made of hardened steel; | cutting part reinforced with wolfram; | highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26528	Goldman Fox scissors Premium- straight/13 cm
26529	Goldman Fox scissors Premium - curved/13 cm

#### Universal scissors Premium

| made of hardened steel; | cutting part reinforced with wolfram; | highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26506	Universal scissors/12 cm

# Surgery



SCISSORS







#### SCISSORS

#### Metzenbaum scissors Premium

| made of hardened steel;

| highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26510	Metzenbaum scissors Premium straight/14 cm
26511	Metzenbaum scissors Premium straight/18 cm
26512	Metzenbaum scissors Premium straight/20 cm
26514	Metzenbaum scissors Premium curved/14 cm
26515	Metzenbaum scissors Premium curved/18 cm
26516	Metzenbaum scissors Premium curved/20 cm



#### Lister scissors Premium

| made of hardened steel;

| highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
38530	Lister scissors Premium/14 cm
38532	Lister scissors Premium/18 cm



#### Kelly scissors Premium

| made of hardened steel;

| highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
38520	Kelly scissors Premium – straight sharp/16 cm
38521	Lister scissors Premium - curved sharp /16 cm

#### Debakey Potts Smith scissors Premium

| made of hardened steel;

| highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26507	Debakey Potts Smith scissors19 cm - angle 25°
26508	Debakey Potts Smith scissors 19 cm - angle 45 $^\circ$
26509	Debakey Potts Smith scissors 19 cm - angle 60 $^\circ$
26517	Debakey Potts Smith scissors 23 cm - angle25°
26518	Debakey Potts Smith scissors 23 cm – angle $45^{\circ}$
26519	Debakey Potts Smith scissors 23 cm - angle $60^\circ$

# Surgery



#### SCISSORS





#### TWEEZERS

#### Adson Forceps

made of hardened steel;highly reliable for hardness and corrosion resistance.

#### Cat. No. Name/Size

26560 Adson froceps Premium - 12 cm



#### Potts Smith Dissecting Forceps Premium

| made of hardened steel;

| highly reliable for hardness and corrosion resistance.

Cat. No.	Name/Size
26565	Potts Smith Dissecting Forceps - 15 cm
26566	Potts Smith Dissecting Forceps - 18 cm
26567	Potts Smith Dissecting Forceps - 20 cm
26568	Potts Smith Dissecting Forceps - 23 cm
26569	Potts Smith Dissecting Forceps - 25 cm



#### Anatomical tweezers

made of stainless steel;ribbed grabbing part.

#### Cat. No. 140710

Cat. No.	Name/Size
140713A	Anatomical tweezer 13 cm
140714A	Anatomical tweezer 14 cm
140715A	Anatomical tweezer 15 cm
140716A	Anatomical tweezer 16 cm
140718A	Anatomical tweezer 18 cm
140726	Anatomical tweezer 25 cm
140727	Anatomical tweezer 30 cm
140730	Anatomical tweezer 30 cm curved

#### Surgical tweezers

| made of stainless steel;

| grabbing part with a characteristic thickening.

Nr kat.	Nazwa/rozmiar
140713CH	Surgical tweezer 13 cm
140714CH	Surgical tweezer 14 cm
140715CH	Surgical tweezer 15 cm
140716CH	Surgical tweezer 16 cm
140718CH	Surgical tweezer 18 cm
140720CH	Surgical tweezer 20 cm

#### Rochester pean clamp

| made of stainless steel;

6 types;

| used for clamping vessels (haemostatic action).

Cat. No.		
Straiht		— Name/Size
140210	140120	Rochester pean clamp/14 cm
140220	140130	Rochester pean clamp/16 cm
140230	140140	Rochester pean clamp/18 cm

#### Surgery



#### TWEEZERS / ROCHESTER PEAN CLAMP







#### Mathieu Standard needle forceps

made of stainless steel;used for holding the needle during surgery.

#### Cat. No. Name/Size

151721SMathieu Standard needle forceps/14 cm151724SMathieu Standard needle forceps/17 cm



#### Mathieu Premium needle forceps

made of hardened steel;needle-grabbing part reinforced with wolfram.

Cat. No.	Name/Size
151721P	Mathieu Premium needle forceps/14 cm
151724P	Mathieu Premium needle forceps/17 cm



#### Mathieu Olsen Premium needle forceps

| made of hardened steel;

- | needle-grabbing part reinforced with wolfram;
- | ntegrated cutting part.

#### Cat. No. Name

151760Mathieu Olsen Premium needle forceps/14 cm151780Mathieu Olsen Premium needle forceps17 cm



#### Kocher intestinal clamp

| made of stainless steel;

| clamp with a special thickened tip for easy grabbing.

Cat. No.		- Norse (ci
Straight		– Name/Size
151700	151703	Kocher intestinal clamp/14 cm
151701	151704	Kocher intestinal clamp/16 cm
151702	151705	Kocher intestinal clamp/18 cm

#### Backhaus clamp

| made of stainless steel;

| clamp used for securing surgical drapes/cloths.

Cat. No.	Name/Size
140580	Backhaus clamp/9 cm
140590	Backhaus clamp/13 cm
140585	Backhaus clamp/10 cm

#### Kilian/Hartmann nasal clamps

| made of stainless steel;

| 14 cm long;

used in nasoscopy.

Cat. No.	Name/Length of endoscopic part
26792	Kilian/Hartmann nasal clamps/15 cm
26778	Kilian/Hartmann nasal clamps/22 cm
26772	Kilian/Hartmann nasal clamps/35 cm
26776	Kilian/Hartmann nasal clamps/50 cm
26764	Kilian/Hartmann nasal clamps/75 cm

# Surgery



#### NEEDLE FORCEPS / CLAMPS







#### CLAMPS

07

#### | made of stainless steel.

Cat. No.	Name/Type/Size
26714	Mosquito forceps – straight/12.5 cm
26872	Mosquito forceps - straight/14 cm
26715	Mosquito forceps - straight/16 cm
26716	Mosquito forceps- curved/12.5 cm
26873	Mosquito forceps – curved/14 cm



#### | made of stainless steel.

Cat. No.	Name/Type/Size
26757	Kelly forceps – straight/14 cm
26717	Kelly forceps - straight/16 cm
26758	Kelly forceps - straight/14 cm
26718	Kelly forceps- curved/16 cm



#### | made of stainless steel.

Cat. No.	Name/Type/Size
26749	Klemmer forceps- straight/14 cm
26709	Klemmer forceps- straight/16 cm
26835	Klemmer forceps – straight/18 cm
26713	Klemmer forceps – straight/20 cm
26836	Klemmer forceps – curved/14 cm
26837	Klemmer forceps – curved/16 cm
26838	Klemmer forceps – curved/18 cm
26839	Klemmer forceps – curved/20 cm

| made of stainless steel.

Cat. No.	Name/Type/Size
26978	Mixter Forceps/ 23 cm





#### Weitlaner 13 cm wound dilator

- | made of stainless steel;
- | used for dilating wounds;
- | 3 to 4 teeth on each arm.

Cat. No.	Name/Type
26858	Weitlaner 13 cm wound dilator 13 cm/Sharp teeth
26859	Weitlaner 13 cm wound dilator 13 cm/Blunt teeth

- | made of stainless steel;
- compatible with all standard blades;
- | size 3 suitable for blades up to size 21;
- | size 4 suitable for blades size 22, 23, 24.

Cat. No.	Name/Size	
141109	Scalpel handle3	
141110	Scalpel handle/4	

super sharp;

| top workmanship.

Cat. No.	Name/Size
141141	Swann Morton scalpel blade nr 22
141151	Swann Morton scalpel blade nr 23
141182	Swann Morton scalpel blade nr 24
141114	Swann Morton scalpel blade nr 10 /100 pcs/
141115	Swann Morton scalpel blade nr 11 /100 szt/
141116A	Swann Morton scalpel blade nr 12 /100 szt/
141117	Swann Morton scalpel blade nr 15 /100 pcs/
141119	Swann Morton scalpel blade nr 21 /100 szt/
141119A	Swann Morton scalpel blade nr 22 A /100 pcs/
141120	Swann Morton scalpel blade nr 20 /100 pcs/
141128	Swann Morton scalpel blade nr 9 /100 pcs/
141129	Swann Morton scalpel blade nr 6 /100 pcs/
141130A	Swann Morton scalpel blade nr 26 /100 pcs/
141132A	Swann Morton scalpel blade nr 22 /100 szt/
141141A	Swann Morton scalpel blade nr 23 /100 pcs/
141142A	Swann Morton scalpel blade nr 15C /100 pcs/
141182A	Swann Morton scalpel blade nr 24 /100 pcs/

Surgery



#### WOUND DILATORS / SCALPEL ACCESSORIES







#### DIATERMO MB 122 electrocauter

#### Cat. No. 30540

#### DIATERMO MB 160 electrocauter

Cat. No. 30541

- DIATERMO MB 122 and MB 160 units are designed for performing mono- and bipolar surgical procedures.
   they cut and coagulate;
- the strong MB 120 is dedicated to gynaecology and orthopaedics;
- | with an OC safety alarm informing about circuit breakage;
- with OVIT safety system alarming in case of excessive energy supply;
- standard version with: disposable 2-button handle, a pedal used as an alternative to the handle, a neutral plate with 2 m cable, a set of short autoclavable 2.4 mm electrodes.

Parameter	MB 122	MB 160
Max cutting power (CUT)	120 W	160 W
Max cauterising power (BLEND)	90 W	120 W
Max forced coagulation power (FORCED COAG)	80 W	100 W
Max soft coagulation power (SOFT COAG )	60 W	80 W
Max bipolar operation power (BIPOLAR)	40 W	60 W
Working frequency	600 kHz	600 kHz
Weight	7 kg	7 kg
Dimensions [mm]	260 x 110 x 265	260 x 110 x 265
Voltage	115-230	50-60
230 V fuses	2 x 3.15 A	2 x 3.15 A

#### Adhesive cloth spray

spray used for surgical cloths;ensures complete cloth adhesion.

Cat. No. 141865





#### Electrocauter 50 W / 80 W

#### Cat. No. 30625 / 30626

- | monopolar units for small surgical procedures;
- 50 W apparatus used in small procedures of cutting and coagulation;
- 80 W apparatus used for cutting and cauterisation as well as coagulation and microcoagulation;
- standard equipment of a disposable handle, pedal, plate, and cable, 6 autoclavable electrodes (diameter 2.4 mm, length 5 cm);
- 2 m electrical wire.

Parameter	50	80
Max cutting force (CUT)	50 W	80 W
Max power (CUT/COAG1)	45 W	-
Max power (CUT/COAG2)	40 W	-
Max cauterising power BLEND	-	60 W
Max coagulation power	40 W	50 W
Max microcoagulation power MICRO	-	20 W
Working frequency	600 kHz	700 kHz
Weight	2,5 kg	2,5 kg
Dimensions [mm]	190 x 85 x 239	190 x 85 x 239
Voltage	50-60	115-230
230 V fuses	2 x 3.15 A	2 x 3.15 A

#### Electrocauter replacement parts

Cat. No.	Name
30566	Metal plate without cable
30573	20 x 15 cm plate without cable
30563	5 m plate cable
30559	2.5 m plate cable
30686	Cutting end 13 cm x 10 mm x 8 mm
30687	Cutting end 13 cm x 15 mm x 8 mm
30688	Cutting end 13 cm x 2 0mm x 10 mm
30547	ORION monopolar handle
30548	STAR monopolar handle
30667	Straight bipolar clamp 15 cm x 6" x 0,5 mm
30654	Straight bipolar clamp 18 cm x 7" x 1 mm
30652	Bent bipolar clamp 18 cm x 7" x 1 mm
30651	Straight bipolar clamp 20 cm x 7 3/4" x 1 mm
30643	122 clamp cable, 160 MB

#### Surgery



#### ELECTROCAUTERS / ELECTROCAUTER ACCESSORIES





#### SURGICAL DRAPES/FIELDS / LININGS

fabric (SMS); | with or without opening; sterile;

disposable.

Cat. No.	Name
141720	Surgical drape without opening 30 x 45 cm
141721	Surgical drape without opening 60 x 90 cm
141722	Surgical drape without opening 90 x 120 cm
141730	Surgical drape with opening 30 x 45 cm
141731	Surgical drape with opening 60 x 90 cm
141732	Surgical drape with opening 90 x 120 cm





#### Operating field

green;

sterile;

- disposable;
- | single sheet in sterile packaging;
- polyethylene.

Cat. No.	Name
141765	Operating field 35 x 40 cm
141770	Operating field 60 x 90 cm
141780	Operating field 90 x 120 cm
141840	Operating field 120 cm x 120 cm

| for clean and comfortable surgery;

- | available in 10 pc packages;
- disposable;
- super-absorbent.

141871_V	Absorbent lining VETPLUS 45 x 60 cm
141872_V	Absorbent lining VETPLUS 60 x 60 cm
141873 V	Absorbent lining VETPLUS 60 x 90 cm





# Care and medical products

**SNAP-IN COLLAR** 

# Oster GOLDEN A5 45W/220V clipper







MORE INFORMATION ON PAGE 149

## Care and medical products

80

## MORE INFORMATION ON PAGE 146

#### COLLARS

#### Standard snap-in collar

ribbed edges;

- | increased strength;
- streamlined shape;
- | 3–5 regulation grades
- | ribbed surface ensures elasticity and resistance against scratches and folding;
- | made of an ecological material.

Cat. No.	Name
273480S	Standard snap-in collar 7,5 cm
273481S	Standard snap-in collar 10,0 cm
273482S	Standard snap-in collar 12,5 cm
273483S	Standard snap-in collar 15,0 cm
273484S	Standard snap-in collar 20,0 cm
273485S	Standard snap-in collar 25,0 cm
273486S	Standard snap-in collar 30,0 cm



#### Super Premium snap-in collar

Soft and elastic collars with a special rubber edge. Secures the environment against damage, ensures greater comfort to the owners. Transparent, snap-in closing system: 100% polypropylene. Collective box contains 5 pcs.

Cat. No.	Name
273480C	Super Premium snap-in collar 7,5 cm
273481C	Super Premium snap-in collar 10 cm
273482C	Super Premium snap-in collar 12,5 cm
273483C	Super Premium snap-in collar 15 cm
273484C	Super Premium snap-in collar 20 cm
273485C	Super Premium snap-in collar 25 cm
273486C	Super Premium snap-in collar 30 cm



#### Claw pincers

dog and cat pincers;

| made of stainless steel.

Cat. No.	
131300	Claw pincers 12,5cm DOGGYMAN
131305	Claw pincers 16,5cm DOGGYMAN
131415	Claw pincers for cats 8,5 cm

#### Catheter for dogs

disposable;

Cat. No.		Size
273400	Catheter for dogs	1,3 x 500 mm
273405	Catheter for dogs	1,8 x 500 mm
273410	Catheter for dogs	2,0 x 500 mm
273420	Catheter for dogs	2,6 x 500 mm
273415	Catheter for dogs	3,3 x 500 mm

#### Catheter for cats

disposable;

Cat. No.	Name	Size
273451	Catheter for cats with stylet	1,0 x 130 mm
273453	Catheter for cats with stylet	1,0 x 130 mm
273461	Catheter for cats with stylet	1,3 x 130 mm
273463	Catheter for cats with stylet	1,3 x 130 mm

## Care and medical products



#### GROOMING SETS / CATHETERS







#### CATHETERS / WEIGHT SCALES

#### Electronic platform vet scale

- | weight range up to 300 kg;
- dimensions: 600 mm x 800 mm;
- | measurement with an accuracy of 100 g;
- professional;
- | overload and low battery function;
- digital LED display..

Cat. No. 80211



#### Digital weight scale, for small animals

| minimum measurable weight : 2 kg | weight range up to 20 kg; | measurement with an accuracy of 10 g; electronic; | LCD display; automatic turn-off switch; | lithium battery attached. Cat. No. 27261

#### Digital weight - large

- | weight range up to 100 kg;
- dimensions: 62 x 83 cm;
- | measurement with an accuracy of 100 g;
- electronic;
- LCD display;
- | automatic turn-off switch;
- | lithium battery attached.

Cat. No. 27259





#### Digital weight - medium

- | weight range up to 100 kg;
- | dimensions: 35.5 x 63.5 cm;
- | measurement with an accuracy of 100 g;
- electronic;
- LCD display;
- automatic turn-off switch;
- | lithium battery attached.

Cat. No. 27258

#### Oster GOLDEN A5 Limited Edition 2


#### Clipper Cryogen Blade

- Powerful two-speed motor. Cooler-running motor operates 20 degrees cooler than before.
- Over 2700 strokes per minute. Gear driven Universal motor has speeds for more efficient clipping.
- Goes through thick, matted hair easier because it's faster than the single speed clipper. Easy to handle.
- Normal speed for general grooming and precise clipping and High speed for 1/3 faster - quick body clipping. Virtually unbreakable housing.
- Detachable blade system for quick, easy blade changes and increased versatility. Smooth, quiet, reliable and cool-running. Powerful gear driven universal motor uses very high rpm & unique gear down assembly to give this all-around clipper great power for its size. Unitary cutting assembly for permanent parts alignment. Perfect for clipping all animals.

Cat. No. 273145

## Care and medical products

#### WEIGHT SCALES





#### Care and medical products 80

#### CLIPPERS

#### Powerful Professional Clipper Oster GOLDEN A6 Comfort

- | Three Speed 3100/3600/4100spm
- Quiet, heavy duty DC motor designed for continuous load
- Cooler running experience
- Extra slim casing with excellent weight distribution
- Securely sealed to keep fur out
- Zinc alloy technology
- | Compatible with all professional A5 detachable style blades
- | High-efficiency, cooler running motor
- | No greasing of internal parts
- 3.66 metre cord, designed to minimize outlet space
- Only weighs 392g
- | Includes size 1 year limited warranty

Cat. No. BRAK KODU



#### Clipper Moser 1245 Max 45

- professional and reliable;
- | minimum noise and vibration;
- Has powerful rotary motor 45 Watts for maximum power;
- | Can operate at two sped;
- | Ergonomically designed;
- Weight: 450 gram
- Cable lenght 3 m

Cat. No. 273177

#### Oster<sup>®</sup> Professional blades

| Oster<sup>®</sup> Professional blades are produced from the finest high-carbon steel and handcrafted to our rigorous standards, providing unparalleled precision and performance. Available in multiple sizes for various cutting lengths and maximum versatility.

Cat. No. 273150



#### Heiniger Saphir Cordless Dog Clipper

- | This is proving to be one of our best selling clippers and is manufactured by Heiniger in Switzerland who are well renowned manufacturers of horse and sheep clippers. The quality is without doubt excellent and the power is truly amazing for a cordless pet clipper. It has a slim grip hand piece and stays cool even when clipping for long periods. Takes all A5 type snap on blades. It will clip for 50-60 minutes on each battery with a recharge time of 45 mins.
- Comes with 2 batteries and a charger which will allow continuous clipping all day.
- Weighs only 410 grams and measures 204 x 50 x 41 mm Clipper sold with full 12 months warranty and comes in a very attractive carry case.

Cat. No. 273144

## Care and medical products



#### CLIPPERS





80

## Care and medical products

#### LASERS

#### MR4 Acti Vet lase

| portable;

- great power (25000 mW);power greater that in most class IV lasers;universal;
- . .
- non-invasive;
- 6 treatment programs;
- | combines blue and red light;
- easy operation;
- | up to 8 hours of power supply;
- | perfect for sports events;
- device properties:
- 470 nm blue light;
- 640 nm red light;
- 870 nm infrared;
- 950 nm pulsating laser;

| treating chronic pain;

- | improved elasticity and endurance;
- accelerated healing;
- | FDA-certified;
- | blue light therapy helps control MRSA (Staphylococcus aureus);
- treats for example injuries of the talocrural joint, the shin, elbows, contortions, pain of the back and the Achilles tendon;
- access with acupuncture probe that reaches the target nervous-muscular point and increases laser efficiency.

Cat. No. 273200

#### TQ SOLO laser

| portable and light;

- easy to use (locate the point of pain, turn on the laser and set the suitable frequency in accordance with the manual, press the 'start' button);
- | 15000 mW;
- | minimises chronic and acute pain;
- 4 types of radiation (red light, infrared, pulsating laser, static magnetic field);
- laser: 905 nm;
- infrared: 875 nm;
- red light: 660 nm;
- | lithium battery for 18 h;
- 3 preset frequencies;
- | manual.

Cat. No. 273201



'The technology of the MR4 laser is a rescue for me and my sportsmen. The knee ligament injury was chronic. After two days the scar was smaller by 50%, and the pain was smaller by 90%'

#### Tom Barbeau

NHL hockey player, coach of the national representation, headmaster of Waterville Valley skiing academy



# Diagnostic imaging

# Zoncare V3 vet Veterianry Color Doppler Ultrasound

Careray CareView Detectors





MORE INFORMATION ON PAGE 173

MORE INFORMATION ON PAGE 160

#### PROTECTION

#### Rubber and surgical X-ray gloves

made from lead vinyl, seamless;
universal size, black;
available in equivanet of: 0.25 or 0.5 mm Pb;
surgical version in blue colour.

Cat. No. 27100



#### Protective X-ray coat

| one-sided, with two protection levels:

- light: 0.25 mm Pb on entire area,
- heavy: 0.5 mm Pb on the front and 0.25 on other parts,
- all one-sided coats protect the front, sides, and shoulder blades of the user, and come in three types: 110, 120, and 130 cm;
- colours: blue and 'leopard'.

Cat. No.	Name
283023	110 cm x 0.5 Pb protective X-ray coat
283024	120 cm x 0.5 Pb protective X-ray coat
283025	130 cm x 0.5 Pb protective X-ray coat
283026	110 cm x 0.25 Pb protective X-ray coat
283027	120 cm x 0.25 Pb protective X-ray coat
283028	130 cm x 0.25 Pb protective X-ray coat

#### Thyroid cover

#### Cat. No. 283029





#### X-Ray Positioning Aid

- | suitable positioning aid for dogs and cats
- | adjustable to all sizes required
- | wood construction

Cat. No. 703510

#### Berlin type positioner

durable construction;

- | transparent plexiglass;
- | the shape of the letter V;
- ideal also as a positioner during chest surgery; different sizes: 5 sizes.

at. No. 703550

#### RTG positioner for birds

comfortable for all sizes of birds; they can also be used as a restraint during operation; contains cervical plates (12/21/30/35 mm).

at. No. 703506

#### X-Ray Positioning Aid – Navicular Stanc

- outstanding elevation of the navicular bone is possible for a 60° oblique view
- | includes a bracket for cassettes
- | plexiglas construction

Cat. No. 703540

POSITIONING









#### SANDBAG

#### Pawpositioner

does not absorb x-rays;

- made of unique blend of high density polyofelins in USA;
   protection during operations or diagnostic imaging;
- durable;
- easy to clean;
- portable;
- can be used in CT and computed tomography, magnetic resonance imaging and myelography;
- 5 sizes;
- | able to contain animals weighing between 2 and 60 kg.

#### Cat. No. 605920-605925



#### Sandbag

- position aid for patients requiring x-rays or surgeryred, durable vinyl cover
- | malleable
- easy to clean and disinfect
- \* Not safe for washing machines

Cat. No. 704515



#### Tubular Sandbag

- positioning aid for animals during X-Ray or surgical procedures
- | cylinder
- | durable, wipe clean plastic cover

Cat. No. 704550



#### Tubular Sandbag

- | for animals during radiographs and surgery
- | durable, easy to clean plastic
- | 1 kg of sand
- colour: Ruby Red

Cat. No. 704551

#### Foam Positioning Aid

| foam positioning aid for X-Ray or surgical procedures.

Cat. No. 605900

#### Digital Portable X-ray System EPX-Seriese

| EPX-F2400 (100kV/40mA)2.4KW

- innovation designed Ultra-Light, High frequency portable x-ray.
- | high Performance EPX-Series Portable X-ray system.
- | most compact size and lightest weight.
- | provide high resolution imaging at lower dosages.
- durable design with solid body.
- | soft Touch Controls with Digital displays.
- | calibrated cassette size indicator dials
- available with Inverted control panel for table use.
- two stage dynamic, Auto line compensation.
- | easy access to DR systems.
- | dual Integrated Laser Pointer.

## Diagnostic imaging

## FOAM POSITIONING / GENERATOR








#### CARERAY CAREVIEW DETECTORS


#### Careray CareView Detectors


- multi-purpose detector CareView is designed to meet a wide-range of requirements of digital radiography applications, e.g. retrofit of existing radiography systems (e.g. CCD-DR, CR and film) integration with fix-mounted systems, and also various mobile and portable DR systems;
- the most advanced Csi direct-deposition technology ensures excellence of image quality at low x-ray dose nad improves operation safety;
- | multiple modes of synchronization provides an universal flexibility of system manufacturers and integrators;
- | unique FAED (Full-Field Automatic Exposure Detection) Detector offers a full-imaging-area coverage and low dose sensitivity for accountable operation in challenging situations;
- | 16 bit grey scale;
- | high DQE, high resolution, exceptional dynamic range and very low noise level;
- | reinforced housing made of carbon fiber and aluminum;
- simple and intuitive software in Polish, containing veterinary measurement tools and hints as to the suggested doses of X-ray radiation;
- detectors available in wired and wireless versions (exceptionally long battery life up to 1400 photos with 10-second breaks between exposures, 7-hour standby time;
- | we offer various sizes of detectors (10 x 12 ", 14 x 17", 17 x 17 ")

Cat. No. 283300









## Diagnostic imaging

CARERAY CAREVIEW DETECTORS

09

#### INTERMEDIATE RADIOGRAPHY SYSTEM

#### DIVARIO T2

The extremely fast desktop CR system with a maximum throughput of 73 cassettes per hour. Easy to operate, the Divario DR-T2 offers superior image quality with a pixel size of up to 50  $\mu$ m, a high-speed mode of 5 pixel/mm and includes dicomDX-R PACS® professional image acquisition software with all the necessary image processing tools. With a maximum processing capacity almost 70% faster than previous CR systems, the Divario CR-T2 will reduce patient wait times and increase the efficiency of the examination process. Designed to be unobtrusive and compact, the Divario CR-T2 is the ideal CR system for all standard x-ray examinations.

#### Includes:

- 1 x cassette 1417 incl. imaging plate for images up to 14" x 17" (35 x 43 cm)
- Resolution to 10 pixels / mm, speed up to 73 plates / hour
- | Divario CR-Flash software package:
- | dicomPACS®DX-R Console Software
- | dicomPACS®DX-R DICOM Send SCU
- | dicomPACS®DX-R DICOM Worklist SCU + dicomPACS® Connect
- | dicomPACS®DX-R DICOM Patient CD
- | dicomPACS®DX-R Cognition Optimised Processing

Cat. No. 283024

#### Cassette sizes for type CC:

| 35 x 43 cm (14" x 17") | 35 x 35 cm (14" x 14") | 24 x 30 cm (10" x 12")

- | 18 x 24 cm (8" x 10")
- | 15 x 30 cm (6" x 12")





OPERATION OF THE ACQUISITION SOFTWARE

#### Job planning



X-ray positioning guide





#### **Generator control**



The generator panel displays all values and settings (kVp. mAs, focus etc.) recommended for a specific examination

## Diagnostic imaging





Preview of the X-ray image and worklist



Preview of the current X-ray image

INTERMEDIATE RADIOGRAPHY SYSTEM

#### DIVARIO T2


#### INTEGRATED PROFESSIONAL VIEWER



Completely integrated dicomPACS® viewer for image diagnosis.



Comprehensive search tools enable the comparison of X-ray examinations of one or more patients.



Useful tools such as the configurable measuring magnifier make diagnosis much easier



An integrated prosthesis documentation module provides preoperative planning (optional).








SPECIAL FUNCTIONS FOR DIGITAL X-RAY IMAGING








TTA measuring tool



TPLO measuring tool




Buchanan's Vertebral Heart Score

164

## Diagnostic imaging

#### INTERMEDIATE RADIOGRAPHY SYSTEM



#### DIRECT RADIOGRAPHY SYSTEM

#### EVS 4343

Robust and safety design against shock and drop

- High definition images by direct deposition CslFast image acquisition time less than 2 sec.
- | Highly reliable and stable genrad.
- Highly reliable and stable genrad

#### Cat. No. A35000007

Model	EVS 4343 (Wired)
Detector Type	Direct deposition CsI or Gadox
Weight	4.5 kg
Active Area	430 X 430 mm
Pixel Pitch	140 µm
Resolution	3,072 X 3,072
A/D Conversion	14 bit
Input Voltage	DC 12V, 5A
Communication	Giga Ethernet
X-ray I/F	Lossless AED / Sync Trigger



#### EVS 2430W

- $\mid\,$  High definition imaging with smallest 76  $\mu m$  pixel
- | Ultimate portability
- | WPCS Power System (Wireless Charging)
- | Smart information by OLED display
- $\mid$  Light weight 1.9 kg (w. battery) and fancy design
- | Various applications available
- (Genrad., Mobile, Neonatal/Pediatric, Equine, Vet. Portable, NDT, etc)

#### Cat. No. A44000000

Model	EVS 2430W (Wireless)
Detector Type	Direct deposition Csl or Gadox
Weight	1.9 kg
Active Area	233.47 x 291.84 mm
Pixel Pitch	76 µm
Resolution	3,072 X 3,840
A/D Conversion	16 bit
Input Voltage	DC 12V, 5A
Communication	Giga Ethernet / IEEE 802.11n (5 GHz)
X-ray I/F	Lossless AED



#### EVS 3643

- | Ultimate portability
- | WPCS support (wireless battery charging)
- | Light weight 2.98 kg and robust design
- | Water proof function for increased reliability [IPX4]

#### Cat. No. A26000005

Model	EVS 3643 (Wireless)
Detector Type	Direct deposition Csl or Gadox
Weight	2.98 kg
Active Area	358 X 430 mm
Pixel Pitch	140 µm
Resolution	2,560 X 3,072
A/D Conversion	14 bit
Input Voltage	DC 12V, 5A
Communication	Giga Ethernet / IEEE 802.11n (5 GHz)
X-ray I/F	Lossless AED / Sync Trigger

#### EVS 2430

- | High definition imaging with smallest 76µm pixel
- | High definition images by direct deposition Csl
- | Highly reliable and stable genrad.
- | Fast image acquisition
- | Embedded x-ray trigger interface
- | Designed to withstand shock and drop

Cat. No. A17000000

Model	EVS 2430 (Wired)
Detector Type	Direct deposition CsI or Gadox
Weight	1.9 kg
Active Area	233.47 x 291.84 mm
Pixel Pitch	76 µm
Resolution	3,072 X 3,840
A/D Conversion	16 bit
Input Voltage	DC 12V, 5A
Communication	Giga Ethernet
X-ray I/F	Lossless AED/ Sync Trigger

## Diagnostic imaging

#### DIRECT RADIOGRAPHY SYSTEM





#### **INNOVATIVE TECHNOLOGIES FOR LIMITLESS PORTABILITY**

#### EVS 3643 with Water Proof [IPX4]

EVS 3643 has the water resistance [IPX4] feature which provides protection against blood, sudden rain or othercontaminants which can aect the detector performance in mobile applications.

EVS 3643 delivers reliable performance in any environments and unexpected situations providing perfect portable solution for x-ray diagnostic imaging.



#### Protection Suits

Available in three dierent designs to match varying user requirements and to protect your Exprimer detector in style from various environmental hazards. With an ergonomic design, you can achieve optimum usability in any x-ray environments. Combined with a tablet mounter, you can heighten your productivity with one compact package in any portable situations.

#### Slim Dual Battery Charger

Simultaneous charging to two batteries is possible with the dual battery charger. With its horizontal and vertical insertion design, this charger supports both EVS 3643 and EVS 2430W batteries.

Supported with 12V, charging on-the-go is possible as it connects with the cigar outlet of any car. The charger comes standards with EVS 3643 and EVS 2430W detectors.



#### PCP (Portable Console PC)

Instanct digital upgrade is possible when PCP is combined with EVS 3643 or EVS 2430W. Eliminate the need for additional control box between the detector and the PC, and acquire digital images instantly with DRTECH's PCP. Inside AP maximizes the portability of Exprimer's portable detectors by

allowing direct communication with tablet PCs and smart phones for image viewing on mobile devices. With a PCP and wireless communications, high definition images can be acquired anywhere and anytime allowing for unlimited outdoor and mobile imaging.

## Diagnostic imaging

DIRECT RADIOGRAPHY SYSTEM






#### **BENEFIT FROM EXCELLENT QUALITY X-RAY IMAGING WITH EXPRIMER**

Exprimer, the innovative digital x-ray solutions suit multiple diagnostic environments with dierent needs.

- | Excellent image quality using direct deposited Csl
- | Ultimate sharpness image by TRUVIEW® ART
- Instant upgrade to digital mobile X-ray system
- Patient dose reduction with reliable Lossless AED
- High resistance to impact and vibration
- Low price fixed grid (120 lines / inch)
- Light weight and durable design for portable applications



Lossless AED' promising 'Dose Reduction' & 'TRUVIEW® ART' promising sharper and more accurate images with '30% higher MTF'



## Intensified image sharpness with directly deposited CsI

Directly deposited CsI can provide clearer images at lowest dispersion compared to conventional CsI and GOS scintillator. High quality images are not a result of any one feature such as smallest pixel size or low electronic noise, but achieved when all components of the detector are optimized to operate in harmony with each corresponding specification

## Better DQE performance in higher spatial frequencies

Exprimer with its well oriented direct deposition columnar structure CsI + TFT has high DQE performance providing outstanding high quality images. It also demonstrates comparably excellent DQE\* performance in high spatial frequency range.

\*DQE measurement condition: RQA-5 (Typical 2.6  $\mu$ Gy)

## Faster image display with high speed operating scheme


Fast image preview and display time of Exprimer products allow for more eective and ecient operation leading to increased productivity. Image preview in both wired and wireless modes are achieved less than 2 seconds and a full image is acquired in 4.5 and 6 seconds respectively.

#### Best image processing solution for fine tuned quality images ECONSOLE1 & TRUVIEW® ART

ECONSOLE1 (UI Software) and TRUVIEW® ART (image sharpening algorithm) provides a perfect solution to increase diagnostic productivity and accuracy. With easy to use convenient user interface design and powerful image processing engine, ECONSOLE1 & TRUVIEW® ART enables more accurate diagnosis with high quality and highly defined images.

## Diagnostic imaging

#### DIRECT RADIOGRAPHY SYSTEM



SCANNER / DENTAL GENERATOR

#### Midmark Digital Dental System

- | Midmark DC is a high frequency veterinary X-ray system designed for digital imaging;
- I the DR sensor has been designed to enable the visualization of both soft tissues as well as the crown and root, which are necessary to provide patients with the highest quality dental care;
- | the cable is reinforced, designed to increase its durability;
- waterproof, hermetic gasket eliminates problems occurring during routine cleaning;
- CMOS technology provides better image quality;
- fast USB connection is designed to take photos 37% faster;
- the rounded corners of the sensor allow for easier positioning in tight places in the mouth to better grasp the roots;
- Progeny software provides excellent image quality, has an intuitive and logical interface, simplified acquisition. It allows you to: save images in dicom format, backup and restore data, 22 practical and useful tools, create your own custom templates, connect to most major practice management software.
- The product is available in a mobile or wall-mounted version, different arm lengths available

Cat. No. 271528







#### Ultrasound Q3Vet

#### ADVANCED DOPPLER MODEL

#### Applied technologies:

DSC (Digital Scan Converter), DBF (Digital Beam Forming), RDA (Real-time Dynamic Aperture), DRA (Dynamic Receiving Apodization), DRF (Dynamic Receiving Focusing), DFS (Dynamic Frequency Scan) 32Ch, THI (Tissue Harmonic Imaging), SNR (Speckle Noise Reduction), 3D/4D Imaging Module, Pure digital time-delay ASICs, Auto Flow Angle calibration, Extreme Resolution adaptive imaging, Transducer Signal Processing

#### Functions:

- one touch Optimization;
- 15" ekran wysokiej rozdzielczości + 8.4" ekran dotykowy;podgrzewacz do żelu;
- auto pomiary IMT;
- B, Dual B, Quad B, M, B/M;
- CFM, PDI, PWR D, THI;
- Duplex/Triple synchronous;
- dwa lata gwarancji;
- dostępne sondy: Convex, Microconvex, Liniowa, Rektalna, 4D volume,

#### Cat. No. ZON\_300

#### V3Vet Veterinary Color Doppler Ultrasound


- | 12.1 inch LED monitor with 90 degree tilting;
- 3 activated transducer connectors, transducer auto-freeze protection;
- build-in Lithium-in battery supports more than 1.5 hours continues work;
- outstanding portable model, compact and streamlined design weighing 7.5 kg;
- | 320G hard disk for large capacity storage;
- professional measurement software package applicable for general and OB diagnosis of ve such as dog, cat, bovine, horse, shepp etc.;
- data port: USB 2.0, Video, S-Video, RS-232, VGA.

Cat. No. ZON\_200

## **Diagnostic imaging**

## ULTRASOUNDS





#### ULTRASOUNDS

#### Ultrasoud i50V

#### Function and configuration:

- | 12.1 inch High resolution medical monitor with 90 degrees angle rotation;
- | double beamform imaging technology;
- outstanding portable model, compact and streamlined design;
- | all-round adjusting knob, afficient and fast operation;
- one key full screen design (FSI) for better visual field
- support multiple measurement software packages;
- optional Trolley, Built-in lithium battery and suitcase, convenient for outdoor diagnosis;
- support multiple data transmission interface;
- auto generated BMP format report, support USB import and export;
- support laser printer module, directly connect the laser printer to output report;
- suppored Transducers: Convex, Linear, Transvaginal, and Micro-convex.

Cat. No. ZON\_100



#### Chison USG Eco 5 Vet

#### Features and technologies:

- | ultra long battery life over 2 hrs;
- | full digital beam-former;
- | colour doppler technology;
- | advanced fast frame rate technology;
- colour filter;
- applications: abdomen, OB&GYN, Urology, Vascular, Small Parts, Pediatrics, Basic Cardiology
- | display modes: B, B/B, 4B, B/M, CFM, PW, trapezoidal;
- i-Image, THI, SRA;
- compound imaging,
- rotatable LED monitor (0-30°)
- USB, VGA;
- Wide range of probes (2.5 11 MHz).

Cat. No. CHIS\_ECO5VET

## Diagnostic imaging

#### ULTRASOUNDS



#### ULTRASOUNDS

#### Chison USG Sonobook 6 Vet

#### Features and technologies:

- | 15 inch HD LED monitor with 120° rotatable;
- Compact and robust laptop design with 5.5 kg for challenging VET environment and easy mobility;
- | 2 hrs battery life in active mode for outdoor VET scanning
- Professional VET software and wide range of VET transducers;
- dedicated VET calculation package, body mark and annotation;
- | FHI for deep penetration of large animals;
- | Q-bean for rapid breath of small animals;
- X-contrast, Q-flow, Q-image, Super Needle;
- Comprehensive cardiac measurement packages;
- Dedicated trolley with triple probe connectors, adapter space, printer socket

Cat. No. CHIS\_SONOBOOK6



#### Chison USG Sonobook 8 Vet

#### Features and technologies:

- | 15 inch HD LED monitor with 120° rotatable;
- Compact and robust laptop design with 5.5 kg for challenging VET environment and easy mobility;
- | 2 hrs battery life in active mode for outdoor VET scanning
- Professional VET software and wide range of VET transducers;
- dedicated VET calculation package, body mark and annotation;
- | FHI for deep penetration of large animals;
- | Q-bean for rapid breath of small animals;
- X-contrast, Q-flow, Q-image, Super Needle;
- Comprehensive cardiac measurement packages;
- Dedicated trolley with triple probe connectors, adapter space, printer socket

Cat. No. CHIS\_SONOBOOK8

## Diagnostic imaging

## 09

#### ULTRASOUNDS



#### ULTRASOUNDS

#### Chison USG EBit 60 VET

#### Features and technologies:

- | 15 inch HD LED monitor , rotatable;
- Compact and robust laptop design with 7.5 kg for challenging VET environment and easy mobility;
- | 2 hrs battery life in active mode for outdoor VET scanning
- | professional VET software ;
- dedicated VET calculation package, body mark and annotation;
- | FHI for deep penetration of large animals;
- | Q-bean for rapid breath of small animals;
- X-contrast, Q-flow, Q-image, Super Needle, Hip Graf;
- Ports: USB, LAN, VGA, CVI, Video, Remote
- | Wide range of probes (15-18 MHz)
- | Dedicated trolley (option)

Cat. No. CHIS\_EBit60VET



#### Chison Ultrasonograf QBit7 VET

#### Features and technologies:

- | 19 inch HD LED monitor , rotatable;
- | professional VET software and wide range of VET probes;
- dedicated VET calculation package, body mark and annotation;
- | FHI for deep penetration of large animals;
- | Q-bean for rapid breath of small animals;
- X-contrast, Q-flow, Q-image, ;
- | advanced cardiac package;
- | Hassle-Free maintenance (Hero Kit);
- | streamlined workflow (6-one-key step);
- everyday efficiency;
- | walk-up, usability, ergomonics and mobility
- | excellent features: Stress ECO and elastography

Cat. No. CHIS\_QBit7VET

## Diagnostic imaging

#### ULTRASOUNDS



#### ULTRASOUNDS

## Setting A New Standard For HD+ Video

equipped with a sophisticated built-in workstation and advanced image processing technology, the new HD-500 image processor provides high definition images for endoscopic examinations. The HD-500 is ideal for practitioners who require extraordinary performance during their daily endoscopic examinations and challenging research.

#### Easy operation

Smart workflow makes the endoscopic examination simpler.

#### HD+ Image Processing

Processing over one megapixel image with advanced technology

Intelligent built-in workstation HD+ image and patient information management makes for a more efficient examination.

Cat. No. SONO\_HD500



#### ENDOMED Arthroscopy

#### FULL STATION WITH TOOLS

#### Included:

- | ENDOcam 100, Digital Video single CCD, without TV-Adapter
- ENDOlite 200, Wysokowydajne źródło światła LED
- ENM 20-191, Monitor 19" LED, Roz. 1280x1024 HD
- ENDOshaver 6000 , Shaver System with 8000 rpm
- ENDOpump 100, Pompa do artroskopii, 10-150mmHg .0.1-2.0
- Endomed Wózek model S-move

#### Tools:

- | Full radius resector, aggressive, Ø3,5mm
- | Meniscus cutter, Ø3,5mm
- Full radius resector, Ø4,2mm
- | Full radius resector, aggressive, Ø4,2mm
- | Meniscus Cutter, Ø4,2mm
- Round Burr, Ø5.5mm







- Arthroscope Ø4mm, View 0°, angle 95°, długość 175mm
- Arthroscope Ø4mm, View 30°, angle 95°, długość 175mm
- Sheath arthroscopy, two stopcocks, Ø4mm, 30°, length 175mm
- | Obturator arthroscopy, pyramidal, Ø4mm, długość 175mm
- | Obturator arthroscopy, blunt, Ø4mm, długość 175mm
- Hook probe 3mm, graduated
- Scalpel knife
- Sickle knife
- Grasping forceps, Ø3,5mm, długość 130 mm
- Hook punch with scoop, 30° right, Ø3,5mm
- Hook punch with scoop, 30° left, Ø3,5mm
- Wide jaw punch medium, Ø3,5mm, straight
- | Wide jaw punch, medium, Ø3,5mm, 90°, left angled
- Wide jaw punch, medium,  $\phi$ 3,5mm, 90°, right angled
- Hook Scissors Straight ,Ø3,5mm, length 130 mm

Cat. No. ENDO ART

#### ENDOSCOPES

#### ENDOMED Combo-3

#### FULL ENDOSCOPIC STATION

#### Included:

truck;

light source / processor C3P13, ENDOMED COMBO-3, DVI;
 1 × LED;

- | 1 x C3G13, ENDOMED COMBO-3, 1050 mm Gastroscope,
- $| \phi$  9.6mm, working channel  $\phi$  2.8 mm;
- | 1 x C3C16, ENDOMED COMBO-3, Colonoscope 1650 mm,
- $| \phi$  11.8mm, working channel  $\phi$  3.2 mm;
- | medical monitor Monitor 19 "LED, Res 1280x1024 HD, DVI-I, medical;
- Vaco SP120, ENDOMED Suction pump 2 liters.

Cat. No. ENDO\_COM



## Endoscope Sonoscape HD-350

- Creates clear images through the use of edge enhancements and reducing image noise reduction.With the CMOS image sensor, HD-350 gives high quality images with low noise levels, which leads to improves clarity and optimized version is not only routine procedures but also in complex interventions.
- | To increase efficiency, patient data management is accessible with the press of a buton. The smart workflow of the built-in workstation of HD-350 helps doctor to edit reports, record videos and manage patient files with ease.
- The waterjet function of the scopes help to maintain a clear endoscopic view by washing away blood or mucous.

Cat. No. SONO\_HD350





## Anaesthetics

# Anesthesia device with DM-6B ventilator

# Mindray 8VET Monitor





MORE INFORMATION ON PAGE 186

## Anaesthetics

10



## MORE INFORMATION ON PAGE 191

#### Matrx<sup>®</sup> VME2

10

Sleek, stylish and featuring a chrome design, the Midmark Matrx<sup>®</sup> VME2 incorporates all veterinary anaesthesia equipment onto one post. A portable stand and 5-leg spider base makes for convenient transportation and easy mobility. Optional add-on components include shelving, a second flowmeter, a second vaporizer, an occlusion valve with safety relief and/or an emergency air intake valve.

Cat. No. 271516



#### Anesthesia device with DM-6B ventilator

- easy to use and move construction on 4 wheels;
- $\mid\,$  fan operating in 3 modes: IPPV; A / C; SIMV.7 "TFT LCD;
- | the screen displays breathing parameters and alarms;
- | supports animals up to 200 kg;
- oxygen flow: from 0 to 5l / minute;
- | tidal volume (Vt) 0.20 ml 1500 ml;
- | frequency 1bpm ~ 100 bpm.

Cat. No. 271599



#### Moduflex Access 2 Veterinary Anaestesia Machine

- 0.1 to 4 LPM O2 flowmeter ;
- rotary float and large scale for easy readability;
- recessed flush valve to prevent accidental activation;
- both the anesthesia machine head and the absorber head are precision machined to reduce the amount of connections which greatly reduces the risk of leaks;
- quick-connect fresh gas outlet for easy and secure switching between breathing circuits;
- cage mount compatibility manifold for easy and quick installation of vaporizer;
- | five wheel spider base;
- 1500 ml CO2 absorber designed for easy and fast replenishing of sodalime, ondensation trap to prevent water build-up in sodalime, easy and efficient gas sealing.;
  5 year warranty.

Cat. No. DISPOMED\_MODUFLEXACCESS2

#### Coaxial Veterinary Anaestesia Machine

- | 0.1 to 4 LPM O2 flowmeter;
- rotary float and large scale for easy readability;
- | recessed flush valve to prevent accidental activation;
- quick connect fresh gas outlet for easy and secure switching between breathing circuits;
- SELECTATECTM compatibility manifold for easy and quick installation of the vaporizer;
- | 1500 ml CO2 absorber with integrated coaxial breathing circuit reduces patient heat loss, helps to humidify the inspired gas, and reduces clutter, condensation trap to prevent water build-up in soda lime. Easy and efficient gas sealing
- momentary closed feature incorporated into the pop-off valve
- convenient anesthesia monitor shelf
- 5-year warranty.

Cat. No. DISPOMED\_COAXIA

## Anaesthetics

ANESTHETIC APPARATUS





#### MANUAL RESUSCITATORS / RESPIRATION BAGS

| disposable self-inflating bag;

- | complete with a transparent disposable mask;
- | superb fitting to patient's airways;
- latex-free;
- capacity: 1600 cc.

Cat. No. 34277

#### Manual resuscitator

silicone;

- transparent;
- delivered with a silicone mask;
- latex-free;
- capacity: 1600 cc.

Cat. No. **34245** 



#### Double-chamber manual resuscitator

professional case; | self expanding internal chamber; capacity: 1900 cc; | fully compatible; reusable.

Cat. No. 34240



| made of black antistatic rubber;

various sizes.

Cat. No.	
271600	Respiration bag 0,5 l
271610	Respiration bag 11
271620	Respiration bag 1,5 l
271630	Respiration bag 2 l
271640	Respiration bag 3 l
271645	Respiration bag 4 l
271646	Respiration bag 5 l



#### Maska anestetyczna

| made of plastic;

- | with a rubber diaphragm;
- | for use with manual resuscitators;
- for dogs and cats..

Cat. No. 271672, 271673

- | made of a high quality material;
- with balloon and connector;
- solid and flexible;
- autoclavable at 134°C.


Cat. No.	Name
272020	3,0 mm rubber endotracheal tube
272021	3,5 mm rubber endotracheal tube
272022	4,0 mm rubber endotracheal tube
272023	4,5 mm rubber endotracheal tube
272024	5,0 mm rubber endotracheal tube
272025	5,5 mm rubber endotracheal tube
272026	6,0 mm rubber endotracheal tube
272027	6,5 mm rubber endotracheal tube
272028	7,0 mm rubber endotracheal tube
272029	7,5 mm rubber endotracheal tube
272030	8,0 mm rubber endotracheal tube
272031	8,5 mm rubber endotracheal tube
272032	9,0 mm rubber endotracheal tube
272033	9,5 mm rubber endotracheal tube
272034	10,0 mm rubber endotracheal tube
272035	11,0 mm rubber endotracheal tube

## Anaesthetics

10

#### ANAESTHETIC MASKS / ENDOTRACHEAL TUBES





#### Silicone endotracheal tube

#### | with balloon and connector;

- solid and flexible;
- | autoclavable at 120°C.

Cat. No.	Name
271580	3,0 mm silicone endotracheal tube
271581	5,0 mm silicone endotracheal tube
271582	6,0 mm silicone endotracheal tube
271583	7,0 mm silicone endotracheal tube
271584	8,0 mm silicone endotracheal tube
271585	10,0 mm silicone endotracheal tube
271586	11,0 mm silicone endotracheal tube
271587	12,0 mm silicone endotracheal tube
271588	14,0 mm silicone endotracheal tube



#### Rurka dotchawicza PVC

#### | with balloon and connector.

Cat. No.	Name
272083	3,5 mm PVC endotracheal tube
272084	4,0 mm PVC endotracheal tube
272085	4,5 mm PVC endotracheal tube
272086	5,0 mm PVC endotracheal tube
272087	5,5 mm PVC endotracheal tube
272088	6,0 mm PVC endotracheal tube
272089	6,5 mm PVC endotracheal tube
272090	7,0 mm PVC endotracheal tube
272091	8,0 mm PVC endotracheal tube
272092	8,5 mm PVC endotracheal tube
272093	9,0 mm PVC endotracheal tube
272094	9,5 mm PVC endotracheal tube
272095	10,0 mm PVC endotracheal tube
272096	10,5 mm PVC endotracheal tube
272097	11,0 mm PVC endotracheal tube



#### Mindray 8VET 8.4" Monitor

8.4" High-resolution LED display with 8 waveforms. Portable, less than 3.2 kg· No-fan design for a quiet care environment. Easy to use for intuitive operation. Up to 4 hours continuous monitoring with large capacity Li-ion battery. Trust for veterinary use, parameters including ECC, RESP, SpO2, NIBP, PR, TEMP, IBP, CO2· 360-degree visible alarm indicator. 120 hours trend and 48 hours waveforms reviewing

#### Standard Configuration

· ECG, RESP, NIBP, Mindray SpO2, TEMP, Li-ion battery

#### Optional

 $\cdot$  Nellcor SpO2  $\cdot$  2-ch IBP CO2 (Sidestream, Microstream or Mainstream)  $\cdot$  Thermal recorder

Cat. No. 240120

#### Zoncare PM-7000 VET

- 12.1 inch LED screen, high brightness, Low power consumption and long lifetime
- Standard 6 parameters configurations: ECG. RESP, NIBP. SP02. PR. TEMP
- Large font display offers a clear view of all vital signs
- Support AC and DC, build-in Lithium-ion battery, can provide the power to the monitor a
- | long time after the power off
- Precise measurement technology, professional veterinary measurement software, and
- suitable for the long time monitoring on dog, cat and other animals
- | Large-capacity storage and playback function

Cat. No. Pm7000 - 240120

#### Mindray PM 60Vet pulse oximeter

- | Suitable for adult, pediatric and neonatal patients
- | High resolution 2.4" color LCD display
- | Two selectable working modes: spot check and continuous monitoring
- | Adjustable audible and visual alarms
- Real time data transfer by infrared technology
- | Batteries to support up to 36 hour continuous working
- Protective cover and carrying case perfect for portable monitoring needs

Cat. No. Mindray23

## Anaesthetics









POMPS

10

#### Mindray Tiva SK-500 II Syringe Pump

- three kinds of modes: Rate Mode, Time Mode, Body Weight Mode;
- double CPU, making the process of injection safer and more reliable;
- accurate and sensitive pressure sensor detecting occlusion pressure accurately;
- injection time, dose, solution can be checked during the process of injection;
- precision step motor, there is no pulsating wave even during low-rate injection;
- excellent injection operation, the accuracy will be ensured even at small-volume injection;
- accept all standard syringes of 5ml, 10ml, 20ml, 30ml and 50ml;
- alarm is given when there is malfunction or operation error, and the machine stops running automatically;
- Save the preset volume automatically when there is unpredictable switch off, ensuring the accuracy of injection volume.

Cat. No. MINDRAY30



#### Mindray Infusion Pump SK-6001

- | Middle Detection and Ultrasound Detectors to double guarantee the correct IV set loading Flow rate mode (SK-600I: 1-600 ml/h)
- anti-bolus system to reduce significantly bolus after occlusion sudden release.
- | reminder alarm in case of forgetting to start the operation.
- high-speed infusion up to 2000ml/h satisfies the special clinical scenario
- | (Optional) drip sensor and fluid heater
- wireless data transmission with an infusion supervision system.

Cat. No. MINDRAY1





# Dentistry

# Highdent Trio Plus Veterinary Dental Unit

# Highdent Quattro Plus Veterinary Dental Unit





## Dentistry

11

MORE INFORMATION ON PAGE 198



#### DENTAL UNITS

11

#### Highdent Intro Veterinary Dental Unit

- low risk of breakage and reduced risk of injury due to its recessed handpieces;
- low risk of leaks due to its innovative piping system configuration;
- 5 year warranty on the veterinary dental unit and 2 years on Piezo;
- | ergonomic Design: adjustable height;
- | Highdent<sup>™</sup> dental units include a muffler to reduce noise from the handpieces;
- | unique tabletop design for instruments;
- | available in different colours.

Cat. No. DISPOMED\_INTRO



#### Highdent Trio Plus Veterinary Dental Unit

- the new and easy to clean shelf give you more space for your instruments;
- | height and handpiece console are adjustable;
- | clear indication of handpieces;
- low risk of breakage and reduced risk of injury due to its recessed handpieces;
- low risk of leaks due to its innovative piping system configuration;
- | automatic handpiece selection eliminates the risk of error;
- | 5 year warranty on the machine and 2 years on Piezo
- Highdent<sup>™</sup> dental machines include a muffler to reduce noise from the handpieces
- water and air flow adjustments for individual handpieces offer better control.

Cat. No. DISPOMED\_TRIO

## Dentistry

#### DENTAL UNITS



#### DENTAL UNITS

11

#### Highdent Quattro Plus Veterinary Dental Unit

- the new and easy to clean shelf give you more space for your instruments;
- | height and handpiece console are adjustable;
- | clear indication of handpieces;
- low risk of breakage and reduced risk of injury due to its recessed handpieces;
- low risk of leaks due to its innovative piping system configuration;
- automatic handpiece selection eliminates the risk of error;
- 5 year warranty on the machine and 2 years on Piezo
- | Highdent<sup>™</sup> dental machines include a muffler to reduce noise from the handpieces
- water and air flow adjustments for individual handpieces offer better control
- | Detachable suction hose for easier cleaning.

Cat. No. DISPOMED\_QUATTRO



#### Scaler Woodpecker UDS-K-LED

- ompatible with EMS's handpiece and tips;
- he handpiece is cast from Titanium alloy to ensure sound mechanic function and durabiltiy;
- assembled with excellent piezo-ceramic, the handpiece has high and long lasting power and generates little heat;
- detachable handpiece, can be autoclaved under the high temperature of 135° C and pressure of 0.22Mpa;
- | five pieces of scaler tips attached, 2pcs of G1, 1 pc of G2, G4, P1;
- automatic frequency tracking ensures that the machine always work on the best frequency;
- the adjustment range of the amplitude is much wider than before;
- | to eliminate the dental calculus, bacterial plaque, and tea scales on the teeth quickly without hurting the gum and the enamel;
- digitally controlled, easy operation and more efficient for scaling.

Cat. No. 271305 \_UDSKLED

#### Skaler Woodpecker UDS-L LED z butelką

- | four fuctions:Scaling.endo,automatic water supply system(optional) and LED;
- auomatic frequency tracking ensures that the machine always works on the best frequency and performs more steadily :
- the handpiece is detachable and can be autoclaved to the high temperature of 135 and pressure of 0.2Mpa;
- digital control of Woodpecker dental ultrasonic scaler ,easy operation and more efficient for scaling;
- the curing light can be used along with BSC-US13 Ultrasonic scaler, utilizing its ultrasonic power to provide electricity for the curing light.

Cat. No. 271305\_UDSLLED

## Dentistry







#### DENTISTRY LUKSATOR / GROOVED HANDLES

#### Dentistry luksator

11

erfect for cutting the fibers of the cervical periodontium and widening the periodontal space;

| do not use to lift a tooth.

Cat. No.	Name
175317	Luksator 1 mm
175318	Luksator 2 mm
175319	Luksator 3 mm
175320	Luksator 4 mm
175321	Luksator 5 mm



#### Mini Friedmann Rongeurs 1.3 mm

| horizontal grooved handles to enhance grip;

- | double spring mechanism allows self-opening action;
- | tip width: 1.3 mm, lenght: 140 mm.

Cat. No. 177314



#### Mini Friedmann Rongeurs 2.0 mm

horizontal grooved handles to enhance grip;

| double spring mechanism allows self-opening action;

| tip width: 2.0 mm, lenght: 120 mm.

Cat. No. 177312

Rodent Molar Cutter with Tungsten Carbide

| suitable for trimming long molars leaving teeth with a

smooth surface;

| 170 mm lenght.


Cat. No. 177316

## Dentistry

11

#### GROOVED HANDLES / PINCERS





## PINCERS / DENTAL ELEVATOR

#### Rodent Molar Cutters

helps to debulk spurs and hooks on molars;180 mm lenght.

Cat. No. **177315** 



#### Fahrenkrug Dental Elevator

dedicated to help extract deciduous teeth;For dogs - 200 mm lenght (175801);

| For dogs and cats - 190 mm (175802).

Cat. No. 175801, 175802



#### Dental Extraction Forceps

diamond coated;

| 110 mm lenght.

Cat. No. 176010

#### Dental Extraction Forceps

premolars straight;

| 110 mm lenght.

Cat. No. 176011

## Dentistry

11

## DENTAL EXTRACION FORCEPS





#### DENTAL EXTRACION FORCEPS

#### Dental Extraction Forceps for dogs

roots;

| 165 mm lenght, straight.

Cat. No. **176016** 



## Dental Extraction Forceps for dogs

| premolar; | 165 mm lenght, diamond coated.

Cat. No. 176015



#### Dental Extraction Forceps for dogs

| premolar; | 165 mm lenght, straight.

Cat. No. 176017

#### Witzel Root Extraction Forceps

50 mm lenght.

Cat. No. 176215

## Dentistry

11

## DENTAL EXTRACION FORCEPS





## DENTAL EXTRACION FORCEPS

#### Dietrich Tartar Removing Forceps

| 165 mm lenght.

Cat. No. 176116



## Dental Forceps for cutting roots

| 150 mm lenght.

Cat. No. **176024** 



Extraction forceps

150 mm lenght.

Cat. No.. 176023

#### Extraction oblique forceps

for cats and small dogs125 mm lenght.

Cat. No. 176025

## Dentistry

11

## DENTAL EXTRACION FORCEPS





PROPHY PASTE / PROPHY CUPS / PERIODONTAL PROBES

#### Prophy paste for animals

200 per box.

Cat. No. Dispomed7



#### Prophy cups

| 100 per box.


Cat. No. Dispomed6



#### Periodontal Probes

| 170 mm lenght; | graduated tip.

Cat. No. 17734011



| 160 mm lenght;

graduated tip.

Cat. No. 17647202

#### Periodontal Probes

160 mm lenght;graduated tip.

Cat. No. 176456

#### Periodontal Probe/Explorer Size 9

| 155 mm lenght.

Cat. No. 176459

## Dentistry

#### PERIODONTAL PROBES







#### PERIODONTAL PROBES / SCALERS

#### Periodontal Probe/Explorer Size 17


| 155 mm lenght; dedicated for cats.

Cat. No. **176467** 



| size 3/4"; suitable for incisor teeth; length: 170 mm.

Cat. No. 17647203



# | 165 mm lenght. Cat. No. **176410**

| 165 mm lenght.

Cat. No. 176412

#### Dental Mirror

| 165 mm lenght.

Cat. No. 176422

| suitable for Chinchillas;

lightweight;


length: 115mm.

Cat. No. 177303

## Dentistry

11

## SCALERS / DENTAL MIRROR / DILATORS







#### DILATORS

11

#### Mouth Gags for Rodent

suitable for medium to large breed rabbits,robust

| length: 110 mm.

Cat. No. **177304** 



#### Mouth Gags for Rodents

| cam ratchet mechanism; | robust; | length: 110 mm.

Cat. No. **177330** 



#### Rodent Cheek Dilator

| XS size; | for guinea pigs and chinchillas; | blade dimension: 16 x 13 mm (L x W).

Cat. No. **177302** 



#### Rodent Cheek Dilator

S size;

for dwarf rabbits;

| blade dimension: 14 x 11 mm (L x W).

Cat. No. 177306

#### Rodent Cheek Dilator

M size;

| for medium ans small breed rabbits;

| blade dimension: 21 x 15 mm (L x W).

Cat. No. 177305

## Dentistry







#### DILATORS

11

#### Rodent Cheek Dilator

| L size;| blade dimension: 35 x 10 mm (L x W).

Cat. No.. **177310** 



#### Cheek Dilator

suitable for small pets;

| blade dimension: 11 x 9 mm.

Cat. No. 177320

## Rodent Cheek Dilator

suitable for chinchillas;blade dimension: 21 x 9 mm.

Cat. No. **177311** 



## Dental Mouth Gag

- designed for oral examination and to facilitate dental procedures;
- stainless steel construction with spring action and rubber inserts to protect teeth;
- | available in three sizes.

Cat. No. 200117, 200110, 200115

## Dentistry







Dentistry

#### JAW MODEL

11

### Dental Model – cat

- | transparent hinged jaw helps visualise the roots of all teeth (incisors, canine, premolars, molars);
- $\mid$  can be separated for closer examination.

Cat. No. **174092** 



### Dental Model - dog

transparent hinged jaw helps visualise the roots of all teeth (incisors, canine, premolars, molars);
can be separated for closer examination.

Cat. No. **174093** 



### Dental Model - rodent

transparent hinged jaw helps visualise the roots of all teeth (incisors, canine, premolars, molars);
can be separated for closer examination.

Cat. No. **174091** 







# Stainless steel cage

# Emaled 560/560 lamp





MORE INFORMATION ON PAGE 235

## Clinic equipment



MORE INFORMATION ON PAGE 223

Clinic equipment

LAMPS

12

#### Saturno lamp

- suitable for minor surgery and first aid;three versions available (mobile on wheels, mountable on
- the wall or on the ceiling); aluminium reflector;
- luminosity from 1 m: 50 000 lux;
- 9 LEDS 2.5 W;
- $\mid\,$  selection of 2 different white light temperatures: 4000  $^{\circ}\text{K}$  and 4500  $^{\circ}\text{K};$
- easy operation;
- stable;
- | minimum room height: 2400 mm (ceiling version).

Cat. No.	Name
49100	Saturno lamp, on wheels
49102	Saturno lamp, wall-mounted
49103	Saturno lamp, ceiling-mounted
30797	Spare autoclavable handle



#### Pentaled 28 LED LIGHT LAMP

- | suitable for precise surgery;
- | for diagnostics and scientific laboratories;
- bright, cold light;
- | colour temperature: K 4.500/5.000;
- number of LED: 28;
- average LED life: 50.000 hours;
- does not strain the eyes, facilitates concentration;
- | polycarbonate corpus;
- | reflector can be rotated on the vertical and horizontal axis;
- | luminosity from 1 m: 120 000 lux.

Cat. No. 49115

#### Pentaled 12 lamp

- onomical shape;
- suitable for precise surgery;
- | for diagnostics and scientific laboratories;
- bright, cold light;
- | does not strain the eyes, facilitates concentration;
- | polycarbonate corpus;
- reflector can be rotated on the vertical and horizontal axis;luminosity from 1 m: 50 000 lux.

Cat. No. 30787



#### Solesud2 lamp

- halogen light;
- | for use in gynaecology;
- aluminium reflector;
- | 2 halogen bulbs made of tempered glass;
- | luminosity from 50 cm: 40 000 lux.
- on wheels;
- adjusted height.

Cat. No. 30747







12

#### Solesud 2

Suitable for gynaecology, examination and headbed. Aluminium reflector is equipped with 2 halogen dichroic bulbs, 20W 12V IRC protected by a tempered glass, hermatically sealed against dirt and fluid spray and with ergonomic handle for movement. Bulbs switch on simultaneously and, in case one breaks, light is always assured. Arm is made with articulated joint and mechanical frictions. A light beam of 180 mm diameter at a distance of 50 cm allows precise and reliable examinations. Generates a light intensity of: 40000 lux at 50 cm.

#### **Technical Specifications**

Type of lighting: halogen · Consumption (Watt): 2x20 · Life of the bulbs (hours): 5000 · Light intensity at 50 cm (lux): 40000 · Light intensity at 1 m (lux): 8000 · Colour temperature ± 5% at 1 m (K): 3000 · Light field size at 50 cm (\$\overline{\pmm}\$): 180 · External diameter of the reflector (\$\overline{\pmm}\$ mm): 104 · Colour: white · Norms: IEC 60601-1 · Marks: CE93/42
· Weight: 6 kg · Operating voltage: 230 V - 50/60 Hz
Electronic transformer integrated in the base with double safety fuse · Adjustable height: 60-160 cm

Cat. No. 30748

#### Emaled lamp

- high-quality light;adjustable luminosity;
- | long life (> 40 000 h);
- | little heating;
- | luminosity: 50 000 lux.

Cat. No.	Name
30752	Emaled 200F lamp, standing
30753	Emaled 200W lamp, wall-mounted
30754	Emaled 200D lamp, ceiling-mounted



#### Emaled 560/560 lamp

- surgical lamp;
- does not heat up the operating field;
- | the smallest and brightest light in this class;
- I luminosity adjustment form 10% to 100% using the control panel or the pilot;
- easy assembly.

Cat. No. 30755

#### Emaled 300M lamp


- high-quality light;
- adjustable luminosity;
- | long life (> 40 000 h);
- | little heating;
- | luminosity: 100 000 lux.

Cat. No. 30756









LAMPS

12

#### PENTALED 30E SCIALYTIC LIGHT – trolley

| A compact size hi-tech product that guarantees great performance levels. It has 30 elliptical reflectors split into 6 modules each containing 5 LEDs, providing 160,000 LUX- selectable 4,500 °K - 5,000 °K. By simply pressing the digital key, the surgeon can choose two different white light temperatures- E-view. Innovative lens allows the surgeon to expand/narrow the illuminated operating field electronically, without altering intensity of light.-E-comfort.The indirect LED's creates 3D shadowless illumination of the operating field, no risk of dazzling-E-white. Only white LED's to avoid risk of perception of colour temperature variations or coloured shadows-E-deep. Centre of light is fitted with extra LEDs module specifically designed to reflect deep light. This means that surgeon can work with perfect 3D lighting.Battery and high definition camera on request. Optional battery group with automatic charger (4 hours autonomy) and high definition integrated camera (for ceiling and wall models only).


Cat. No. 30780



#### PRIMALED LIGHT - trolley

LED life: 50.000 hours- Light intensity: 80.000 lux at 50 cm, 20.000 lux at 1 m· Colour temperature: 4.500-5.000 °K· Colour rendering index: 95· Illuminated field at 50 cm:
Ø 140 mm· External diameter of the reflector: 190 mm· Temperature increase in the operating area at 1 m: 1°C· Adjustable height with castors: 600 mm - max 1600 mm· Power supply: 81-240V - 50/60 Hz· Electric absorption: 21VA· Weight: 10 kg (trolley included)

Cat. No. 30777



#### Instrument table

painted;

- | railing on three sides;
- white PVC shelves;
- dimensions: 50 x 60 x 80 cm.

Cat. No. 27921

#### GIMA 2 instrument table

| cart with aluminium body;

2 white shelves;

2 side handles;

dimensions: 60 x 40 x 77;

weight: 15 kg.

Cat. No. 27429

#### MAYO instrument table

| stainless steel tray, on five wheels;

| dimensions: 68 x 46 cm;

- 5 wheels;
- | adjustable height from 70 to 100 cm;

weight: 9.3 kg.

Cat. No. 27453

#### Mayo-U-Shape instrument table

| stainless steel tray, on four wheels;

| internal dimensions: 75 x 45 cm;

adjustable height from 90 to 115 cm.

Nr kat. 45830

### Clinic equipment

#### INSTRUMENT TABLES

12









12

### Clinic equipment

#### INSTRUMENT TABLES

#### Instrument tables Medicazione

stainless steel cart;

- | 1 handle;
- | 2 shelves with 3.5 cm bends on four sides;
- dimensions: 60 x 40 x 81 cm;
- weight: 9 kg.

Cat. No. 27433



### Instrument tables INOX

steel cart;

- 4 plastic drawers;
- 2 shelves with 3.5 cm bends on four sides;
- dimensions: 70 x 50 x 82 cm.

Cat. No. 27430



#### TER Instrument table

| stainless steel on four rubber wheels;

- dimensions: 70 x 45 x 82 cm;
- 3 shelves;
- | weight 13 kg.

Cat. No. 27431

#### Medical cabine

- double door, locked;
- gray;
- | 32 shelves (16 on each side);
- | 5 internal shelves.

Cat. No. 27915



#### Medical cabinet

- enamelled design;
- | tempered glass;
- l door;
- 3 adjustable shelves;
- dimensions: 53 x 36 x 144 cm;
- weight: 26 kg.

Cat. No. 27901

#### Medical cabinet

enamelled design;

- | tempered glass;
- 4 door;
- 3 adjustable shelves;
- | dimensions: 75 x 38 x 169 cm;
- weight: 42 kg.

Cat. No. 27903

#### Medical cabinet


enamelled design;

- tempered glass;
- 2 door;
- 4 adjustable shelves;
- dimensions: 75 x 38 x 169 cm;

weight: 49 kg.

Cat. No. 27905

### Clinic equipment



12







STOOLS

12

#### Ergonomically designed comfortable stool

Ergonomically designed comfortable stool with saddleshaped seat and height adjustment with a gas piston. Height adjustable 52 to 67 cm.

Cat. No.	Name
45250	Black
45252	Blue
45251	White



#### Height adjustable stool with backrest

padded simulated leather seat Ø 33 cm with backrest
height adjustable 53 to 68 cm.

Cat. No.	
45240	black
45241	white
45242	blue
45243	green
45244	belge

#### Height adjustable stool

| Padded simulated leather seat  $\phi$  33 cm with backrest | Height adjustable 53 to 68 cm

Cat. No.	Name
45230	Black
45231	White
45232	Blue
45233	Blue lagoon
45234	Beige

#### Stainless steel stool

| plastic feet; | height adjustment: 47-59 cm.

#### 27950

Stainless steel AISI 304 seat Ø 32 cm.Height adjustable: 47-59 cm. Plastic base with feet.

#### 27951

Stainless steel AISI 304 seat  $\phi$  32 cm.Height adjustable: 51-63 cm. Plastic base with 5 castors.

	Cat. No.	Name
	27950	Steel stool
	27951	Steel stool on wheels



#### Ergo stool with backrest

ergonomically designed comfortable stool with saddleshaped seat and height adjustment with a gas piston;

height adjustable 62 to 78 cm.

Cat. No.	
45260	black
45262	blue
45261	white

## Clinic equipment



STOOLS





#### V Top table

stainless steel table;

- | made of two 762 mm panels arranged into a V-shape;
- | max height: 1041.4 mm;
- | min. height: 774.7 mm;
- | maximum load: 136 kg;
- option of heating the top form 21.2°C to 36.2°C;
- rotatable by 360;
- rails on both sides;
- option of mounting wheels for easier transportation.

Cat. No. 280919



#### Flat Top table

- stainless steel table;
- | max height: 996.95 mm;
- | min. height: 730.25 mm;
- | maximum load: 136 kg;
- | rounded corners;
- rotatable by 360 and almost upright positioning;
- rail on all sides;
- option of heating the top form 21.2°C to 36.2°C;
- option of mounting wheels for easier transportation;
- | perfect positioning.

Cat. No. 280920



#### Stół chirurgiczny Vet-Tables Dispomed

- the electric column table offers stability with a minimum footprint;
- | 5 year warranty;
- | manufactured in Canada by Dispomed;
- | lift capacity of 300 pounds (136 Kg);
- | all tray tops are made of 304 grade stainless steel;
- | maximum height 49" (124.46 cm);
- minimum height of 31.25" (78.74 cm);
- silent motor;
- | electric components are water resistant;
- | manual tilt;
- | weight: 245 lbs (111.13 Kg);
- base size : Length 23,5" (59,69 cm) / Width 20,375" (51,75 cm);
- | available in different configurations

Cat. No.. Dispomed10-14


#### Vet-Tables Scissor Exam and Surgery Table Dispomed

- | lift capacity 300 pounds (136 Kg);
- Canadian made;
- | rapid rail tie down system lets you safely position animals;
- the mobile parts of the table are on bearings;
- | all tray tops are made of 304 grade stainless steel;
- | minimum height 11" (27.9 cm);
- | maximum height 40" (101.6 cm);
- slide cleaners prevent the accumulation of hair in the lifting rails;
- | full aluminium machined structure for a greater robustness;
- silent motor stops automatically when the maximum height is reached;
- base size : Length 42" (106,68cm) / Width 17,5"(43,18cm);
  5 year warranty.

Cat. No. Dispomed15-20







12

#### Moving table for lifting animals

| perfect for lifting wound animals and animals after surgery;

- saves time and strain;
- maximum load: 90.7 kg;
- | Top Steel with elevated edges;
- | long life and quiet operation;
- four dirt-resistant bumpers;
- | full charging in 240 minutes;
- | four wheels for easy manoeuvring;
- min. height: 273mm;
- | max. height: 950 mm.

Cat. No. 280921



#### VET Lift table

- | perfect diagnostic and surgical table;
- | functional, attractive design, and top quality materials at a low price;
- | 14 point of fastening the patient, located under the top;
- height adjusted with foot pedals;
- height from 54 cm to 100 cm;
- stable base;
- | stainless steel top, 130 x 60 cm;
- with drain;
- | carrying capacity: 120 kg;
- 230 V, 50 Hz;
- | with an angling system;
- | wheels with lock can be installed.

Cat. No. 280916



#### MultiVet table

- Compatible with most portable X-ray devices
- | Integrated tripod
- | Deck translucent
- | Self-centering cassette feeder
- The tripod and cassette tray are on the guides and cover more than 80% of the table length
- Overall height: 207 cm
- Height: 64cm, width: 60cm, length: 140cm

Cat. No. 80355

#### MultiVet2 table

- The X-ray tableMultiVet offers a wide range of application possibilities. Whether for vertical, horizontal or inclined X-ray images, the precise positioning tool of the MultiVet makes the adjustment of the X-ray generator very comfortable and easy.
- With the addditional module for X-ray examinations on the standing animal you will obtain optimum examination results during abdominal screenings.
- Same as the CombiVet S table, the MultiVet X-ray table has a built in cassette drawer that moves simultaneously with the X-ray unit along the table.

Cat. No. 80399

#### Operation table, electric

| Operation table lifting and tilting

- | 120 x 48 cm table top with rims and drain hole
- Adjustable height: 780-1120 mm

#### ACCESSORIES:

- | Instrument tray for operation table Code: 40000149
- 2 standard + 2 wheels with break for operation table Code: 40000825
- Rubber mat for table Code: 40000017
- Infusion holder Code: 45909370

Cat. No. 80300









#### CAGES

12

#### Plastic cage

precise design;

- made from durable plastic;
- | resistant to staining and disinfectants;
- easy assembly;
- | multiple cages can be connected;
- one cage can be split into two smaller ones;
- can be mounted on a platform;
- | two sizes: middle and large;
- | removable floor (optional).

Cat. No. 280923



#### Cat cage

- | available as a single unit;
- movable;
- 4 sizes;
- | trapdoor lockable in the open and closed position;
- | stainless steel or acrylic door;
- | laminate or acrylic back;
- | made of damp and abrasion resistant board;
- | can be additionally equipped with drawers etc.

Cat. No. 280925

#### Isolation cage

- | perfect for infected patients;
- easily cleaned;
- elevated floor;
- slanted draining system;
- | a second cage can be mounted on one base;
- | two-layer floor minimises contact with faeces.

Cat. No. 280924




#### Stainless steel cage

- | made form top-quality material;
- | stainless steel fixtures;
- very durable;
- no sharp edges;
- safe lock;
- | very durable hinge;
- | quick and easy removal of door for cleaning, without the use of tools;
- | easy and quick fitting;
- | multiple cages can be connected;
- aesthetic look;
- various sizes.

Cat. No. 280922







12 Clinic equipment

#### CAGE ACCESSORIES

#### Cage cover

provides isolation during treatment;
resistant to stains;
blue;
machine-washable;
4 sizes.

Cat. No. 280927



#### IV infusion hanger

| comfortable height adjustment with one hand;

| easily attachable to cage;

available in 2 sizes.

Cat. No. 280930

#### Infusion pump clip

easy assembly;

| stable support for the pump;

adjustable.

Cat. No. 280931

#### Bowl

easily cleaned;

| resistant to indentation;

| made of stainless steel;

easily mounted handles available.

Cat. No. 280928

#### <u>Clipboards</u>

made of stainless steel;can be attached to door;resistant to fractures and breaking.

Cat. No. 280929



#### Cage divider

made of stainless steel;easy removal.

Cat. No. 280932

## Clinic equipment

CAGE ACCESSORIES

12







### NOTATKI


#### CAGE ACCESSORIES

#### Cage panel

12

- | made of stainless steel;
- aesthetic look;
- | no sharp edges;
- | all panels are 717.55 mm deep.

Cat. No. 280933



#### Oxygen therapy door

isolate the cage;

- | facilitate the use of oxygen therapy in animals;
- | made of acrylic glass;
- thickness: 6.4 mm;
- | temperature adjustment by adding ice or hot water to the chamber;
- a thermo-higrometer mounted on the door, monitoring temperature and humidity without opening the door.

#### Cat. No. 280934

#### Mobile platform

| protects cages from dampness and drafts;

- | polyurethane shelf;
- wheels rotate by 360 $^{\circ}$ ;
- | with safety mechanisms;
- various sizes.

Cat. No. 280935





| protects cages from dampness and drafts;

adjusted feet;

various sizes.

Cat. No. 280936









_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	
_	





POLNET Sp. z o.o. i Wspólnicy Spółka Komandytowa

ul. Sowia 13 B 62-080 Tarnowo Podgórne, Poland tel.: + 48 61 657 67 00, fax: + 48 61 657 67 67 office@polnet.pl

Veterinary Department: + 48 61 657 67 34; + 48 61 657 67 21 weterynaria@polnet.pl, www. polnet-wet. pl

www.polnet.pl







